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Table 1. Sample Analysis Results											
	Milk Powder						Liquid Milk				
	#1				#2	#3	#4				#5 (mg/L)
	Detected (mg/Kg)	Added (mg/Kg)	Found (mg/Kg)	Recovery (%)	Detected (mg/Kg)	Detected (mg/Kg)	Detected (mg/L)	Added (mg/L)	Found (mg/L)	Recovery (%)	
Melamine	n.a.	40	35	88	262	67	1.8	40	3.6	90	22

Note: \* Injections were made for each sample

\*\* Found = Measured value of spiked sample - Measured value of sample

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Sample #	Melamine Concentration using the RP-PLC method (mg/Kg) <sup>1</sup>		Melamine Concentration using the Acclaim Mixed-Mode WCX-1 Column (mg/Kg) <sup>2</sup>	
	Prepared following the description in this AN	Prepared following the description in Ref. 5	Prepared following the description in this AN	Prepared following the description in Ref. 5
#1	n.a.	n.a.	n.a.	n.a.
#2	274	229	262	229
#3	67	68	67	68
#6	2.2	2.3	2.4	2.8
#7	n.a.	n.a.	n.a.	n.a.
#8	n.a.	n.a.	n.a.	n.a.

Note: 1. Using the chromatographic conditions described in Ref. 5

2. Using the chromatographic conditions described in this AN

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1. Filigenzi M.S., Tor E.R, Poppenga R.H., Aston L.A and Puschner B., The determination of melamine in muscle tissue by liquid chromatography/tandem mass spectrometry. *Rapid Commun. Mass Spectrom.* **2007**,21(24), 4027-4032.
2. Ding T., Xu J., Li J., Shen C., Wu B., Chen H. and Li S., Determination of melamine residue in

plant origin protein powders using high performance liquid chromatography-diode array detection and high performance liquid chromatography-electrospray ionization tandem mass spectrometry. *Se Pu*. **2008**, *26(1)*, 6-9.

3. Kim B, Perkins LB, Bushway RJ, Nesbit S, Fan T, Sheridan R, Greene V., Determination of melamine in pet food by enzyme immunoassay, high-performance liquid chromatography with diode array detection, and ultra-performance liquid chromatography with tandem mass spectrometry. *JAOAC Int*. **2008**, *91(2)*, 408-413.

4. Leo Wang, Stacy M. Henday, Xiaodong Liu, Mark Tracy and William C. Schnute, Simultaneous Determination of Melamine and Cyanuric Acid Using LC-MS with the Acclaim Mixed-Mode WAX-1 Column and Mass Spectrometric Detection. Dionex Corporation LPN 1991-01 10/07.

5. Determination of melamine in raw milk and dairy products, GB/T 22388 2008.

6. Rapid determination of melamine in raw milk (High performance liquid chromatography method), GB/T22400 2008.

7. Dionex Corporation. *Product Manual for the Acclaim Mixed-Mode WCX-1 Column*.