
HAYER & BOECKER

DIE DRAHTWEBER

HAYER ZD-T25

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1. Introduction

Only install and start the test sieve shaker after having read the following operating instructions.

The operating instructions give important information about the secure, appropriate and economic operation of the test sieve shaker. They should be respected in order to avoid any dangers as well as to maximize the reliability and durability of the machine.

The safety instructions and regulations of prevention of accidents, as well as the different protective recommendations of the technical associations and employer's liability insurance associations, have to be strictly observed when using the machine.

**All safety instructions within these operating instructions have to be followed!
Always keep the operating instructions with the machine.**

If used properly, the construction of the machine guarantees that it runs and can be operated without any risk to persons.

The guarantee does not apply to damages caused by improper use due to a non-respect of the operating instructions.

All pictures and drawings in these operating instructions are general illustrations. They do not reflect the construction in detail.

The technical data and dimensions are without obligation, claims cannot be derived.

We reserve technical changes of the machine described below, without changing the operating instructions.

All rights on these operating instructions, the added drawings and other documents remain to HAVER & BOECKER. So are the rights on disposition, copying and passing on, including trademark registering.

2. Transport

The HAVER ZD-T25 is delivered completely assembled, except the guide rods.

If not otherwise agreed, it is packed according to the packing directives „HPE“, as defined by „Bundesverband Holzmittel, Paletten, Exportverpackungen e.V.“ (the national association of wooden instruments, pallets, export packaging) and „Verein Deutscher Maschinenbauanstalten“ (the society of German engineering shops)

Observe the current valid national and international directives for accident prevention. Damage caused by improper transport will invalidate any guarantee or replacement.

Warning


Improper lifting may cause bodily harm or material damage. Only lift the machine with the suitable equipment and by qualified persons!


Do not lift the machine with the guide rods!

3. Quantity delivered


3.1 HAVER ZD-T25

*Basic machine HAVER ZD-T25
only use for dry sieving*


cover with inspection glass

1 pair of guide rods, clamping and counter nuts


*1 spanner SW 27 DIN 894 WGB
to fix the guide rods*

1 connection cable, 3 m long, with EURO-plug and appliance plug

3.2 Clamping systems (according to model)

3.2.1 Classic

1. Nut with handlebars
2. Sieve cover
3. Guide rod


Elevation adjustment:

By turning the nut it moves up and down the guide rod.

Tensioning:

The sieve cover has to be installed and both nuts must be screwed down **simultaneously**.

3.2.2 HAVER TwinNut

1. Nut HAVER TwinNut
2. Guide rod
3. Sieve cover

Elevation adjustment:

Pull apart the two parts of the nut and slip it on the guide rod up to the sieve cover. Reclose it.

Tensioning:

Tighten the two HAVER TwinNuts **simultaneously**.

4. Technical data for HAVER ZD-T25

Test sieve shaker with horizontal rotary sieve motion

Operation voltage:	230 Volt or 110 Volt with transformer
Frequency:	50 Hertz / 60 Hertz
Power consumption:	0,18 kW
Max. rotation speed:	300 rpm
Amplitude:	Circular motion with 10 mm radius
Control:	all-electronic, rotary speed, rotation angle, interval and pause time can be set and be stored as sieve program
Weight:	89 kg, without test sieves
Test sieves:	Ø 250 mm
Max. weight of full sieve tower:	15 kg
Max. weight of sieve material:	5 kg
Max. number of sieves:	9
Dimensions:	W 575 x D 583 x H 915 mm
Emitted sound:	The emitted sound of the machine in action has been determined with a comparative sound measurement according to DIN 45635-01 K11. The equivalent continuous sound level A is: $L_{pAeq} \leq 70$ (dB) The intensity of sound can be reduced when using HAVER Sound absorbing Cupboards.

5. Application references

The machine has been constructed according to the actual technical standards and the recognized safety specifications. Although, there might be dangers for the user or third persons or impairments of the machine and other goods.

Only use the machine in perfect technical constitution. Follow the operating instructions conscious of safety and danger. Repair defects at once by authorized and qualified persons!

The machine has exclusively been developed for grain sieving to determine particle distribution. Any other or further application is not as agreed upon. The manufacturer/supplier will assume no responsibility for damages as a result of improper use. The risk is on the user's side.



Only use the HAVER ZD-T25 for dry sieving

During the sieving analysis the movement of the sieving material may cause unbalance. It is recommended to place the machine on a very stable table or on a concrete floor. The machine has to be assured against fall.


6. Description HAVER ZD-T25

1. Nut Haver TwinNut
2. Inspection glass
3. Cover
4. Test sieve set
5. Guide rod
6. Vibration body
7. Rubber buffer
8. Front panel with control board
9. Sieve pan
10. Plug connection with fuse box
11. Counter nut
12. Connection cable


Never use the HAVER ZD-T25 in damp rooms! It is not suitable for wet sieving!

7. Assembly

1. Put the machine on a **very stable, vibration free table or on the floor (concrete ground)**. **Do not use a rubber support!** Placed on a table the machine has to be assured against fall.
2. **The machine has to be fixed at its position.** Uneven ground can be adjusted by turning the rubber buffers and then fixing them with the counter nut.
3. Screw the guide rods (3) into the drilled holes of the vibration body (1) and fix the counter nuts (2) by the spanner.

1. Put on the test sieve set (1) and fill the sieving material in the upper test sieve.
2. Put on the cover (2) and fix it **simultaneously** with the two tension nuts (3).

Please notice item 8, Test sieving according to ISO 2591-1!

Connect the mains plug, pay attention to the power supply!

8. Test Sieving according to ISO 2591-1

During the sieve analysis, a certain quantity of the test product is divided into fractions by one or several test sieves. The difficulty of the sieve analysis is to determine the separation limit of the separation procedure which is carried out by sieving. It is therefore recommended to determine the analytical separation limit.

Details of the basic principals for sieve analysis are included in the standard sheets.

- DIN 66 165, part 1 - Basic standard - edition 1987
- DIN 66 165, part 2 - Execution - edition 1987

These sheets are available at the Beuth-Verlag, Berlin.

8.1 Largest particle to be permitted on a test sieve


To avoid damage to the test sieve, the size of the largest particle in the charge should not exceed $10w^{0.7}$ mm, where w is the nominal aperture width in mm.

Example	
Nominal aperture size, w mm	Size of largest particle, apx. mm
0,045	1
0,25	4
1	10
4	25

9. Recommended bulk volume for test sieves with 200 mm Ø

1	2	3
Nominal aperture width, w	Bulk volume of material	
μm	Volume of charge app. cm^3	Max. volume of residue app. cm^3
25	22	11
32	26	13
45	30	15
63	35	17
90	42	21
125	50	25
180	60	30
250	70	35
355	80	40
500	100	50
710	120	60
mm		
1	140	70
1,4	160	80
2	200	100
2,8	240	120
4	350	175
5,6	400	200
8	500	250
11,2	800	400
16	1000	500
22,4	1600	800

10. Starting and operating

Before connecting with the power supply please verify that the motor voltage corresponds to your supply voltage.

Please notice the power supply!

Starting the sieving

1. Switch on the machine with the main switch - it shines red
2. The last set points appear with an indication in front of

Angle:	180 – 1800°
Speed:	150 – 300 U/min
Cycles:	1 – 10
Pause:	0 – 10 sec.

3. Choose the set value by the key ENTER, the value flashes.
4. By the left arrow key the value is set downwards, by the right arrow key upwards.
5. Press ENTER, the value is stored.
6. The other values can be chosen and set in the same way with the arrow keys up / down.
7. Press the key RUN / STOP to start the sieving cycle.
8. Press the key RUN / STOP again to interrupt the sieving.

ESC

Interruption of an action

SAVE

By pressing this key 10 set points (angle, speed, cycles, pause) can be stored and protected from power failure. Insert the set points as described above and press the key SAVE. The display shows 10 times P for program space. One of the letters flashes. Choose the storage place by the arrow keys UP / DOWN and RIGHT / LEFT. Press the key ENTER. The indication changes from P to P1 up to P10. The set values are stored.

LOAD

By pressing this key the storage place (P1 – P10) can be selected. Press the key LOAD. Choose the storage place by the arrow keys UP / DOWN and RIGHT / LEFT. Press the key RUN / STOP to start the sieving cycle.

12. Guarantee

All functions of the sieving machine are checked and adjusted before its delivery.

The machine guarantee is for **24 months**.

Test sieves are excluded from guarantee!

We are only liable for defects resulting from flaws or faults in material. The guarantee does not apply to damages caused by transport. Repairs resulting by improper handling or wear due to operation are excluded from guarantee. The guarantee is void, if interventions such as variations or improper repairs are effected which are outside our control.

13. Maintenance and repair

Please note that HAVER Test Sieve Shakers only should be repaired, serviced and checked by authorized electro-experts following the national standards (e.f. according to VBG 4). Improper repairs may result in considerable dangers for the user.

If a part of the machine and its equipment is changed or replaced, this part has to be checked again according to EN 60 204-1!

In any case of failure we recommend to contact our service:


HAVER & BOECKER
PARTICLE SIZE ANALYSIS

Ennigerloher Straße 64
59302 OELDE, Germany
Phone: 02522-30363 (GERMANY)
+49-2522-30330 (ABROAD)
Fax: +49-2522-30404
E-Mail: PA@haverboecker.com

The HAVER HAVER ZD-T25 should be checked and serviced at least every 2 years at our works!

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Cigarette—Determination of structure of blending cut tobacco

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