Application Note

UV VIS

Fields of Application / Industry:

- Chemistry / Polymer Industry
- Clinical Chemistry / Medicine / Hygiene / Health Care
- Cosmetics
- Electronics
- Energy
- Environment / Water / Waste
- Food / Agriculture
- Geology / Mining
- **Material Analysis**
- Metallurgy / Galvanization
- Pharmacy
- Refineries / Petrochemistry
- Semi-Conductor Technology
- Others

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Color location determination of glass tube samples

General

Three glass tube samples of sort FIOLAX® (1x tube, 1x flattened tube segment and 1x tube plate) have been analyzed concerning their color. The glass tube was additionally analyzed as half shell.

Because of their chemical resistance, neutrality, temperature resistance and stability this sort of glass is mainly used for package for products of biotechnology and pharmacy. Furthermore the tube shape is resistant against high inner pressure, so it is mainly used for ampullae and tubes for pills and accommodates protection concerning quality assurance of drugs.

For determination of color location the SPECORD® 200 PLUS with solid sample holder was used.

The aim was to test the applicability of the instrument and if the bend of the samples can be disregarded when sample is positioned directly in front of detector.

For the color location determination, VIS transmittance spectra were scanned and then evaluated by the color determination of WinASPECT PLUS software. This module offers the possibility of calculating different color indexes with consideration of different angles of incidence and light models.

The chromaticity coordinates x and y characterize the color in a two-dimensional Cartesian space; in this respect, the saturation must be specified and given separately. The CIE coordinates characterize the color in a three-dimensional Cartesian space, while taking the brightness into account. The brightness of a color tone is represented by the value "L" in the CIE lab system. In this respect, the value "a" stands for the red-green portion and the value "b" for the yellow-blue portion.

Solid sample holder

Figure 1 shows the holder for solid samples. This accessory is suitable for solid samples with a diameter from 20mm and a thickness up to 25mm and is placed on the mount for cells of SPECORD directly in front of detector.

Fig. 1 Solid sample holder for SPECORD®

The measurements were carried out with the following parameter settings:

| | | | |
|--------------------------|------------------------------|---------------------|-------------|
| Title: | Color location determination | Range: | 190 -1100nm |
| Correction: | Reference | Step width: | 1nm |
| Display: | Transmittance | Speed: | 50nm/s |
| Lamp changeover: | 320nm | Accessories: | none |
| Measurement mode: | Scan mode | | |

The reference measurement was conducted against air. For sample measurement the glasses were placed into the solid holder and fixed with the clip.

Results/Evaluation

The figure 2 shows the transmittance spectra of all glass samples. The accordant color values are given below. The calculated CIE L*a*b* coordinates base on an light model C (most corresponding to natural light) and an observation angle of 2°.


Fig. 2: Transmittance spectra of the glass samples

Color coordinates of glass samples:

Glass plate

X= 20.6493 Y= 17.0683 Z= 2.5527
 x= 0.5128 y= 0.4238
 L= 48.3459 a= 20.1055 b= 55.2485
 Cab= 58.7931 hab= 1.2218

Tube

X= 43.7184 Y= 43.5559 Z= 20.4771
 x= 0.4057 y= 0.4042
 L= 71.9307 a= 2.9439 b= 40.1180
 Cab= 40.2258 hab= 1.4975

flattened tube segment

X= 44.6442 Y= 43.1233 Z= 22.3619
 x= 0.4054 y= 0.3916
 L= 71.6386 a= 6.8803 b= 36.2937
 Cab= 36.9401 hab= 1.3834

Tube half shell

X= 63.3266 Y= 64.7305 Z= 47.6702
 x= 0.3604 y= 0.3684
 L= 84.3447 a= -0.3512 b= 25.2509
 Cab= 25.2533 hab= -1.5569

The described method with the SPECORD® 200 PLUS with solid sample holder and was ideally suited for determining the color location of such samples. Taken spectra of glass plate, tube and tube half shell are different distinctly. Transmittance spectra and color coordinates of the tube and the flattened tube segment are conditionally comparable. The basic color coordinates X; Y and Z as well as the chromaticity coordinates x and y with the brightness value L are in a good accordance. of a color tone is represented by the value "L" in the CIE lab system. There are only differences in the red-green portion as well as the yellow-blue portion and the chroma.

The characterization of the color location can be effected for all color coordinates specified. In this respect, the color discrimination is always guaranteed. The bend of the tube effects the color coordinates. Hence for color determination this effect can not be disregarded although the sample was positioned directly in front of detector.

Printout and further use permitted with a reference to the source.

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© 2010 Analytik Jena AG
Analytik Jena AG