



Vacuum Pumps for moist Gases



Series LABOPORT® SD N 820.3 FT.40.18 Vacuum Pump

LABOPORT® SD Self-drying Vacuum Pump for moist Gases

Technical features:

- Pure transferring and evacuation
- Highly compatible with vapours and condensation
- Chemically-resistant
- Therefore suitable for highly aggressive or corrosive gases and vapours
- Maintenance-free
- Environmentally friendly
- Gastight, leakage rate approx. 6×10^{-3} mbar x l/s, not tested in serial production.

The chemically-resistant series N 820.3 FT.40.18 diaphragm pump is a twin-head unit with an integrated KNF self-drying system. There is a wide range of applications for this pump in laboratories, especially whenever clean vacuum is required and moist gases must be pumped down. Examples include vacuum-drying of vacuum drying chambers (for drying or heat-treating substances and components) or steam sterilizers for sterilizing instruments, vessels, filters and textiles. Due to its high resistance to aggressive media, this pump can be used universally.

The heart of these very compact pumps are KNF structured diaphragms (PTFE-coated). These patented diaphragms were stress-optimized using the Finite Elements method. As a result, we were able to make the pumps smaller while increasing the service life of the diaphragm.

The KNF self-drying system allows condensed liquid to be blown out of the pump heads at high speed during evacuation. The vacuum in the recipient remains constant. The drying cycle can be adjusted to the requirements of the individual process using three variables. After drying, the pump reaches a better vacuum and is able to evacuate significantly faster compared with pumps without a drying system.

Material in contact with the pumped media

Type/Order No.	Pump head	Diaphragm	Valves
N 820.3 FT.40.18	PTFE	PTFE-coated	FFPM

Technical data:

	N 820.3 FT.40.18
Delivery (l/min) ¹⁾	20
Ultimate vacuum (mbar abs.)	10
Operating pressure (bar g)	1
Connectors for tube (mm)	ID 10
Permissible gas and ambient temperature	+5...+40 °C
Mains	230V/50Hz
Motor protection	IP 44
Power P ₁	120 W
Operating current	0.7 A
Weight	9.6 kg
Dimensions	
LxHxW (mm)	312/220/177
With thermal switch and power fuse	

Motors with other voltages and frequencies on request.

¹⁾ at atm. pressure

Spare parts


Description	Order No.
Spares kit	057358




Vacuum Pumps for moist Gases

Dimensions and performance characteristics


Pump down time for 10 l receiver



Performance characteristics



Dimensions (mm)





Vacuum Pumps for moist Gases



Series LABOPORT® SD N 840.3 FT.40.18 Vacuum Pump

The chemically-resistant series N 840.3 FT.40.18 diaphragm pump is a twin-head unit with an integrated KNF self-drying system. There is a wide range of applications for this pump in laboratories, especially whenever clean vacuum is required and moist gases must be pumped down. Examples include vacuum-drying of vacuum drying chambers (for drying or heat-treating substances and components) or steam sterilizers for sterilizing instruments, vessels, filters and textiles. Due to its high resistance to aggressive media, this pump can be used universally.

The heart of these very compact pumps are KNF structured diaphragms (PTFE-coated). These patented diaphragms were stress-optimized using the Finite Elements method. As a result, we were able to make the pumps smaller while increasing the service life of the diaphragm.

The KNF self-drying system allows condensed liquid to be blown out of the pump heads at high speed during evacuation. The vacuum in the recipient remains constant. The drying cycle can be adjusted to the requirements of the individual process using three variables. After drying, the pump reaches a better vacuum and is able to evacuate significantly faster compared with pumps without a drying system.

LABOPORT® SD Self-drying Vacuum Pump for moist Gases

Technical features:

- Pure transferring and evacuation
- Highly compatible with vapours and condensation
- Chemically-resistant
- Therefore suitable for highly aggressive or corrosive gases and vapours
- Maintenance-free
- Environmentally friendly
- Gastight, leakage rate approx. 6×10^{-3} mbar x l/s, not tested in serial production.

Material in contact with the pumped media

Type/Order No.	Pump head	Diaphragm	Valves
N 840.3 FT.40.18	PTFE	PTFE-coated	FFPM

Technical data:

	N 840.3 FT.40.18
Delivery (l/min) ¹⁾	34
Ultimate vacuum (mbar abs.)	10
Operating pressure (bar g)	1
Connectors for tube (mm)	ID 10
Permissible gas and ambient temperature	+5...+40 °C
Mains	230V/50Hz
Motor protection	IP 44
Power P ₁	245 W
Operating current	1.5 A
Weight	12.9 kg
Dimensions	
LxHxW (mm)	341/239/189
With thermal switch and power fuse	

Motors with other voltages and frequencies on request.

¹⁾ at atm. pressure

Spare parts


Description	Order No.
Spares kit	057359




Vacuum Pumps for moist Gases

Dimensions and performance characteristics


Pump down time for 10 l receiver



Performance characteristics



Dimensions (mm)





The KNF-Group worldwide

Germany

KNF Neuberger GmbH
Alter Weg 3
D-79112 Freiburg
Tel. +49 7664 5909 0
Fax +49 7664 5909 99
info@knf.de
www.knf.de

Benelux Netherlands

KNF Verder B.V.
Utrechtseweg 4a
NL-3451 GG Vleuten
Tel. +31 30 677 92 40
Fax +31 30 677 92 47
info@knf-verder.nl
www.knf-verder.nl

Benelux Belgium and Luxembourg

KNF Verder N.V.
Kontichsesteenweg 17
B-2630 Aartselaar
Tel. +32 3 871 96 24
Fax +32 3 871 96 28
info@knf.be
www.knf.be

China

KNF Neuberger Trading
(Shanghai) Co., Ltd.
No. 36 Lane 1000
Zhang Heng Road
Shanghai 201203
Tel. +86 21 6859 6566
Fax +86 21 6859 6568
info@knf.com.cn
www.knf.com.cn

France,

Marocco, Tunisia, Algeria

KNF Neuberger S.A.
4, Bld. d'Alsace Z.I.
F-68128 Village-Neuf
Tel. +33 389 70 35 00
Fax +33 389 69 92 52
info@knf.fr
www.knf.fr

United Kingdom

KNF Neuberger U.K., Ltd.
Avenue 2, Station Lane
Industrial Estate
Witney, Oxfordshire OX28 4FA
Tel. +44 1993 77 83 73
Fax +44 1993 77 51 48
info@knf.co.uk
www.knf.co.uk

India

KNF Pumps + Systems (India) Pvt. Ltd.
RAJIV GANDHI INFOTECH PARK,
Phase 1,
Ganga Estate, Survey No. 152/2/2,
Above AXIS BANK,
Hinjewadi, Pune 411 057
Tel. +91 20 640 13 923
+91 20 640 08 923
Fax +91 20 229 33 923
info@knfpumps.in
www.knfpumps.in

Italy

KNF ITALIA S.r.l.
Via Flumendosa, 10
I-20132 Milano
Tel. +39 02 272 03 860
Fax +39 02 272 03 848
info@knf.it
www.knf.it

Japan

KNF Japan Co.Ltd.
Across Shinkawa Bldg. Annex 3F,
1-16-14 Shinkawa, Chuo-ku,
Tokyo 104-0033
Tel. +81 3 3551 7931
Fax +81 3 3551 7932
info@knf.co.jp
www.knf.co.jp

Korea

KNF Neuberger Ltd.
Woosan Bldg.Rm#202
336-4, Hwikyung-Dong
Dongdaemun-Ku
130-090, Seoul
Tel. +82 2 959 0255
Fax +82 2 959 0254
knf@knfkorea.com
www.knfkorea.com

Sweden, Finland, Denmark, Norway

KNF Neuberger AB
Mejerivägen 4
P.O. Box 44060
SE-10073 Stockholm
Tel. +46 8 744 51 13
Fax +46 8 744 51 17
info@knf.se
www.knf.se

Switzerland Sales

KNF Neuberger AG
Stockenstrasse 6
CH-8362 Balternwil
Tel. +41 71 973 99 30
Fax +41 71 973 99 31
knf@knf.ch
www.knf.ch

Switzerland

Production Liquid Pumps
KNF FLODOS AG
Wassermatte 2
CH-6210 Sursee
Tel. +41 41 925 00 25
Fax +41 41 925 00 37
info@knf-flodos.ch
www.knf-flodos.ch

Switzerland

Production Micro Pumps
KNF Micro AG
Zelglimatte 1b
CH-6260 Reiden
Tel. +41 62 787 88 88
Fax +41 62 787 88 99
info@knf-micro.ch
www.knf-micro.ch

Taiwan

KNF Neuberger Ltd.
9F-2, No. 24, Lane 123, Section 6,
Min Chuan East Road
Neihu District
Taipei City, 11490
Tel. +886 2 2794 1011
Fax +886 2 8792 1648
knftwn@knftwn.com.tw
www.knftwn.com.tw

USA, Canada, South America

KNF Neuberger, INC.
Two Black Forest Road
Trenton, New Jersey 08691-1810
Tel. +1 609 890 8600
Fax +1 609 890 8323
knfusa@knf.com
www.knfusa.com
Latin America
Tel. +1 609 649 1010
gb@knf.com

Diaphragm Pumps + Systems



KNF Neuberger GmbH



KNF Neuberger, INC.



KNF Micro AG



KNF Neuberger S.A.



KNF FLODOS AG

The KNF-Production worldwide

Germany

KNF Neuberger GmbH
Freiburg-Munzingen

France

KNF Neuberger S.A.
Village-Neuf

Switzerland

KNF FLODOS AG
Sursee

Switzerland

KNF Micro AG
Reiden

USA

KNF Neuberger, INC.
Trenton - New Jersey



Vacuum Pumps for moist Gases



Series LABOPORT® SD N 842.3 FT.40.18 Vacuum Pump

LABOPORT® SD Self-drying Vacuum Pump for moist Gases

Technical features:

- Pure transferring and evacuation
- Highly compatible with vapours and condensation
- Chemically-resistant
- Therefore suitable for highly aggressive or corrosive gases and vapours
- Maintenance-free
- Environmentally friendly
- Gastight, leakage rate approx. 6×10^{-3} mbar x l/s, not tested in serial production.

The chemically-resistant series N 842.3 FT.40.18 diaphragm pump is a twin-head unit with an integrated KNF self-drying system. There is a wide range of applications for this pump in laboratories, especially whenever clean vacuum is required and moist gases must be pumped down. Examples include vacuum-drying of vacuum drying chambers (for drying or heat-treating substances and components) or steam sterilizers for sterilizing instruments, vessels, filters and textiles. Due to its high resistance to aggressive media, this pump can be used universally.

The heart of these very compact pumps are KNF structured diaphragms (PTFE-coated). These patented diaphragms were stress-optimized using the Finite Elements method. As a result, we were able to make the pumps smaller while increasing the service life of the diaphragm.

The KNF self-drying system allows condensed liquid to be blown out of the pump heads at high speed during evacuation. The vacuum in the recipient remains constant. The drying cycle can be adjusted to the requirements of the individual process using three variables. After drying, the pump reaches a better vacuum and is able to evacuate significantly faster compared with pumps without a drying system.

Material in contact with the pumped media

Type/Order No.	Pump head	Diaphragm	Valves
N 842.3 FT.40.18	PTFE	PTFE-coated	FFPM

Technical data:

	N 842.3 FT.40.18
Delivery (l/min) ¹⁾	34
Ultimate vacuum (mbar abs.)	4
Operating pressure (bar g)	1
Connectors for tube (mm)	ID 10
Permissible gas and ambient temperature	+5...+40 °C
Mains	230V/50Hz
Motor protection	IP 44
Power P ₁	245 W
Operating current	1.5 A
Weight	13.7 kg
Dimensions	
LxHxW (mm)	341/242/189
With thermal switch and power fuse	

Motors with other voltages and frequencies on request.

¹⁾ at atm. pressure

Spare parts


Description	Order No.
Spares kit	057359




Vacuum Pumps for moist Gases

Dimensions and performance characteristics


Pump down time for 10 l receiver



Performance characteristics



Dimensions (mm)





The KNF-Group worldwide

Germany

KNF Neuberger GmbH
Alter Weg 3
D-79112 Freiburg
Tel. +49 7664 5909 0
Fax +49 7664 5909 99
info@knf.de
www.knf.de

Benelux Netherlands

KNF Verder B.V.
Utrechtseweg 4a
NL-3451 GG Vleuten
Tel. +31 30 677 92 40
Fax +31 30 677 92 47
info@knf-verder.nl
www.knf-verder.nl

Benelux Belgium and Luxembourg

KNF Verder N.V.
Kontichsesteenweg 17
B-2630 Aartselaar
Tel. +32 3 871 96 24
Fax +32 3 871 96 28
info@knf.be
www.knf.be

China

KNF Neuberger Trading
(Shanghai) Co., Ltd.
No. 36 Lane 1000
Zhang Heng Road
Shanghai 201203
Tel. +86 21 6859 6566
Fax +86 21 6859 6568
info@knf.com.cn
www.knf.com.cn

France,

Marocco, Tunisia, Algeria

KNF Neuberger S.A.
4, Bld. d'Alsace Z.I.
F-68128 Village-Neuf
Tel. +33 389 70 35 00
Fax +33 389 69 92 52
info@knf.fr
www.knf.fr

United Kingdom

KNF Neuberger U.K., Ltd.
Avenue 2, Station Lane
Industrial Estate
Witney, Oxfordshire OX28 4FA
Tel. +44 1993 77 83 73
Fax +44 1993 77 51 48
info@knf.co.uk
www.knf.co.uk

India

KNF Pumps + Systems (India) Pvt. Ltd.
RAJIV GANDHI INFOTECH PARK,
Phase 1,
Ganga Estate, Survey No. 152/2/2,
Above AXIS BANK,
Hinjewadi, Pune 411 057
Tel. +91 20 640 13 923
+91 20 640 08 923
Fax +91 20 229 33 923
info@knfpumps.in
www.knfpumps.in

Italy

KNF ITALIA S.r.l.
Via Flumendosa, 10
I-20132 Milano
Tel. +39 02 272 03 860
Fax +39 02 272 03 848
info@knf.it
www.knf.it

Japan

KNF Japan Co.Ltd.
Across Shinkawa Bldg. Annex 3F,
1-16-14 Shinkawa, Chuo-ku,
Tokyo 104-0033
Tel. +81 3 3551 7931
Fax +81 3 3551 7932
info@knf.co.jp
www.knf.co.jp

Korea

KNF Neuberger Ltd.
Woosan Bldg.Rm#202
336-4, Hwikyung-Dong
Dongdaemun-Ku
130-090, Seoul
Tel. +82 2 959 0255
Fax +82 2 959 0254
knf@knfkorea.com
www.knfkorea.com

Sweden, Finland, Denmark, Norway

KNF Neuberger AB
Mejerivägen 4
P.O. Box 44060
SE-10073 Stockholm
Tel. +46 8 744 51 13
Fax +46 8 744 51 17
info@knf.se
www.knf.se

Switzerland Sales

KNF Neuberger AG
Stockenstrasse 6
CH-8362 Balternwil
Tel. +41 71 973 99 30
Fax +41 71 973 99 31
knf@knf.ch
www.knf.ch

Switzerland

Production Liquid Pumps
KNF FLODOS AG
Wassermatte 2
CH-6210 Sursee
Tel. +41 41 925 00 25
Fax +41 41 925 00 37
info@knf-flodos.ch
www.knf-flodos.ch

Switzerland

Production Micro Pumps
KNF Micro AG
Zelglimatte 1b
CH-6260 Reiden
Tel. +41 62 787 88 88
Fax +41 62 787 88 99
info@knf-micro.ch
www.knf-micro.ch

Taiwan

KNF Neuberger Ltd.
9F-2, No. 24, Lane 123, Section 6,
Min Chuan East Road
Neihu District
Taipei City, 11490
Tel. +886 2 2794 1011
Fax +886 2 8792 1648
knftwn@knftwn.com.tw
www.knftwn.com.tw

USA, Canada, South America

KNF Neuberger, INC.
Two Black Forest Road
Trenton, New Jersey 08691-1810
Tel. +1 609 890 8600
Fax +1 609 890 8323
knfusa@knf.com
www.knfusa.com
Latin America
Tel. +1 609 649 1010
gb@knf.com

Diaphragm Pumps + Systems



KNF Neuberger GmbH



KNF Neuberger, INC.



KNF Micro AG



KNF Neuberger S.A.



KNF FLODOS AG

The KNF-Production worldwide

Germany

KNF Neuberger GmbH
Freiburg-Munzingen

France

KNF Neuberger S.A.
Village-Neuf

Switzerland

KNF FLODOS AG
Sursee

Switzerland

KNF Micro AG
Reiden

USA

KNF Neuberger, INC.
Trenton - New Jersey



Vacuum Pumps for moist Gases



Self drying Vacuum Pump for moist Gases

Technical features:

- Pure transferring and evacuation
- Highly compatible with vapours and condensation
- Chemically-resistant
- Therefore suitable for highly aggressive or corrosive gases and vapours
- Maintenance-free
- Environmentally friendly
- Gastight, leakage rate approx. 6×10^{-3} mbar x l/s, not tested in serial production.

Series N 860.3 FT.40.18 Vacuum Pump

The chemically-resistant series N 860.3 FT.40.18 diaphragm pump is a twin-head unit with an integrated KNF self-drying system. There is a wide range of applications for this pump in laboratories, especially whenever clean vacuum is required and moist gases must be pumped down. Examples include vacuum-drying of vacuum drying chambers (for drying or heat-treating substances and components) or steam sterilizers for sterilizing instruments, vessels, filters and textiles. Due to its high resistance to aggressive media, this pump can be used universally.

The heart of these very compact pumps are KNF structured diaphragms (PTFE-coated). These patented diaphragms were stress-optimized using the Finite Elements method. As a result, we were able to make the pumps smaller while increasing the service life of the diaphragm.

The KNF self-drying system allows condensed liquid to be blown out of the pump heads at high speed during evacuation. The vacuum in the recipient remains constant. The drying cycle can be adjusted to the requirements of the individual process using three variables. After drying, the pump reaches a better vacuum and is able to evacuate significantly faster compared with pumps without a drying system.

Material in contact with the pumped media

Type/Order No.	Pump head	Diaphragm	Valves
N 860.3 FT.40.18	PTFE	PTFE-coated	FFPM

Technical data:

	N 860.3 FT.40.18
Delivery (l/min) ¹⁾	60
Ultimate vacuum (mbar abs.)	4
Operating pressure (bar g)	1
Connectors for tube (mm)	ID 12
Permissible gas and ambient temperature	+5...+40 °C
Mains	230V/50Hz
Motor protection	IP 54
Power P ₁	220 W
Operating current	1.6 A
Weight	14.8 kg
Dimensions	
LxHxW (mm)	331/278.5/291
With thermal switch and power fuse	

Motors with other voltages and frequencies on request.

¹⁾ at atm. pressure

Spare parts


Description	Order No.
Spares kit	047499




Vacuum Pumps for moist Gases

Dimensions and performance characteristics


Pump down time for 20 l receiver



Performance characteristics



Dimensions (mm)





The KNF-Group worldwide

Germany

KNF Neuberger GmbH
Alter Weg 3
D-79112 Freiburg
Tel. +49 7664 5909 0
Fax +49 7664 5909 99
info@knf.de
www.knf.de

Benelux Netherlands

KNF Verder B.V.
Utrechtseweg 4a
NL-3451 GG Vleuten
Tel. +31 30 677 92 40
Fax +31 30 677 92 47
info@knf-verder.nl
www.knf-verder.nl

Benelux Belgium and Luxembourg

KNF Verder N.V.
Kontichsesteenweg 17
B-2630 Aartselaar
Tel. +32 3 871 96 24
Fax +32 3 871 96 28
info@knf.be
www.knf.be

China

KNF Neuberger Trading
(Shanghai) Co., Ltd.
No. 36 Lane 1000
Zhang Heng Road
Shanghai 201203
Tel. +86 21 6859 6566
Fax +86 21 6859 6568
info@knf.com.cn
www.knf.com.cn

France,

Marocco, Tunisia, Algeria

KNF Neuberger S.A.
4, Bld. d'Alsace Z.I.
F-68128 Village-Neuf
Tel. +33 389 70 35 00
Fax +33 389 69 92 52
info@knf.fr
www.knf.fr

United Kingdom

KNF Neuberger U.K., Ltd.
Avenue 2, Station Lane
Industrial Estate
Witney, Oxfordshire OX28 4FA
Tel. +44 1993 77 83 73
Fax +44 1993 77 51 48
info@knf.co.uk
www.knf.co.uk

India

KNF Pumps + Systems (India) Pvt. Ltd.
RAJIV GANDHI INFOTECH PARK,
Phase 1,
Ganga Estate, Survey No. 152/2/2,
Above AXIS BANK,
Hinjewadi, Pune 411 057
Tel. +91 20 640 13 923
+91 20 640 08 923
Fax +91 20 229 33 923
info@knfpumps.in
www.knfpumps.in

Italy

KNF ITALIA S.r.l.
Via Flumendosa, 10
I-20132 Milano
Tel. +39 02 272 03 860
Fax +39 02 272 03 848
info@knf.it
www.knf.it

Japan

KNF Japan Co.Ltd.
Across Shinkawa Bldg. Annex 3F,
1-16-14 Shinkawa, Chuo-ku,
Tokyo 104-0033
Tel. +81 3 3551 7931
Fax +81 3 3551 7932
info@knf.co.jp
www.knf.co.jp

Korea

KNF Neuberger Ltd.
Woosan Bldg.Rm#202
336-4, Hwikyung-Dong
Dongdaemun-Ku
130-090, Seoul
Tel. +82 2 959 0255
Fax +82 2 959 0254
knf@knfkorea.com
www.knfkorea.com

Sweden, Finland, Denmark, Norway

KNF Neuberger AB
Mejerivägen 4
P.O. Box 44060
SE-10073 Stockholm
Tel. +46 8 744 51 13
Fax +46 8 744 51 17
info@knf.se
www.knf.se

Switzerland Sales

KNF Neuberger AG
Stockenstrasse 6
CH-8362 Balterswil
Tel. +41 71 973 99 30
Fax +41 71 973 99 31
knf@knf.ch
www.knf.ch

Switzerland

Production Liquid Pumps
KNF FLODOS AG
Wassermatte 2
CH-6210 Sursee
Tel. +41 41 925 00 25
Fax +41 41 925 00 37
info@knf-flodos.ch
www.knf-flodos.ch

Switzerland

Production Micro Pumps
KNF Micro AG
Zelglimatte 1b
CH-6260 Reiden
Tel. +41 62 787 88 88
Fax +41 62 787 88 99
info@knf-micro.ch
www.knf-micro.ch

Taiwan

KNF Neuberger Ltd.
9F-2, No. 24, Lane 123, Section 6,
Min Chuan East Road
Neihu District
Taipei City, 11490
Tel. +886 2 2794 1011
Fax +886 2 8792 1648
knftwn@knftwn.com.tw
www.knftwn.com.tw

USA, Canada, South America

KNF Neuberger, INC.
Two Black Forest Road
Trenton, New Jersey 08691-1810
Tel. +1 609 890 8600
Fax +1 609 890 8323
knfusa@knf.com
www.knfusa.com
Latin America
Tel. +1 609 649 1010
gb@knf.com

Diaphragm Pumps + Systems



KNF Neuberger GmbH



KNF Neuberger, INC.



KNF Micro AG



KNF Neuberger S.A.



KNF FLODOS AG

The KNF-Production worldwide

Germany

KNF Neuberger GmbH
Freiburg-Munzingen

France

KNF Neuberger S.A.
Village-Neuf

Switzerland

KNF FLODOS AG
Sursee

Switzerland

KNF Micro AG
Reiden

USA

KNF Neuberger, INC.
Trenton - New Jersey