Software/hardware for e-beam lithography

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
Nanomaker system can be delivered both complete and module wise in accordance with the requirements:

- - available
- - emulation

NanoMaker
Complete deliverable

NanoMaker Workbench
To prepare data on offline PC

NanoMaker Writer
To design, perform exposure and image acquisition

	NanoMaker Complete deliverable	NanoMaker Workbench To prepare data on offline PC	NanoMaker Writer To design, perform exposure and image acquisition
Graphics Editor (specialized CAD system)	■	■	■
Import of lithographic structures and images from: *.DXF, *.CSF, *.GDS, TIFF Image Files (*.TIF), Bitmap Image Files (*. BMP)	■	■	■
Export of lithographic structures and images to: *.DXF, *.CSF, *.GDS, *. ELM	■	■	■
Recommended Parameters Reference Table (specialized Database)	■	■	■
Postprocessing (Negative, Union, Frame, Shrink, Erase, Overlaps Out, Dividing, Dose Stratification)	■	■	■
Proximity Effect Correction (including 3D structures)	■	■	■
Resist Development Simulation	■	■	■
Dose Curves using (for 3D structuring)	■	■	■
Exposure	■	■	■
Image Acquisition and Processing	■	■	■
Stitching and Alignment of Exposure fields	■	■	■
Job Processing	■	■	■
Active Compensation for Distortion and Dynamic delays	■	■	■
User Guide in electronic form	■	■	■
Pattern Generator + Drivers	■	■	■
Stage and Microscope Control Drivers (optional)	■	■	■

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
Output XY DACs and input ADC voltages can be tuned for arbitrary intervals in ± 10.0 V range.


Time characteristics of the Pattern Generator

Minimal exposure dwell time	50 ns (20 MHz)
Dwell Time resolution	10 ns (100 MHz)
On-board image accumulation (per accumulation cycle)	50 ns (20 MHz)
Minimal image acquisition time (time of XY DACs addressing + transferring of ADC signal)	0.8 μ s (implying 12 accumulation cycles)


Alignment for successive e-beam lithographies

Contact pads were made by photolithography then e-beam lithography was used to make mesa structure (rings). Another e-beam lithography layer was then incorporated to produce the metal contacts.

3D Proximity Correction and 3D Structuring

Kinofilm optical device produced by 3D lithography. By using 3D proximity correction, objects with arbitrary 3D shape can be created during a single exposure.

Large scale device production in combination with high resolution interferometer stage


NanoMaker's unique features of exposure fields stitching based on dynamic compensation of beam position in combination with an interferometer stage, demonstrate excellent results in production of large-scale devices (up to a few centimeters).

A fragment of device for X-ray optics and spectroscopy (by the courtesy of A. Firsov, BESSY II)


Rainbow holograms (OVD) writing with a SEM

Kinematic rainbow hologram produced at Bessy II using LEO 1560 SEM controlled by NanoMaker (by courtesy of A. Firsov). 400 nm PMMA resist covered by 50 nm of AL. Sample size is 2 cm by 1 cm

Gas Assisted Focused Electron Beam Induced Deposition

The NanoMaker can be successfully used for controlled deposition of thin film coatings on substrates by decomposing of metal-organic vapors induced by focused electron beam irradiation (so called Gas Assisted EBL). Photo demonstrates deposited thin Pt lines of 10, 20 and 30 nm width. By courtesy of Dr. L. Rotkina, UPenn.

Ultra High Resolution

This smallest world map produced by EBL. It has been produced by Proxy-Writer (the MS DOS predecessor of NanoMaker) in connection with a scanning electron microscope (UHV FE-STEM HB 501). The substrate was a SiN membrane covered with the electron sensitive material AIF3.