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Suitable for: Vulcanized rubber, Non-foam plastic, Carbon true density, True density of rock or coal, Pitch coke true density, Abrasive materials true density, Pigment density, Solid and half-solid asphalt density, Cement powder, Refractory materials, Ceramic, Natural stone, Water-coal-slurry density, Carbon materials true density by using boiling method, Soil true density, Ion resin wet true density, Powder true density research lab.

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Principle: According to the standards of GB/T533-A, GB/T1033.1, GB/T6155, GB/T23561, YB/T5300, JB/T7984.3, GB/T 1713, GB/T8929, GB/T1713, GB/T208, GB/T5071, QB/T1010, GB/T9966, GB/T18856, GB/T24203, GB/T8330, SL-237, and by adopting the immersed volume displacement method of Archimedean principle, and cooperated with pycnometer, it can show the measuring result directly.

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Notices:

- Under immersed method, choosing the liquid which can not dissolve but can moisture the sample surface is important.
- For ceramic materials, such as feldspar, quartz and ceramic products, distilled water can be used as liquid medium.
- For cement material, kerosene, xylene or other organic liquid can be used as testing medium.
- For inorganic powder, organic liquid also can be chosen.
- By using porcelain mortar to rub powder and filter it by 240 mesh sieve, then put the powder in the weighing bottle, and dry it by moisture meter under 105°C, take out and cool it to room temperature in desiccator.