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## LC-3060 Minisize HPLC Pump (10mL/50mL/100mL)

### Core Technology

- ARM7 core 32-bit embedded microprocessor;
- Electronic pressure pulsation suppression technology;
- High-precision DC servo motor.

### Application

- Solvent delivery pump for analytical HPLC system
- Solvent delivery pump for semi-preparative HPLC system
- Sampling pump for preparative HPLC system
- Solvent delivery pump for ion chromatography
- Other demand for high pressure and stable solvent delivery

### Features

- Ultra-compact design
- Expandable to binary high pressure gradient system
- Expandable to quaternary low pressure gradient system
- No mechanical buffer (electronic pressure pulsation suppression)
- Pump head changeable in accordance with requirement (analytical/semi-preparative)
- Changeable to Peek pump head when needed
- External control realized via RS-232/485 interface

### Other

- LC-3060 High Pressure Solvent Delivery pump adopts electronic pressure pulsation suppression technology displacing traditional mechanical buffer, efficiently controlling the fluctuation of the flow rate and further improving the reliability of the instrument, and at the same time reducing the dead volume of solvent delivery system to the minimum.
- High precision DC servo motor and encoder are used to provide strong guarantee to the implementation of electronic pulsation suppression technique. The volume/ weight ratio of the pump parts reduces a lot compared to that of the step motor. Noise has almost been eliminated.
- The transfer efficiency of mobile phase under different column dumping (On/off of check-valve; compressibility of mobile phase, etc..) is automatically adjusted and compensated by corresponding parameters, to ensure the stability of flow rate in the

chromatograph system.

- The flow range can be adjusted by changing the pump head and corresponding parameters, from 10ml analytical type to 50ml semi-preparative type.
- It can be controlled by chromatographic workstation via RS-232 port.

## Technical Specifications

### 10mL

<b>Flow Rate Range</b>	0.01~9.99 mL/min
<b>Delivery System</b>	Double-piston serial reciprocating pump
<b>Flow Rate Accuracy</b>	±0.06% (1mL/min, secondary water, room temperature)
<b>Flow Rate Reproducibility</b>	< 0.03% (1mL/min, secondary water, room temperature)
<b>Pressure Range</b>	0~43 MPa (1mL/min, secondary water, room temperature)
<b>Pressure Pulsation</b>	<0.1MPa (1mL/min, secondary water, 15MPa)
<b>Control</b>	RS-232/485 port, PC control
<b>System Protection</b>	Maximum pressure protection setting
<b>Display</b>	LCD 2×8 characters
<b>Power supply</b>	DC 24V, 3A
<b>Dimensions/Weight</b>	130×110×210 mm (W×H×D) /2.8Kg

### 50mL

<b>Flow Rate Range</b>	0.01~49.99 mL/min
<b>Delivery System</b>	Double-piston serial reciprocating pump
<b>Flow Rate Accuracy</b>	±0.2% (1mL/min, secondary water, room temperature)
<b>Flow Rate Reproducibility</b>	< 0.1% (1mL/min, secondary water, room temperature)
<b>Pressure Range</b>	0~32 MPa (1mL/min, secondary water, room temperature)
<b>Pressure Pulsation</b>	<0.1MPa (1mL/min, secondary water, 15MPa)
<b>Control</b>	RS-232/485 port, PC control
<b>System Protection</b>	Maximum pressure protection setting
<b>Display</b>	LCD 2×8 characters
<b>Power supply</b>	DC 24V, 3A
<b>Dimensions/Weight</b>	130×110×210 mm (W×H×D) /2.8Kg

### 100mL

<b>Flow Rate Range</b>	0.01~99.99 mL/min
<b>Delivery System</b>	Double-piston parallel reciprocating pump
<b>Flow Rate Accuracy</b>	±1% (10mL/min, secondary water, room temperature)
<b>Flow Rate Reproducibility</b>	< 0.5% (10mL/min, secondary water, room temperature)
<b>Pressure Range</b>	0~15 MPa (1mL/min, secondary water, room temperature)