用于制药产品密封件的清洗、硅化、灭菌、干燥设备  
Washing, siliconizing, sterilizing and drying machine for pharmaceutical industry

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外形尺寸(mm) Outer Dimension			篮框数 Basket NO.	有效容量 Useful Volume	料斗总容量 Total Volume	注射用水耗量 WFI Water Consumption	洁净蒸汽耗量 Saturated Steam Consumption	软化冷却水耗量 Cold Softened water Consumption	工业蒸汽耗量 Industrial Steam Consumption	压缩空气耗量 Instrument Compressed air Consumption	电源 Power		重量 Weight	
宽度 Width	深度 Depth	高度 Height	Set	Lt	Lt	Kg/cycle	Kg/cycle	Kg/cycle	Kg/cycle	N Lt	控制电源 Control Power AC220V/50HZ/KW	动力电源 Main Power 3-AC380V/50HZ/KW	Kg	
LST80	2800	1330	2600	8	80	240	1000	30	250	40	180	1	31	2370
LST120	3000	1550	2600	12	120	360	1200	60	300	70	180	1	35	3000
LST180	3200	1700	2800	18	180	540	1500	120	350	100	180	1	38	4000
LST240	3400	1930	3000	20	240	720	2000	200	400	120	180	1	43	4800

用于制药产品密封件的清洗、硅化、干燥设备  
Washing, siliconizing, and drying machine for pharmaceutical industry

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外形尺寸(mm) Outer Dimension			篮框数 Basket NO.	有效容量 Useful Volume	料斗总容量 Total Volume	注射用水耗量 WFI Water Consumption	工业蒸汽耗量 Industrial Steam Consumption	压缩空气耗量 Instrument Compressed air Consumption	电源 Power		重量 Weight	
宽度 Width	深度 Depth	高度 Height	Set	Lt	Lt	Kg/cycle	Kg/cycle	N Lt	控制电源 Control Power AC220V/50HZ/KW	动力电源 Main Power 3-AC380V/50HZ/KW	Kg	
LT80	2600	1330	2600	8	80	240	1000	40	180	1	31	2370
LT120	2600	1550	2600	12	120	360	1200	70	180	1	35	3000
LT180	2800	1700	2800	18	180	540	1500	100	180	1	38	3700
LT240	3000	1930	3000	20	240	720	2000	120	180	1	43	4300

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
Equipment is designed for fully treatment of sealing devices for pharmaceutical products, like rubber and plastic stoppers, droppers, metal rings and similar components.

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Unique Separated I Basket Design: Products are distributed into each separated basket by physical effect.

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
- Maximum reduction of friction:
  - \* Avoid over volume stack of products at each process.
  - \* Decrease mechanical strike of products.
- Maximum water saving, to reduce production cost.
- More even siliconizing: with special silicone oil measuring pump and filling system.
- Better sterilization: surfaces of all the rubber stopper products can be sterilized.

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
Products filled at each separated basket, so to achieve better treatment effect. Washing is realized by whole basket drum rotating (each individual basket is fixed on the hollow main shaft), the fluid bed effect is realized by circulation pump and spray nozzle so to eliminate the suspending particles.

- Sterilizing under temperature of 121°C
- Drying with 130°C hot air passed HEPA filters. (temperature can be set at HMI)
- Cooling by 20°C sterilized air passed HEPA filters.
- Loading and unloading door has reliable safety inter-lock devise.
- Basket door open and close automatically
- Process door conform to sterilizing requirements.
- Automatic loading.
- Unloading door connect with isolator.
- Multiple ways for automatic unloading
- Equipment has connections for validation tests

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- Machine is consisted of chamber, baskets, hollow main shaft, loading and unloading door, baskets auto open and close devise, circulation pump system, detergent measuring and filling system, silicone oil measuring and filling system, hot air drying system, pipes system, and electrical controlling system.
- The parts contacting with material such as chamber, basket, hollow main shaft, loading/unloading door, baskets auto open and close devise, pipes, etc are made of SS AISI316L.
- The inner chamber is round angle structure.
- Inner chamber is adopted with TIG welding.
- All the parts, inner surface of pipes inside chamber are polished, with roughness less than 0.5u.
- The outer surface of chamber is adopted ceramic for insulation. The outside surface is SS plate structure, the special design avoid the possibility of heat transfer.
- Unloading door has connections with sterile area
- Maximum temperature of equipment outer surface is 40°C
- Automatic loading
- Multiple ways for automatic unloading

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Series “LT” and “LST” equipped with complete controlling system, including touchable screen, PLC, recorder and printer. The controlling system is for controls of equipment and processing parameters, such as F0 value, pressure and water level, etc, it has simple HMI, and can define circulation process sequences, and can be operated by hand.

Program design conform to GAMP5, and FDA 21 CFR part 11 requirements.

In addition, one set of complete processing document is provided for assisting operator realizes proper operations.

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- Washing and rinsing  
Including: 5 phases of spray, overflow, drain, fluid bath and CIP auto washing
- Siliconizing
- Sterilizing
- Drying, cooling  
Time of washing, rinsing and each process can be adjusted according to the appropriate treating cycle of the products. The entire process is designed with maximum flexibility.