

水銀相對密度、濃度測試儀

Mercury Relative Density and Concentration Tester

KBD-600BM / 1200BM

適用於：水銀回收廠、生產廠的水銀相對密度、濃度測試。

Suitable for: Mercury recycle, Mercury relative density and concentration testing of manufacturer.

原理：根據GB/T21862、ASTM D1475、ISO2811-2、DIN53217-3之規定，採用阿基米德的浸漬Gamma球體積置換法。

Principle: GB/T21862, ASTM D1475, ISO2811-2, DIN53217-3 standards; Immersed Gamma volume displacement method of Archimedean principle.

| | |
|---------------------------|---|
| 型號 Model | KBD-600BM/1200BM |
| 測量範圍 Measuring range | 0.01g~600g/1200g |
| 密度精度 Density precision | 0.001g/cm ³ |
| 測試種類 Value showed | 直接讀取水銀的比重、濃度 Can show mercury density and concentration directly |



功能：

- 使用旋轉下降裝備可輕易將Gamma球浸入水銀液體中。
- 可任意調整Gamma球液面浸沫深度。可輕易取下Gamma球更換。
- **經醫學家證實水銀是有毒的液體，如吸入汞蒸氣，都可引起汞中毒。請為測試人員完善處理好排氣系統。**
- 任何重量的Gamma球砝碼皆可為其標準值，不受體積的限制，操作方便。
- 水銀比重資料很依賴溫度，可配合恒溫水槽得更精確的資料。
- 能補償因空氣浮力所造成的誤差。
- 可經由所測溶液的比重值直接顯示溶液的濃度值。

Functions:

- By using this machine, the density and concentration of mercury can be showed directly.
- By using spinning falling equipment, the Gamma ball can be easily dipped into the mercury liquid.
- The liquid immersed depth can be adjusted freely, and the Gamma ball can be replaced easily.
- **Mercury is poisonous liquid by medical scientists confirmed, if inhaled mercury vapor, it can cause mercury poisoning. Please supply good exhaust system for testers.**
- Any weight of Gamma ball can be as the standard value, without the limit of volume, easy to operate.
- Mercury density data depends on the temperature, the machine can cooperate with the thermostatic water tank to get more precise data.
- Can compensate the deviation caused by air buoyancy.
- Can directly show the solution concentration value according to the specific gravity of measuring solution.

水銀的密度和溫度的關係表：

The relational table of mercury density and temperature:

| T溫度(°C) | d密度(g/cm ³) | T溫度(°C) | d密度(g/cm ³) | T溫度(°C) | d密度(g/cm ³) |
|---------|-------------------------|---------|-------------------------|---------|-------------------------|
| 0 | 13.628 | 100 | 13.385 | 250 | 13.026 |
| 20 | 13.579 | 150 | 13.264 | 315.5 | 12.847 |
| 50 | 13.506 | 200 | 13.145 | | |