## EnviroLyzer-FE<sup>®</sup>

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**Interferences:** Metal ions like Lead [(Pb)<sup>2+</sup>] > 10 mg/l, Zinc [(Zn)<sup>2+</sup>] > 2 mg/l, Nickel [(Ni)<sup>2+</sup>] > 2 mg/l. Strong oxidizing agents, Cyanide [(CN)<sup>-</sup>], Nitrite [(NO<sub>2</sub>)<sup>-</sup>] and Phosphate [P] (poly phosphate more than orthophosphate), Chromium [Cr], Zinc [Zn] in concentrations exceeding 10 times that of Iron. Iron [Fe] and Copper [Cu] in excess of 5 mg/l, Nickel [Ni] in excess of 2 mg/l. Bismuth [Bi], Cadmium [Cd], Mercury [Hg], Molybdate and Silver [Ag] precipitate Phenanthroline.

Poly phosphate must be absent!

Large amounts of color and turbidity interferes.

Fats, Oil, Proteins, Surfactants and Tar.

#### Note:

-Initial boiling with acid solution converts Polyphosphate to Orthophosphate and removes Cyanide and Nitrite.

-Adding excess Hydroxylamine eliminates errors caused by excessive concentrations of strong oxidizing reagents.

-In the presence of interfering metal ions, use a large excess of Phenanthroline to replace that complexed by the interfering metals. Extraction methodology may be used when excessive concentrations of metals are present.

**Summary:** Iron is a compound with the symbol Fe. Iron exists from oxidation state -2 to +6, although +2 and +3 are the most common. It forms binary compounds with the halogens and the chalcogens. Elemental Iron is reactive; it oxidizes in air to give Iron oxides, also known as rust. The rusting of Iron and Iron alloys is undesirable, and has a major economic impact.

Ferric refers to iron-containing materials or compounds. In chemistry the term is reserved for iron with an oxidation number of +3, also denoted iron(III) or Fe<sup>3+</sup>. On the other hand, ferrous refers to iron with oxidation number of +2, denoted iron(II) or Fe<sup>2+</sup>. Iron(III) is usually the most stable form of iron in air, as illustrated by the pervasiveness of rust, an insoluble iron(III)-containing material.

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CONSUMABLES FOR 1 YEAR OPERATION: 一年内所需备品备件:

| Item                                    | Qty | Reference | Description                                    |
|---|-----|-----------|--|
| <b>EnviroLyzer® colorimetric 正常备品备件</b> |     |           |  |
| 1                                       | 1   | A0020290  | EPDM Duckbill for micropump 50µl (2pcs) 单向阀    |
| 2                                       | 1   | B0011500  | Tubing size 14 Norprene (15,2 m) (sample) 样品管  |
| 3                                       | 1   | C0010010  | Magnetic stirrer bar 13mmX3.0mm,COL 30mm 电磁搅拌子 |
| 4                                       | 1   | O0001710  | tubing pinch-valve 1/8"OD x 1/16" ID x 15m 管夹阀 |
| <b>稀释备品备件</b>                           |     |           |  |
| 5                                       | 1   | B0011600  | Tubing size 16 Norprene (15,2m) (dilution) 稀释管 |
| 6                                       | 1   | I0000300  | Valve/24000/6000/1000 (dispenser valve) 分配器阀门  |
| 7                                       | 1   | I0000700  | Syringe/6000 – 25 ml (range dependent) 分配器注射器  |

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