35 mm micro cast film line

•

10 times faster, 10 times cheaper R&D





This micro cast film line is based on Xplore's proprietary technology. In combination with our micro-compounder, it is now possible to obtain cast film with just 5 - 10 g of material in 5 - 10 minutes. In this way material development time is drastically reduced and considerably less R&D budget is needed.

A representative test film is now obtained in a few minutes with only a fraction of material used compared to conventional film development on laboratory equipment.

This film line is especially useful when only a small amount of material is available, or when expensive additives are used, like nano-tubes, active pharmaceutical ingredients, etc..

With our continuous feed kit, which consists of a mountable barrel hopper and continuous feed screws, it is now also possible to process larger samples into extremely long film via continuous extrusion.

Film characteristics

The micro cast film line is able to produce cast film and sheet with 5 - 10 g samples. The integrated air knife enables to process film up to a width of ca. 35 mm.

The cast film drawing speed is between 100 and 5000 mm/min and the cast roll of the first calender can be adjusted with increments of 1 mm/min. The winding roll of the second calender is torque controlled (its diameter will grow during the test run), and the winding torque can be adjusted with increments of 1 Nmm.

In combination with our micro-compounder, which is equipped with throughput control, Xplore enables you to produce a uniform film thickness without the use of a gear pump.

Latest results show that films as thin as 5 μm can be obtained without any holes.



Technical Specifications

- Cast film die: width 35 mm, slit thickness between 0.2, 0.4 or 0.6 mm
- I Calender with a speed controlled cast film roll and an adjustable guide roll (Speed controlled cast film roll can be adjusted between $100-5000\ \text{mm/min}$)
- I Calender with a torque controlled winding roll (Torque controlled winding roll can be adjusted in 230 steps with increments of I Nmm)
- Overall dimensions (I \times w \times h): 43 \times 25 \times 20 cm
- Weight: ca. 10 kg

Controls

- Control box with integrated operating display
- Supply voltage: 230 or 115 Volts

Optionally

- Custom made slit die thickness
- Surface finish of slit die can be custom manufactured
- Continuous feed kit to use with one of our micro compounders:
- Cooled hopper
- Continuous feed screws
- Additional calender(s) with a roll(s) with elastic surface finish for processing of slippery film
- Heated film die, width 35 mm (including control box)

Xplore Instruments BV

Brightlands Campus Urmonderbaan 22 6167 RD, Geleen The Netherlands

Tel: +31 46 4763388 Fax: +31 46 4763389

info@xplore-together.com www.xplore-together.com

Trade Register: 60040114