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Virtually all materials interact more or less strongly with humidity. The knowledge of moisture induced changes of the material properties is a key issue in decisions related to processing, packaging, storage and shelf life of a product.

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

The equilibrium water content of a sample plotted as a function of the relative humidity, at constant temperature, is called the water vapour sorption isotherm.

The sorption isotherm enables to draw conclusions on fundamental material properties and provides information on material handling with respect to moisture.

From the sorption kinetics curve of a material, moisture induced structural changes are determined that are often accompanied by sudden, temporary changes in moisture content.

微晶纤维素粉末的吸附等温线和吸附动力学
Sorption isotherm and sorption kinetics of microcrystalline cellulose powder

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Multi-sample instrument

The SPS sorption test systems enable the measurement of up to 23 samples simultaneously. This gives a very high analysis throughput compared with a single-sample instrument.

Sample material in small dishes is placed in a circular arrangement on a sample tray. By continuously turning and lowering of the tray, the dishes are placed onto the balance and weighed at fixed time intervals.

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Variable sample size

Sample dishes of different size enable the handling of samples up to 50 mm in diameter. A large dish gives a large sample surface area for sorption to take place, resulting in shorter measurement cycles.

High dynamic weight range

Sample weights from a few milligrams up to several grams are handled by the SPS system. In a test run, the full weighing range up to the maximum load is available for all samples independent from the sample weight.

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For materials with very low water uptake, high accuracy is achieved by using a larger volume distributed in a thin layer using a large dish.

All different kinds of samples from fine powders, granules, tablets and capsules, papers and foils, even cotton and fibers and also heavy and bulky samples like wood and construction materials are handled by the SPS instruments.

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SPS Versions and Technical Data

| | High versatility with the SPSx-1 μ series | |
|----------------------|--|---|
| SPS instrument | SPSx-1 μ Advance | SPSx-1 μ High Load |
| | High resolution and reproducibility for applications from pharmaceutical research up to material testing | High resolution and very high load range for applications in the food and chemical industry, including permeability tests |
| Number of samples | 11/23 | |
| Sample size, max | Ø 33 mm/Ø 50 mm | |
| Load range, min/max | < 10 mg - 22 g | < 30 mg - 22 g/220 g (dual range) |
| Resolution | 1 μ g | 1 μ g/10 μ g |
| Repeatability | $\pm 2 \mu$ g | $\pm 5 \mu$ g/ $\pm 20 \mu$ g |
| Temperature | 5 °C to 60 °C | |
| Accuracy | ± 0.1 K time/ ± 0.5 K uniformity | |
| Relative Humidity | 0 % RH to 98 % RH* | |
| Accuracy | ± 0.6 % RH (0 ... 98 % RH) at 23 ± 5 °C | |
| Hardware accessories | CCD video camera Raman spectroscopy | CCD video camera Raman spectroscopy Permeability Kit |
| Software option | CFR21 part 11 software package | |

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* Note: The full humidity range is only available at chamber temperatures slightly above room temperature.

| | Solutions for high demanding applications | |
|----------------------|---|--|
| SPS instrument | SPS23-100n | SPS11-10 μ |
| | Highest resolution and sensitivity for applications in pharmaceutical research and preformulation | High load range for applications in building physics, quality control and for permeability tests |
| Number of samples | 23 | 11 |
| Sample size, max | \varnothing 18 mm | \varnothing 50 mm |
| Load range, min/max | < 10 mg - 2000 mg | < 50 mg - 111 g/220 g (dual range) |
| Resolution | 0.1 μ g | 10 μ g/100 μ g |
| Repeatability | \pm 1.5 μ g | \pm 20 μ g/ \pm 100 μ g |
| Temperature | 5 °C to 40 °C | 5 °C to 60 °C |
| Accuracy | \pm 0.1 K time/ \pm 0.5 K uniformity | |
| Relative Humidity | 0 % RH to 98 % RH* | |
| Accuracy | \pm 0.6 % RH (0 ... 98 % RH) at 23 \pm 5 °C | |
| Hardware accessories | CCD video camera Raman spectroscopy | CCD video camera Raman spectroscopy Permeability Kit |
| Software option | CFR21 part 11 software package | |

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CCD video camera

In-situ documentation of moisture induced visual changes during a sorption measurement with an adapted high-resolution CCD camera.

Permeability Kit

Determination of the permeability of foils and films with the SPS moisture sorption instruments according the EN ISO 7783-1 for the determination of water-vapor transmission rate part 1: "Dish method for free films".

Raman spectroscopy

Sensor fusion combining water vapour sorption and Raman spectroscopy opens new and exciting perspectives for the analysis of solid materials.

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Sample preparation

Samples are simply filled into dishes that are easily placed on the sample tray. The sound construction and easy to access measurement chamber render sample preparation a comfortable procedure.

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
Step-by-step procedure

- A tabular input enables quick and easy setup of measuring cycles even for very complex water sorption studies
- Programmable sorption / desorption cycles
- Optional initial sample drying, or equilibration before starting the actual measurement
- Templates are available to specify standard measurement procedures
- User-friendly and intuitive interface with active user guidance prevent operating errors

SPS 4.2.0 New Project: 20130307_8CK302

| Date | Time | Elapsed Time | Cycle | Step | Sample | r.H. Set [%] | r.H. Read [%] | Temp. Set [°C] | Temp. Read [°C] | State |
|----------|------------|--------------|-------|------|--------|--------------|---------------|----------------|-----------------|---------|
| 3/7/2015 | 5:39:40 PM | 00:00:14 | 1 | 0 | | 0.0 | 2.9 | 25.0 | 25.0 | Running |

Setup | Climate | Samples | Measurement | Results

SPS

Logout

Dr. Petra Expert Lab Group Leader

Close

- Water
- Door
- Autosampler
- Peltier
- Column

On

Off

Climate

On

Off

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Temperature 25.0 25.0

Humidity 0.0 2.9

Set

Measurement

Project

New

Open

Data

- Climate
- Tare
- Weight

Measure

Start

Next Step

Stop

Pause

On

Next Sample

First Setting

Last Setting

Other

°C

% r.H.

Logging

Sample 1: Batch changed to CRM 0459
 Sample 7: Batch changed to CRM 0459
 Sample 1: Initial Moisture Content changed to 1.23 %
 Sample 7: Initial Moisture Content changed to 1.23 %
 Sat Mar 07 16:28:34 2015: Dr. Petra Expert Lab Group Leader
 Door opened
 Sat Mar 07 16:29:47 2015: Dr. Petra Expert Lab Group Leader
 Door closed
 Sat Mar 07 16:30:34 2015: Dr. Petra Expert Lab Group Leader
 Tare
 Sat Mar 07 16:36:29 2015: Dr. Petra Expert Lab Group Leader
 Confirmed sample changes:
 Sample 0: Tare changed to 1151.193 mg
 Sample 1: Tare changed to 1171.741 mg
 Sample 7: Tare changed to 1176.672 mg
 Sat Mar 07 16:37:29 2015: Dr. Petra Expert Lab Group Leader
 Climate set. Nominal 25.0°C 2.0 %r.H. Actual 25.0°C 20.0 %r.H.
 Sat Mar 07 16:37:45 2015: Dr. Petra Expert Lab Group Leader
 Door opened
 Sat Mar 07 16:38:38 2015: Dr. Petra Expert Lab Group Leader
 Door closed
 Sat Mar 07 17:07:30 2015:
 Logged out automatically
 Sat Mar 07 17:32:41 2015: Dr. Petra Expert Lab Group Leader
 Checked in
 Sat Mar 07 17:32:41 2015: Dr. Petra Expert Lab Group Leader
 Weight
 Sat Mar 07 17:38:39 2015: Dr. Petra Expert Lab Group Leader
 Confirmed sample changes:
 Sample 1: Wet Weight changed to 145.392 mg
 Sample 7: Wet Weight changed to 211.240 mg
 Sat Mar 07 17:39:26 2015: Dr. Petra Expert Lab Group Leader
 Measurement start

选配了CFR21的SPS软件
 SPS Software with CFR21 option

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Fully automated process

As soon as the “Start” button is pressed, experiments run unattended by automated measurement and data acquisition. The transparent glass cover enables the direct observation of the samples during the test.

Optimized gas flow

Permanent circulation of the atmosphere with an efficient gas flow along the sample surface and precise humidity regulation enable a fast equilibration of the samples to the respective climatic conditions, resulting in shorter measurement cycles.

Balance stability

An internal reference ensures an excellent drift compensation and baseline stability. Calibration of the balance is fully automated.

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Data format

Acquired data is saved in an electronically secured data format. Data import into an Excel® template offers all the statistical and visual capabilities of this standard office tool.


The graphics and tables in the Excel® template include:

- Dynamic change of the sample mass as a function of time and climatic conditions
- Equilibrium water uptake at different relative humidity levels


All graphics and calculations in Excel® are fully user configurable.

Data security

The standard SPS software is upgradable with an optional SPS CFR 21 package in order to meet the data security requirements for electronic records in accordance with CFR21 part 11.

| No. | Time [h] | Total Weight [mg] | Net Weight [mg] (Sample Tara) |
|-----|----------|-------------------|-------------------------------|
| 1 | 0,05 | 1115,448 | 136,824 |
| 2 | 0,19 | 1115,350 | 136,770 |
| 3 | 0,29 | 1115,330 | 136,754 |
| 4 | 0,41 | 1115,387 | 136,741 |
| 5 | 0,53 | 1115,399 | 136,745 |
| 6 | 0,64 | 1115,396 | 136,738 |
| 7 | 0,76 | 1115,387 | 136,728 |
| 8 | 0,88 | 1115,385 | 136,726 |
| 9 | 0,99 | 1115,377 | 136,718 |
| 10 | 1,11 | 1115,378 | 136,710 |

lowest weight 136,617
Data for Equilibrium Curve

| | RH Set [%] | Temp. Set [°C] | Temp. Read [°C] | Equilibrium (1) | Mol ratio | net weight at Equilibrium | |
|--|------------|----------------|-----------------|-----------------|-----------|---------------------------|---------|
| | 0,0 | 11,5 | 25,0 | 24,7 | 0 | -0,05847 | 136,617 |
| | 0,0 | 4,9 | 25,0 | 25,1 | 0 | -0,09792 | 136,62 |
| | 0,0 | 3,6 | 25,0 | 24,9 | 0 | -0,10958 | 136,61 |
| | 0,0 | 2,5 | 25,0 | 25,0 | 0 | -0,11907 | 136,62 |
| | 0,0 | 2,1 | 25,0 | 24,9 | 0 | -0,11611 | 136,63 |
| | 0,0 | 1,7 | 25,0 | 25,0 | 0 | -0,12125 | 136,62 |
| | 0,0 | 1,5 | 25,0 | 25,0 | 0 | -0,12857 | 136,62 |
| | 0,0 | 1,3 | 25,0 | 24,9 | 0 | -0,13005 | 136,62 |
| | 0,0 | 1,2 | 25,0 | 25,0 | 0 | -0,13688 | 136,63 |
| | 0,0 | 1,1 | 25,0 | 25,0 | 0 | -0,14171 | 136,64 |
| | 0,0 | 1,1 | 25,0 | 25,0 | 0 | -0,14785 | 136,64 |
| | 0,0 | 1,0 | 25,0 | 25,0 | 0 | -0,15399 | 136,63 |

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
The SPS instruments are designed for a multitude of industrial and scientific applications from different disciplines: pharmaceuticals, food science, cosmetics, plastics industry, chemistry, textile industry, biology, mineralogy, etc.

Dynamic vapour sorption is applied to measure

- Quantity and kinetics of water vapour sorption/desorption
- Equilibrium moisture content


of a variety of solid materials like powders, granules, pellets, tablets, films or packages.

食品原料的水蒸气吸附
Water vapour sorption of food materials

- Marshmallows
- Potato chips
- Basmati rice
- Linseed
- Set relative humidity

纺织材料的吸附动力学
Sorption kinetics of textile materials

- Sheep wool, raw
- Cotton, raw
- Cotton, white purified
- Set relative humidity

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
Specific applications are:


- Hydrate formation and moisture induced structural transformations (crystallisation)
- Thermodynamic parameters of materials such as heat of sorption and desorption
- Determination of porosity differences
- Diffusion and permeability studies (films and packages)
- Quality assurance: determination of subtle differences between different batches of a product

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
Crystallisation of amorphous lactose in skimmed milk powder

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Sorption isotherm of a hydrate forming amino acid

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Water uptake of three batches of cellulose powder at 70 % and 80 % RH, net weight: 100 mg

Perfect repeatability of results within the same batch
Precise determination of subtle differences between batches

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Visual documentation of the crystallisation of hard candy at high relative humidity

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High sample throughput

Measurement of up to 23 samples simultaneously at the same climatic conditions enables the direct comparison of the samples.

High dynamic weight range

Balance resolution of $0.1\mu\text{g}$ for samples up to 2 g (SPS23-100n) and $10\mu\text{g}$ for samples $> 100\text{g}$ (SPS11-10 μ) enables the analysis of small samples as well as the detection of small changes within large samples.

Base line correction

An internal reference weight ensures an excellent drift compensation over the whole weighing process.

Versatility

Optional equipment like video camera, Raman spectroscopy and permeability measurement of foils and packages is available for the SPS instruments.

Sustainable design

Even the first commercial available instruments from 1999 run with the latest SPS software version and are compatible with the newest hardware.