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Jan Philip Plog, Thermo Fisher Scientific, Material Characterization, Karlsruhe, Germany

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Sample	Portland Cement/g	Water/g	Fine Gravel/g
1	275	75	0
2	275	75	125
3	275	75	75

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1. See Wikipedia, Concrete, <http://en.wikipedia/wiki/Concrete> (as of Jan. 13, 2014, 16:50 CET)
2. ASTM C 143 - Standard Test Method for Slump of Hydraulic-Cement Concrete
3. EN 12350-2 - Testing fresh concrete. Slump test
4. N. Roussel, Correlation between Yield Stress and Slump: Comparison between Numerical Simulations and Concrete Rheometers Results, Materials and Structures, May 2006, Volume 39, Issue 4, pp 501-509
5. Dzuy NQ, Boger DV. 1985. Direct yield stress measurement with the vane method. J Rheol 29:335-47
6. Chong Hu, François De Larrard, Odd E. Gjovr, Rheological testing and modelling of fresh high performance concrete, Materials and Structures, January/February 1995, Volume 28, Issue 1, pp 1-7

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Francesca De Vito, Fabian Meyer and Fritz Soergel, Thermo Fisher Scientific, Material Characterization, Karlsruhe, Germany

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Raspberry jam	τ_0 in Pa	t in s
With seeds	590	100
Sieved	967	104
Sieved (measurement 2)	977	104

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Peanut butter	τ_0 in Pa	t in s
Measurement 1 (blue)	1818	52.1
Measurement 2 (red)	1785	52.1
Measurement 3 (green)	1829	53.2

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1. Sun A., Gunasekaran S. "Yield stress in foods-measurements and applications" *Intl. J. Food Prop.* 12(2009) 70-101
2. Steffe J. F. "Rheological methods in food process engineering". 2nd Ed.; Freeman Press, East Lansing(1996)
3. Rao M. A. "Rheology of fluid and semisolid foods – Principles and applications". 2nd Ed., Springer, New York (2007)
4. Oldörp K. "Applied food rheology – Using fast speed control and axial measurements". Thermo Scientific Application Note V238 (2009)
5. Schramm G. "A practical approach to rheology and rheometry". Thermo Scientific, Karlsruhe (2004)
6. DIN Report 143 "Modern rheological test methods – Part 1: Determination of the yield point – Fundamentals and comparative testing methods". (2005)
7. De Vito F., Meyer F., Soergel F. "Investigation of the effect of fatcontent on the yield stress of mayonnaise measured with HAAKEViscotester iQ and vane rotor". Thermo Scientific Application NoteV271 (2014)
8. Javanmard M., Endan J. "A survey on rheological properties of fruit jams". *Int. J. Chem. Eng. Appl.* 1 (2010) 31-37
9. Fischer A. "Qualitätskontrolle mit Hilfe der HAAKE RheoWin Software". Thermo Scientific Application Note V223 (2007)

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热线电话：
800 810 5118
400 650 5118
sales.cad@thermofisher.com
www.thermofisher.com