Monolith® NT.Automated 分子相互作用分析

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
1. Amaning K, Rak A. et al., (2013) The use of virtual screening and differential scanning fluorimetry for the rapid identification of fragments active against MEK1. *Bioorganic & medicinal chemistry letters* 23:3620-3626
2. Breitsprecher D, et al., (2014) Case Study on Automated Screening Project of a Fragment Library against MEK1 kinase, Application Note NT021

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"One mission of the Structural Research group within the Lead Identification and Optimization department is to quantitatively validate and characterize interactions of small molecules as well as new biological entities with protein targets. We employ various traditional biophysical methods such as SPR, Thermal Shift, ITC, NMR and X-ray crystallography. Most recently we included Microscale Thermophoresis (MST) in our standard project support workflow and extended its application from affinity determinations to fragment screening approaches.

MST is a versatile and valuable tool which we quickly adapted in our repertoire of methods. The impressive advantages of MST, namely the low sample consumption, the broad application range, and swift assay development make it a unique biophysical method. The measurement in free solution without the need of surface coupling saves time and avoids a potential source for false positive or negative results. Our Structural Research group now also added the label-free version of the Monolith to our MST instrumentation portfolio, which gives us the possibility to choose to measure with high selectivity and sensitivity (NT.115). The label-free version (NT.LabelFree) allows us to measure without any additional sample modification depending on the need of the particular assay. In some cases, LabelFree MST allowed us to perform assays with otherwise "very ill" behaved proteins which were not amenable to any other biophysical technique. Generally, we find very good consistency between quantitative MST measurements and results stemming from other biophysical methods."

Dr. Nicolas Basse, Department of Structural Biology, UCB-Celltech, UK
英国 UCB 优时比制药

"As part of our drug discovery projects we use Microscale Thermophoresis (MST) as an orthogonal method to measure the binding affinity of compounds to their protein target and apply this fragment hit ID through to lead optimisation. MST complements our biophysical platform and has correlated well with other more established technologies. Because it uses small amounts of protein, MST has proved to be particularly useful to look at molecular interaction involving proteins that are difficult to express or purify. MST requires a relatively short time to setup new assays and is a powerful technique for buffer optimisation. Using the NT.115 MST instrument we have successfully measured small molecule-protein and protein-protein interactions in complex media. Finally, LabelFree MST is one of the few true label-free/immobilisation-free instruments capable of measuring molecular interactions. We have appreciated the professionalism and support from NanoTemper Technologies."

Dr. Alexey Rak, Structural Biology & Biophysics, Sanofi R&D, France
法国赛诺菲集团

"We routinely assess interaction affinity for both small molecule and biologics projects, with NanoTemper Technologies' Microscale Thermophoresis being the most recent addition to the pool of instruments we use to carry out these measurements. It has proved a valuable tool for characterising small molecule-protein and protein-protein interactions, as well as for the study of protein aggregation concentration determination. There is very good agreement with other technologies such as Surface Plasmon Resonance (SPR) and Isothermal Titration Calorimetry (ITC), and we are particularly appreciative of this new technology because of the extremely low protein consumption and relatively short time required for the assay setup.

Dr. Timothy Sharpe, Head of the Biophysics Facility, Biozentrum, University of Basel, Switzerland
瑞士巴塞尔大学

"We have used our NanoTemper NT.115 MicroScale Thermophoresis (MST) instrument extensively in the last one and a half years to study many different types of interaction. Where we have been able to make comparisons, results from MST agree well with those from other established techniques (ITC, fluorescence intensity and anisotropy, SPR). However, we particularly appreciate two distinguishing aspects of MST. Firstly, MST can measure interpretable signal changes in circumstances where many other techniques struggle: i.e. small unlabelled molecules binding to large labelled molecules, where binding is tight (nM Kd) and one or both partners aggregate at micromolar concentrations. Secondly, one can often gain extra information about the interaction partners from MST titrations, particularly for systems that have changes in their conformational or oligomeric state upon binding. In several cases, multi-phasic titrations have inspired experiments with other techniques to characterize behaviours that can't be explained by the simplest binding models."

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1. Ascher, D. B., Wielens, J., Nero, T. L., Doughty, L., Morton, C. J., and Parker, M. W. (2014) Potent hepatitis C inhibitors bind directly to NS5A and reduce its affinity for RNA. *Sci. Rep.* 4
2. Jerabek-Willemsen, M., André, T., Wanner, R., Roth, H. M., Duhr, S., Baaske, P., and Breitsprecher, D. (2014) MicroScale Thermophoresis: Interaction analysis and beyond. *Journal of Molecular Structure*
3. Liebner, R., Mathaes, R., Meyer, M., Hey, T., Winter, G., and Besheer, A. (2014) Protein HESylation for half-life extension: Synthesis, characterization and pharmacokinetics of HESylated anakinra. *European Journal of Pharmaceutics and Biopharmaceutics*
4. Barandun, L. J., Immekus, F., Kohler, P. C., Ritschel, T., Heine, A., Orlando, P., Klebe, G., and Diederich, F. (2013) High-affinity inhibitors of *Zymomonas mobilis* tRNA-guanine transglycosylase through convergent optimization. *Acta Crystallographica Section D* 69, 1798-1807
5. Xiong, X., Coombs, P. J., Martin, S. R., Liu, J., Xiao, H., McCauley, J. W., Locher, K., Walker, P. A., Collins, P. J., Kawaoka, Y., Skehel, J. J., and Gamblin, S. J. (2013) Receptor binding by a ferret-transmissible H5 avian influenza virus. *Nature* 497, 392-396
6. Seidel, S. A., Wienken, C. J., Geissler, S., Jerabek-Willemsen, M., Duhr, S., Reiter, A., Trauner, D., Braun, D., and Baaske, P. (2012) Label-free microscale thermophoresis discriminates sites and affinity of protein-ligand binding. *Angew Chem Int Ed Engl* 51, 10656-10659
7. Ramakrishnan, M., Alves De Melo, F., Kinsey, B. M., Ladbury, J. E., Kosten, T. R., and Orson, F. M. (2012) Probing cocaine-antibody interactions in buffer and human serum. *PloS one* 7, e40518
8. Lippok, S., Seidel, S. A., Duhr, S., Uhland, K., Holthoff, H. P., Jenne, D., and Braun, D. (2012) Direct detection of antibody concentration and affinity in human serum using microscale thermophoresis. *Analytical chemistry* 84, 3523-3530
9. McLaughlin, S. H. (2011) Binding of the geldanamycin derivative 17-DMAG to Hsp90 measured with fluorescence label and label-free. *Application Note NT001*

Europe / International

NanoTemper Technologies GmbH
Floessergasse 4
81369 Munich
Germany

Phone +49 (0) 89 4522895 0
Fax +49 (0) 89 4522895 60

info@nanotemper-technologies.com

USA / Canada

NanoTemper Technologies, Inc.
400 Oyster Point Blvd. Suite 336
South San Francisco
CA 94080, USA

Phone + 1 650 763 1658
Fax + 1 650 350 4390

info@nanotemper-technologies.com

Great China Area

NanoTemper Technologies Ltd.
Room0305, Air China Plaza,
No.36 Xiaoyun Road, Chaoyang District
Beijing, 100027

Phone +86 (10) - 84462100
+86 (10) - 84462300

info@nanotemper-technologies.com

Latin America

NanoTemper Technologies GmbH
Center Alameda Santos
Alameda Santos, 200
São Paulo -SP 01418-000, Brazil

Phone: +56 (2) 2581-4664
or +55 (11) 3587-1408

info@nanotemper-technologies.com

MicroScale Thermophoresis™ is a trademark.
NanoTemper® and Monolith® are registered trademarks.
NanoTemper® and Monolith® are registered in the U.S. Patent and Trademark Office.