

CO₂ adsorption into zeolite 13X with PCTPro

Introduction:

The rising level of CO₂ in atmosphere has been linked to global warming, threatening life as we know it. To mitigate the global warming, R&D is being directed towards understanding the relevant phenomena and foster innovation in the field of CO₂ capture and sequestration (CCS). Due to their well-controlled pore structure and size, zeolites have been primary candidates in the gas separation (e.g. CO₂ capture) in industry. Knowledge about the CO₂ sorption properties of zeolites (adsorption capacity, pressure regimes and kinetics) is essential to the design of advanced materials capable of capturing CO₂ in industrial settings. Among zeolites, 13X is known for its relatively high CO₂ capacity. This application note highlights precision measurements of the absorptive properties of a zeolite 13X over a wide range of temperatures.




Figure 1. Structure of a type X zeolite

Experimental

CO₂ adsorption into zeolite 13X was measured at various temperatures using a PCTPro E Sievert's apparatus which was developed for the study of sorption of a variety of gases from vacuum up to 200 bar and from liquid He to 500 °C. Temperature correction was done by measuring the apparent free gas volume at temperature.

Results and discussion

The PCT isotherms for CO₂ adsorption into zeolite 13X are shown in Figure 2. The zeolite capacity decreases with temperature reflecting the physisorption nature of the adsorption isotherms. The data are in good agreement with literature. For example, the CO₂ capacity at 30 °C 20 bar is 5.7 moles/kg (5.0-6.4 moles/kg in the literature). The PCTPro E&E is well-suited for the detailed characterization of materials used in CCS (adsorption of CO₂ onto different solid sorbents). The ease of use and the temperature and pressure range are ideal for this type of materials application.

Literature

1. D Bonenfant, M Kharoune, P Niquette et al. Sci. Technol. Adv. Mater. 9 (2008) 013007
2. R Siriwardane, M Shen, E Fisher, et al. NETL report, www.netl.doe.gov

www.setaram.com – sales@setaram.com




Figure 2. PCT- isotherms for zeolite 13X at 30, 50, 80, 120 and 180 °C



PCTPro-E&E Sievert's apparatus