

Automated Sample Preparation of Whole Blood for Therapeutic Drug Monitoring and Diagnostics by LC-MS using a Commercial Autosampler

Christian Berchtold¹, Irene Wegner¹, Timm Hettich¹, Reto Bolliger², Guenter Boehm², Goetz Schlotterbeck¹

¹ University of Applied Sciences and Arts Northwestern Switzerland FHNW, Institute for Chemistry and Bioanalytics, Gründenstrasse 40, 4132 Muttenz

² CTC Analytics AG; Industriestrasse 20, CH-4222 Zwingen


Introduction

Automated sample preparation reduces the costs per sample and avoids sample handling errors. This is especially important in therapeutic drug monitoring or diagnostics based on blood samples. The use of robots is well established in these fields since a long time. Usually expensive and highly specialized pipetting robots are used for a high number of samples. However, most of these systems are not designed with a direct interface for LC-MS applications.

In this project the parameters necessary to automatically prepare whole blood samples for on-line LC-MS applications in the field of diagnostics and TDM have been investigated. A strategy and the most important parameters are shown for the optimization of a PAL RTC autosampler for the preparation of whole blood.

Materials and Methods


An integrated system consisting of a PAL RTC autosampler equipped with a 1 mL syringe, a LC-MS injection tool, a vortexer and a centrifuge was used (see below).



Tested laboratory blood was supplied by blood donation and applied into 1.5 ml standard LC-vials with magnetic caps for vial transport. Methanol, acetonitrile, 0.1 mol/L ZnSO₄ solutions in methanol and a mix of isopropanol/water have been tested to achieve protein precipitation. Blood volumes in the range of 50 to 200 µL were investigated at centrifugation times of 1 to 20 minutes (3000 g).

Sample Preparation Procedure



Since the PAL operates with standard vials its necessary to modify the sample preparation steps accordingly




Results - Precipitation and Extraction

Precipitating and extracting depends on the solvent used and is always dependent on the metabolites analysed. The target is the best compromise between efficiency of extraction, pellet size and clear supernatant.


Precipitated and centrifuged



Pellet height depending on the duration of centrifugation and solvent used for 50 µL blood




Pellet height for 0.1 mol/L ZnSO₄ solution in MeOH at 10 min centrifugation. Blue is the high point and red the low one of the pellet.



Total ion chromatogram (TIC) of the solvent after precipitation. Only minor differences are observed.

Needle penetration parameters



To avoid clogging of the needle should not penetrate below 8 mm up to 100 µL. A total sample volume of 200 µL blood (per 1.5 mL vial) causes a pellet which is of acceptable height (10 mm maximum).

Results - Therapeutic Drug Monitoring (TDM)

Diclophenac


Diclophenac is a very typical pain killer. Although it is not very important for general TDM it is a very typical molecule for small molecule with active properties.

Method: Eluent A: Water 0.1% Formic acid / Eluent B: Acetonitril 0.1% Formic acid // Flow: 0.5 ml/min // Runtime: 5 min

Column: Waters Xbridge C18 3x50 mm 3.5 µm
MS Parameter: Mode: negative Gradient:

Fragmentor	Parent ion	Collision	Product ion
50 V	294	50	214

Time	%A
0	70
0.2	70
2	0
3.9	0
4	70



Example 300 ng/ml Diclophenac in whole blood (spiked)

Blood Volume: 50 µL

Injection volume: 2 µL

Precipitation solvent: Acetonitril

Vortex time: 90s Centrifugation: 10 min

Carbamazepine


Carbamazepine is a representative of an anti-epileptic drug and very important to monitor in blood to set a proper and individual dosage. [1]

Method: Eluent A: Water 0.1% Formic acid / Eluent B: Acetonitril 0.1% Formic acid // Flow: 0.7 ml/min // Runtime: 5 min

Column: Waters Xbridge C18 3x50 mm 3.5 µm
MS Parameter: Mode: positive Gradient:

Fragmentor	Parent ion	Collision	Product ion
135 V	237	18	194

Time	%A
0	80
3	0
4	0
5	0
4.1	80



Example: 1 µg/ml carbamazepine in whole blood


Blood Volume: 50 µL / Injection volume: 2 µL

Precipitation solvent: Acetonitril

Vortex time: 90s Centrifugation: 10 min

Immunosuppressiva

Immunosuppressive are measured by the use of an online SPE method, which uses a short column for pre-cleaning and a second column for separation. The method is based on Seger et al. [2]



Example: Cyclosporine: 1100 µg/L, tacrolimus 34.8 µg/L, everolimus 32.7 µg/L, sirolimus, 39.2 µg/L of whole blood reference standard by chromsystems (No 0081)

Results - Diagnostics

Metabolites, which are of diagnostic value, are shown as a typical example. The method is comparable to those used in newborn screening based on dried blood spots [3]

Eluent A: Water 0.1% Formic acid, 2 mM ammonium fluoride

Eluent B: Methanol 0.1% Formic acid, 2 mM ammonium fluoride Isocratic 35% A

Flow: 0.6 ml/min Runtime: 5 min

Column: Waters Xbridge C18 3x50 mm 3.5 µm

Blood Volume: 50 µL

Injection volume: 2 µL



Crashing solvent: Acetonitril

Vortex time: 90 s

Centrifugation time: 3 min

MS Parameter: Mode: Positive

	Fragmentor	Parent ion	Collision	Product ion
Cortisol	120 V	263.2	25	121.1
C4-carnitine	120 V	232.1	20	85.1
C3-carnitine	100 V	218.1	20	85.1
C2-carnitine	100 V	204.1	20	85.1
Arginin	100 V	175.1	20	70
Carnitin	100 V	162.1	25	85.1
Leucin	80 V	132.1	10	86
Valin	80 V	118.1	10	72
Prolin	100 V	116.1	0	116.1



Example: Analysed donated laboratory blood

References

- [1] Linder et al. Bioanalysis (2015) 4 (16), 2013-2039
- [2] Seger et al. Nature Protocols (2009) 4 (4), 526-534
- [3] Fingerhut et al. Rapid Communications in Mass Spectrometry (2014) 28 (8), 965-973

