

NORHOF N2 Cold Gas Supply System #855 Series

- Controllable N2 Cold Gas flow up to 60 liter / minute;
- Several N2 Cold Gas control possibilities:
 - Pressure control;
 - PID temperature control with one temperature sensor;
 - Software monitoring

Norhof manufactures N2 cold gas supply systems. Nitrogen (LN2) is used as the cooling medium and is taken from a storage vessel (Dewar) with low pressure (max. 300 mBar) and delivered (pumped) through a line to the application in a gaseous cold flow.



Working principle

Norhof 855 series pump, mounted on a 50 Liter Dewar

The pressure above the liquid level inside the Dewar is built by heating a small amount of liquid in the bottom of the Dewar. With only up to 300 mBar of overpressure, the cold gas will gentle flow out of the outlet pipe. In these systems a pressure less storage Dewar is used as a reservoir for LN2. On the Dewar our unique cold gas pump is mounted. Inside the housing of the pump a microprocessor is used to control the various actions, depending upon the mode selected. The selection of a particular mode is made by a 16-position switch, also located inside the pump housing. By controlling the pressure at time of cold gas transport, the flow is controlled and thus any temperature between ambient and -160°C can be realized.


Important: cold N2 gas supply may have its application advantages, however, important to know it has much lower heat removal capacity compared with liquid Nitrogen

Any #855 series model will work without modifications as:

- autonomous stand-alone unit;
- direct remote controlled unit;
- (controlled by your existing PID controller, PLC, computer with A/D conversion card, etc.);
- will operate on 24 V, DC or with our power supply (115V or 230 V);
- a system with almost no installation time required.
- no need for a pressurized supply of LN2;
- no need for a cryogenic solenoid valve ;
- no need for additional control instruments.

Software

To display sensor temperature, vessel pressure, status of LED's on the pump etc. our Norhof Monitoring software is included with #855 series pump. This software works under Windows '98 - 2000 - ME -NT - Vista - W7 – Windows 10





LN2 Microdosing Systems

855 series Technical Specifications

Static evaporation rate	< 0,5 liters per day		
Flow rate	adjustable from 0.01mBar to 270 mBar (by potmeter on pump)		
	== up to 60 liter / minute		
	(with a fill line of 2 meter and 6.6 mm ID)		
Maximum working pressure	< 300 mBar		
Reaction time	+/- 1 minutes for cooling down the fill line (with 2 meters fill line)		
Power connection	115V / 230V AC with supplied power supply or 24 Volt DC		
Power consumption	average 125 Watts, during pumping 250 watts		
Storage container volume	50 Liter	100 Liter	
Outside dimensions (diameter)	500 mm	500 mm	
Height dimensions	875 mm	1235 mm	
Weight (empty, full)	17 / 57,5 kg	32 / 113 kg	
Standard fill line	9.6 mm OD, 6.6 mm ID PTFE tube, with 32mm foam insulation		
System includes	Dewar, pump, fill line 2.00 m, power supply, cables, 2 level sensors,		
	PC software.		
Working modes	#855, working mode 1 to 7 for control		
External control	#855, 5 volt signals for ON, OFF and RS232 signals for ON, OFF		
PC software	Monitor software, to monitor pump behavior, and for some working		
	modes to adjust some parameters.		
Alarms/warning	Dewar empty, Dewar 5 liters LN2 left, broken sensor(s), frozen		
acoustical/ visual / mechanical	alarm, mechanical overpressure protection valve.		
Options	Transport trolley 5 wheels (10 cm height)		
	Stand for pump (when Dewar is refilled)		
	Custom built adaptor to fixate sensor(s) on application		

Summary of #855 series

Model:	#855
autonomous operation	V
remote control by TTL signal	V
remote control by +24V signal	V
display system status on PC screen	٧

855 series advantages:

- there is no LN2 valve required; that implies no unnecessary heat input
- there is no additional control unit required; ۵ which adds to a clean and elegant setup
- the pump is software driven and many control modes are already built-in; ۵ subzero temperature control, flow control by 0-5V input, control by RS232 line, etc.
- temperature sensor(s) are plugged directly into the pump housing; ۵ and not into a separate control box
- the system can deliver N2 cold gas with a flow optimized for the application; without noise, vibration, excessive waste, etc.
- the variable flow feature makes subzero temperature control extremely easy; ۵ see the website under SOFTWARE > sample drivers > freeze curve
- the system is prepared to be connected to a PC; perfect for monitoring and data logging or remote control
- P.E.D. 99/36/EC (Pressure European Directive) for pressurized vessels does not apply for ٨ this system;

The maximum possible pressure is lower than 300mBar. Therefore this system is allowed to be used inside the lab, near your working place, without danger.



奥然科技 上海

地址:上海市徐汇区漕溪路250号银海大厦A1012室 电话: 021-6451 6293 6451 6297 6451 6299

奥然科技 北京

地址:北京市海淀区北三环西路甲30号双天大厦423室 电话:010-6840 0176 6840 0177

奥然科技 广州

地址:广州市天河区五山路141号尚德大厦A栋1707室 电话:020-8551 6682 8551 6679

奥然科技 南京

地址:南京市中央路417号先锋广场1036室 电话:025-5850 1866 5850 1966