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
Figure 2: Breakthrough curves determined for an 80%/20% N₂/O₂ mixture at 10 bar and 25°C for 13X zeolite. The concentration is the mass spectrometer signal, in arbitrary units, for the peaks in the mass spectra at m/e = 28 (N₂) and 32 (O₂), as determined from the peak heights.

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Metal-organic framework with optimally selective xenon adsorption and separation

Debasis Banerjee¹, Cory M. Simon², Anna M. Plonka³, Radha K. Motkuri⁴, Jian Liu⁴, Xianyin Chen⁵,

Berend Smit^{2,6}, John B. Parise^{3,5,7}, Maciej Haranczyk^{8,9} & Praveen K. Thallapally¹



Figure 3 | Single column breakthrough experiments using SBMOF-1 at room temperature and 1 atm. Column is initially purged with He. (a) Inlet is a dry gas mixture with 400 p.p.m. Xe and 40 p.p.m. Kr balanced with air. (b) Inlet is the same gas mixture as in (a) with 42% relative humidity. Note that the Xe breakthrough time is only marginally decreased in the presence of water.

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Selective CO₂ capture from flue gas using a metalorganic framework—a fixed bed study
Jian Liu†, Jian Tian†, Praveen K.

Thallapally†,*, and B. Peter McGrail†

†Energy and Environment Directorate, Pacific Northwest National Laboratory, Richland, WA 99352, USA.

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
NEW PRODUCT

ABR AUTOMATED BREAKTHROUGH ANALYZER

Hidden Isochema announces
the launch of the all-new

ABR Automated Breakthrough Analyzer.

KEY FEATURES

- ▶ Fully automated breakthrough analyzer
- ▶ Optimised for research-scale sample sizes (bed sizes from 2 cc)
- ▶ Ultralow dead volume automated switching (purge gas /process mixture)
- ▶ Options for gas-gas, gas-vapor, vapor-vapor configurations
- ▶ Optional high pressure (50 bar) configuration
- ▶ Fully integrated close-coupled mass spectrometer

For more information
please contact our Sales
and Application team via
info@hiddenisochema.com

The **ABR** is a dedicated breakthrough analyzer, fully automated and supplied with an integrated close-coupled mass spectrometer.

The ABR is available in a range of configurations to suit research-scale samples, with bed volumes from 2 cc to 20 cc. Up to 6 gas inlets are available as well as a dedicated purge stream. Flow rates are selected to suit the specific applications, and the ABR includes an ultra-low dead volume switching valve. Options include an upgrade for operation at pressures to 50 bar, and an integrated vapor generator module for gas-vapor operation.

The ABR is designed to meet the needs of researchers wishing to characterise the gas separation performance of novel materials such as metal-organic frameworks (MOFs), zeolitic imidazolate frameworks (ZIFs) and covalent organic frameworks (COFs) without the time or expense of synthesising larger quantities of material. The breakthrough data obtained is complementary to the adsorption-desorption isotherms measured with our iGA, iGAsorp, XEMIS and IMI sorption analyzers.

Hidden Isochema ABR System: Academic References

Listed below are peer-reviewed publications featuring data measured using Hidden Isochema ABR breakthrough analyzers. All data uses Hidden Isochema ABR breakthrough analyzers with integrated Hidden DSMS dynamic sampling mass spectrometers. In all cases, data was measured by end users at their laboratory.

- 1. Porous organic cages for sulfur hexafluoride separation**
T. Hasell, M. Miklitz, A. Stephenson, M. A. Little, S. Y. Chong, R. Clowes, L. Chen, D. Tribello, K. E. Jelfs, and A. I. Cooper
Journal of American Chemical Society (2016)
DOI: 10.1021/jacs.5b11797
- 2. Metal-organic framework with optimally selective xenon adsorption and separation**
D. Banerjee, C. M. Simon, A. M. Plonka, R. K. Motkuri, J. Liu, X. Chen, B. Smit, J. B. Parise, M. Haranczyk and P. K. Thallapally
Nature Communications (2016)
DOI: 10.1038/ncomms11831
- 3. Metal-organic frameworks for removal of Xe and Kr from nuclear fuel reprocessing plants**
J. Liu, P. K. Thallapally and D. Strachan
Langmuir, 2012, 28 (31), pp 11584–11589
DOI: 10.1021/la301870n
- 4. Selective CO₂ Capture from Flue Gas Using Metal–Organic Frameworks—A Fixed Bed Study**
J. Liu, P. K. Thallapally and B. P. McGrail
J. Phys. Chem. C, 2012, 116 (17), pp 9575–9581
DOI: 10.1021/jp300961j
- 5. Supramolecular binding and separation of hydrocarbons within a functionalized porous metal–organic framework**
S. Yang, A. J. Ramirez-Cuesta, R. Newby, V. Garcia-Sakai, P. Manuel, S. K. Callear, S. I. Campbell, C. C. Tang and M. Schröder
Nature Chemistry 7, 121–129 (2015)
DOI: 10.1038/nchem.2114
- 6. Separation of rare gases and chiral molecules by selective binding in porous organic cages**
L. Chen, P. S. Reiss, S. Y. Chong, D. Holden, K. E. Jelfs, T. Hasell, M. A. Little, A. Kewley, M. E. Briggs, A. Stephenson, K. M. Thomas, J. A. Armstrong, J. Bell, J. Busto, R. Noel, J. Liu, D. M. Strachan, P. K. Thallapally and A. I. Cooper
Nature Materials 13, 954–960 (2014)
DOI: 10.1038/nmat4035