V-COMPACT 、 V-PAD 旋转粘度计

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



The screenshot shows the Fungilab Data Boss Help window. The left pane contains a table of contents with the following structure:

- Contenido | Índice
 - Introduction
 - About
 - System requirements
 - Registering Fungilab Data Boss
 - Contacting Fungilab
 - Getting started
 - Software Installation
 - Hardware installation
 - Starting Fungilab Data Boss
 - Using Fungilab Data Boss
 - Login Window
 - Main window
 - Plot Window
 - Samples Window
 - Comments Window
 - New Experiment Window
 - Load Queue Window
 - Open Experiment Window
 - Spindles Window
 - Math models Window
 - Password Window
 - Users Window
 - Edit User Window
 - Logs Window
 - Troubleshooting Fungilab Data Boss
 - Login failure
 - Connection error
 - Database error

About Fungilab Data Boss

Fungilab Data Boss is a software for Windows designed to gather data from Fungilab viscosimeters. Its features include:

- Customizable graphical user interface
- Numerical and graphical representation of gathered data
- Background execution of data gathers (real-time plotting of data)
- Fast design of experiments
- Export data and plot to several formats
- Mathematical analysis of collected data
- Users' privacy
- Log reporting
- Viscosimeter plug and play support

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The FDB software developed by Fungilab has been designed to work with the Premium series viscometer. It offers the possibility to program the viscometer for step, ramp and multi-step curves, and it is an important tool to study the viscosity behaviour of different liquids.

Control software

Establishing viscosity programs, documenting the procedure and the results in real time are some of the options offered by the Fungilab Data Boss software developed by FUNGILAB.

The software has been designed to be very intuitive. The option menu provides a clear view of the program. The help menu provides charts and descriptive texts that offer a clear explanation of the working options.

Features

- > Complete viscometer control.
- > Easy to use. All programs eliminate user errors when programming the instrument to collect data.
- > Provides instantaneous viscosity flow curves when performing new experiments, with definable parameters.
- > Clear view of program options using flanges.
- > Definable graphics and zoom function.
- > Different types of experiments can be programmed: simple curves, ramps and multi-step.
- > All experiments are recorded in different databases to be able to consult them anytime.
- > Experiment documentation with name, number and additional data.
- > In order to compare different flow curves, experiments can be plotted simultaneously.
- > Over 12 different charts can be obtained

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