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Newport Scientific, Inc. *formerly Aminco*

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
Newport Scientific, Inc. formerly Aminco

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
## SCF EXTRACTION DRYING AND CLEANING SYSTEM --- NP Series

### ---FEATURES

- Designed for easy cleaning and filling of vessels
- Flow rate controlled by back pressure regulator
- Digital temperature control with
  - Extraction vessel
  - Separation vessel
  - Back pressure regulator
- Digital pressure readings with output
  - Extraction vessel
  - Separation vessel
- Separation at atmospheric and controlled pressure @ 10,000 psi (back pressure regulator)
- Proven components
- Can Pump gas or liquid CO<sub>2</sub> without any change required
- Valve controlled flow rate
- Digital totalizer
- Flow meter with output (higher flow rate available)
- Non-contaminating diaphragm compressor
- Complete flexibility in a single system permits modifications readily, when necessary
- Recorder output available

### OVERVIEW

- Supercritical Fluid (SCF) extraction is receiving considerable attention as a means of separation in the chemical, food, pharmaceutical, and energy industries. The process is the latest technique used to extract oils from natural products, organic pollutants from wastewater, aromatic isomers from mixtures, low-molecular weight materials from polymers and light components from coal.
- Several firms are employing the technique in various commercial applications, such as coffee decaffeination or botanical extraction. Newport Scientific offers the Supercritical Fluid Extraction System is a complete turnkey for extractions up to 10,000 psi and separations up to 10,000 psi.

### SYSTEM OPERATION

Gas or fluid from your commercial gas bottle passes through a filter to the compressor. The compressor stroke rate is 58 RPM. The extraction pressure and flow rate is controlled by a back pressure regulator as is the separation pressure. The extractions and separation temperatures are controlled and indicated independently. A second atmospheric separation may be performed in a glass vessel or cold trap before the gas stream passes through the flow rate indicator and flow totalizer and is vented to the atmosphere.

### SPECIFICATIONS

#### Standard Configuration

**Temperature and Pressure Display:** Independent LED displays.

**Temperature Range:** Ambient to 250°C.

**Temperature Precision:** +/- 0.5°C.

**Operating Pressure:** 10,000 psi upper pressure limit. Front keypad control, with LED display. "Constant pressure" mode of operation.

**Flow Rates:** 15 L/Hour liquid CO<sub>2</sub> (+/- 2% accuracy), Larger rate available.


**Over Pressure Safeguard:** High / Low pressure alarms and rupture disc protection.

**Electrical Driven Diaphragm Compressor:** Double ended diaphragm compressor. Dual aluminum heads, furon seals and sapphire pistons, integrated thermoelectric cooling, cam-driven pump mechanism with single stepper motor drive, dual ball and seat check valves (ruby ball, sapphire seat), stainless steel fluid path, prime-purge valve, and pressure transducer. The pump's constant pressure mode features a selectable pressure set point. Flow rate will auto-adjust to maintain pressure.

**Restrictor Valve:** Heated up to 250°C; resistant to blockage.

#### Data Acquisition System (DAQ)

The compact OM-Daq-USB-2401 is used for portable data acquisition applications, with accessory probes (shown L to R), KMQSS-125U-12, TMQSS-125U-12, JMQSS-125U-12.

#### Extraction Vessel Sizes:

- 50, 100, 500 mL
- 1, 2, 5, 10, 20, 50 L

Vessels come with 5 micron fritts and CO<sub>2</sub> flow direction controller, interchangeable or parallel operation.

#### Separation Vessel Sizes:

- 50, 100, 300, 500 mL
- 1, 2, 5, 10, 20, 50 L

Vessels come with 5 micron fritts and are interchangeable or parallel operation.

**Dimensions:** depending on the configuration

**Weight :** depending on the configuration

#### Configuration Options

**Co-solvent Addition:** Doping module or direct, on-line metered addition.

**Parallel Operation Vessels:** 50mL to 50L (with 5 micron Fritts).

**Sample Baskets:** S/S mesh, with lids.

**Cold Trap:** for receiving small molecular materials

**Gas Totalizer:** for measuring total and on-time flow of CO<sub>2</sub> gas

**Power Requirements:** US 110 VAC, single phase, 10 amps (International Power, 220 VAC Available)