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Series KDL 5 Laboratory Plants

modular laboratory plants • flexible for all distillation tasks • single or multi-stage

Our modular KDL 5 plants offer unmatched flexibility, with several hundred plants distilling a huge variety of substances all over the world. The modular concept includes a large range of proven components, which are continuously improved. Each project carefully selects the components to ensure a laboratory distillation plant that will run reliably and is easy to operate.


The main component of all KDL 5 plants is the Short Path Distillator with an internal condenser made of borosilicate glass, or

the Thin Film Evaporator with an external condenser. Both types are equipped with an evaporator surface area of 0.05 m².

The secondary components required for continuous operation are compatible to these distillators and evaporators. For the distillation of highly viscous or even solid (at ambient temperature) products, the components which are in direct contact with the products are also available in heatable versions.

Two-stage laboratory plant, series KDL 5, with Thin Film Evaporator and Short Path Distillator

The modular KDL 5 laboratory plants can also be assembled using combinations of Thin Film Evaporators and Short Path Distillators – i. e. for multi-stage plants.

In the two-stage plant shown above, a Thin Film Evaporator strips off the light volatile fraction (solvent) before the main fraction is separated from the heavy boiling residue in the subsequent Short Path Distillator stage.

Components of the KDL 5 modular laboratory plant

Main features of the KDL 5 laboratory plants

Main components	Short Path Distillator Thin Film Evaporator
Evaporation surface area	0.05 m ²
Heating range: Evaporation area	up to 350°C
Residue discharge	up to 250°C
Distillate discharge	up to 200°C
Typical feed rates	0.5 - 1.5 kg/h
Operating pressure	
Short Path Distillator	from 0.001 mbar
Thin Film Evaporator	from 1 mbar
Feed systems	Dosing vessel Gear pump Dosing pump
Discharge systems for distillate and residue	Glass bulb Cut measuring vessel Sample collecting carousel Gear pump
Vacuum systems	Rotary vane pump Oil diffusion pump
Heating devices and heating / cooling devices	Large selection, depending on temperature range and required capacity
Racks	Table rack Mobile floor rack

KDL 5 laboratory plant - Table rack version

KDL 5 laboratory plant - Mobile floor rack version

KDL 5 CADi Laboratory Plant (Computer Aided Distillation)


automated distillation series • up to 5 distillate and residue samples • suitable e.g. for distilling mineral oil residues

Routine distillation processes are performed in some areas of process industry, such as determining distillation curves where the cut ratio (the ratio of the distillate quantity to the feed quantity) is shown as a function of the distillation temperature. Distillation curves are important for process management in the mineral oil industry, where Short Path Distillation plants are widely used for the distillation of (high boiling) vacuum residues.

KDL 5 CADi for automated routine distillations

The KDL 5 CADi Short Path Distillation lab plant is ideal for automated distillation processes. A PLC system controls the distillation sequence, with data input and output for up to 5 sets of distillation parameters done by a PC.

KDL 5 CADi system operation via mouse and monitor

The process sequence of the automatic program includes:

- Positioning the sample collecting carousels to “intermediate runs”
- Heating the distillator up to the operating temperature for the first cut
- Positioning the sample collecting carousels to “samples No. 1”
- Collecting distillate and residue in sample cups for a set period of time
- Positioning the sample collecting carousels to “intermediate runs”
- Heating the distillator to the operating temperature for the second cut, etc.

After termination of the automatic sequence, the sample collecting carousels contain 5 distillate samples and 5 residue samples that are used for analytical evaluations. The distillation parameters are recorded by the plant operating system and printed out.

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UIC GmbH