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# The Instructions of TLC Plate Cutter

Name: TLC Plate Cutter


Cat. No.: CS-E0019

Material: Carbide scribe, Aluminum board, PVC plastic board

Application: Scoring and cutting glass TLC plate is a routine in many chemical laboratories, either to economize on plate usage or cut the plate precisely.

## I The Cutting Operation (Size 20\*20 cm as an example)

1. Place a TLC plate on the plate holder with the glass surface up, the plate must be placed firmly at the top edge of the plate holder and cutter body (leave no gap).
2. Figure1-1: carbide scribe mounted in movable plastic head; Figure1-2: There is a small white marker at the bottom, which corresponds to the position of the carbide scribe. Move the cutter head to start position (top of the plate cutter) without pressing down.
3. Align the scale with the white marker according to the required size. Push the cutter head down lightly until it stops at the lowest position and move it slowly along the plate. While scoring the glass surface, place one hand firmly on the plate to prevent it from moving.
4. After longitudinal cutting, take out the entire plastic plate, rotate 90 degrees, the horizontal scale is facing the position of Figure 1-2, Repeat step 3.
5. Remove the scored glass plate from the plate holder and snap it by hand.
6. Scoring and cutting glass TLC plate of other specifications, Repeat step1-5;

- Figure1 -

## I Replacement of the carbide scribe

It is necessary to replace the carbide scribe, once there is a noticeable deterioration in cutting efficiency.

- Figure2 -

- Figure3 -

1. Loosen two screws (Figure 1) on the slider and remove the cutter head.
2. Remove three screws of the cutter head, then replace the scribe (Figure 2);
3. When a new cutter head is reassembled, check the small spring is positioned correctly between the cutter head and the slider, then tighten both screws on the slider (Figure 3).

If you need to replace the cutter head during use, you can contact ChemScene to purchase and we will also provide installation instructions.

## I Note

1. Glass may have sharp edges! Please handle with care.
2. If the total thickness of coating and glass exceeds 2 mm, place the plate on clean paper instead of the plate holder
3. To prolong the life time of the scribe, the end position of scoring must be set 1-2 mm shorter than the glass edge by adjusting screw (Figure1-3).