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

## **HPR-Series SCF Extraction, Reaction, Cleaning and Drying System**

**HPR-Series SCF Extraction, Reaction  
Cleaning and Drying System**

### **Bench Top SCF units for Universities and Industry**

- Sample vessel Volume from 5mL to 8000mL
- Operating pressure up to 10,000 psi (68.9Mpa)
- Operating Temperature up to 350°C
- Safety Rapture Disc Assembly
- PID Control of Pressure and Temperature
- PID Control of Flow Rate for SFT-10 pump
- Integrated Fluid Preheater and Flow Meter
- Thermoelectrically Cooled Pump Technology
- Various Extract Collection Options
- Optional Co-Solvent Addition Modules
- Optional Vessel Windows (up to 200°C)
- Optional Magnetic Drive Mixing
- Optional Data and Drawing Acquisition System

The SFT-HPR-Series Supercritical Fluid Extraction (SFE), Reaction (SFR), Cleaning (SFC), and Drying (SFD) systems are entry level system which possesses many features typically found in more costly the equipment. It may be used for a variety of applications from routine analytical work to basic process development.

The SFT HPR-Series was developed for people who want to investigate the feasibility of applying supercritical fluid techniques to a wide variety of analytical and processing problems. In addition to its numerous industrial uses, the SFT HPR-Series is well suited to the needs of colleges and universities. It is so affordable that it may be incorporated into teaching laboratories. At the same time, it is capable enough to be used for serious research.

The SFT HPR-Series accommodates 5 ml to 8000 ml extraction vessels. It may be operated at pressures up to 10,000 psi. (68.9 MPa) and at temperatures ranging from ambient to 350°C. It is possible to have up to four vessels installed. The wide range of vessel volumes available makes the SFT HPR-Series well suited to both analytical

scale SF applications and basic process development work. With a 8000 ml vessel, the SFT HPR-SERIES can extract very low levels of key components from materials and process larger amounts of bulk material than would be possible with smaller, analytical scale SF equipment. Inside the SFT HPR-Series oven, a pre-heater ensures that the temperature of the fluid reaching the extractions vessel is controlled precisely. This is essential to obtaining accurate, repeatable results.

The SFT HPR-Series utilizes a high performance, dual piston pump which produces the high pressures required for supercritical fluid work. The system has built-in safety precautions to prevent accidental over-temperature or overpressure conditions. As an additional safety back up, a rupture disc assembly provides mechanical protection against accidental over pressurization of the system.

Manually operated valves ensure long term, maintenance free operation. An integrated program logic controller monitors and adjusts the fluid flow rate to achieve and maintain a desired pressure set point. A PID temperature controller monitors and maintains the precise fluid

temperature inside the high pressure vessel.

The SFT HPR-Series utilizes the latest variable restrictor valve (back pressure regulator) technology, providing precise control over the flow rate of the expanding gas. This BPR was graduated for metering and setting flow rate precisely. This is essential for obtaining highly reproducible results. Flow rates can range from 0.01 to 25 ml/min (0.008 to 18 grams/ min) of liquid CO<sub>2</sub> under typical operating conditions for 500mL less vessels, while 0.1 to 2500 ml/min (0.08 to 2000 grams/ min) for 500mL greater vessels. As carbon dioxide is the most commonly used solvent, the SFT HPR-Series, with some modification, allows the user flexibility to work with a variety of supercritical fluids.

Extract collection options include:

- solid phase extraction (SPE) cartridges
- solvent filled vessels
- fractional cyclone separators
- and EPA sample vials.

Optional co-solvent addition modules, Vessel Windows (up to 200°C), Magnetic Drive Mixing are available for the SFT HPR-Series.

## HPR-Series SCF System Specifications

### Standard Configuration

**Temperature and Pressure Display:** Independent LED displays.

**Temperature Range:** Ambient to 350°C.

**Temperature Precision:** +/- 0.5°C.

**Operating Pressure:** 10,000 psi upper pressure limit.

Front keypad control, with LED display. "Constant pressure" mode of operation.

**Flow Rates at 10,000psi pressure difference:**

- 0.01 to 25 ml/min (0.008 to 18 grams/ min) of liquid CO<sub>2</sub> under typical operating conditions for 500mL less vessels, liquid CO<sub>2</sub> (+/- 2% accuracy). 0 psi flow shift achieved.
- 0.1 to 2500 ml/min (0.08 to 2000 grams/ min) of liquid CO<sub>2</sub> under typical operating conditions for 500mL less vessels, liquid CO<sub>2</sub> (+/- 2% accuracy)

**Over Pressure Safeguard:** High / Low pressure alarms and rupture disc protection.

**Electrical Driven Diaphragm Compressor:** Double ended diaphragm compressor. Dual aluminum heads, furon seals and sapphire pistons, integrated thermoelectric cooling, cam-driven pump mechanism with single stepper motor drive, dual ball and seat check valves (ruby ball, sapphire seat), stainless steel fluid path, prime-purge valve, and pressure transducer. The pump's constant pressure mode features a selectable pressure set point. Flow rate will auto-adjust to maintain pressure. With Peltier electrical cooling at -4°C.

**Restrictor Valve:** Heated up to 250°C; resistant to blockage, with graduated micrometer.

**Preheater:** Improves temperature consistency of the fluid by heating the fluid before it reaches the main pressure vessel.

### Extraction Vessel Sizes:

- 50, 100, 500 mL
- 1, 2, 5, 8 L

Vessels come with 5 micron fritts and CO<sub>2</sub> flow direction controller, interchangeable or parallel operation.

### Vessel Sizes for Parallel and Series Operation:

- 10, 25 50, 100, 500 mL
- 1, 2, 5, 8 L

Vessels come with 5 micron fritts and are interchangeable or parallel operation.

**Dimensions:** depending on the configuration

**Weight :** depending on the configuration

**Compatible with RoHS standard**

### Configuration Options

**Co-solvent Addition:** Doping module or direct, on-line metered addition.

**Parallel Operation Vessels:** 10mL to 8L (with 5 micron Fritts).

**Sample Baskets:** S/S mesh, with lids.

**Cold Trap:** for receiving small molecular materials

**Gas Totalizer:** for measuring total and on-time flow of CO<sub>2</sub> gas

**Mixer:** Magnetic Driven Mixing System

**Software:** Computer Data Acquisition System

**Power Requirements:** US 110 VAC, single phase, 10 amps (International Power, 220 VAC Available)

**Liquid Gas Supply:** Liquid CO<sub>2</sub> Cylinder with Dip Tube.