Quanterix™
The Science of Precision Health

Quanterix SP-X™ Imaging and Analysis System

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


Quanterix
The Science of Precision Health

Simoa™ CorPlex™ Human Cytokine Panel 1
Validation Report for SP-X Item # 85-0329

Table of Contents


1. Introduction	Page 3
2. Standard Curve Ranges	Page 3
3. Standard Curve Plots	Page 4
4. Sensitivity and Dynamic Range	Page 4
5. Inter-Run Precision	Page 5
6. Intra-Run Precision	Page 6
7. Lot-to-Lot Precision	Page 6
8. Spike & Recovery	Page 7
9. Linearity of Dilution	Page 8
10. Matched Serum/Plasma Sample Testing	Page 9
11. Comparison to NIBSC Standards	Page 12
12. Freeze/Thaw Stability	Page 13
13. Single-Plex Correlation	Page 14
14. Specificity - Unrelated Antigens	Page 14
15. Cross-Reactivity in Assay Buffer	Page 15
16. Cross-Reactivity of Individual Detection Antibodies in Sample Matrix	Page 17
17. Admixture Linearity	Page 18

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Simoa® Planar Homebrew 工作流程

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Simoa® Planar Homebrew 启动套装 (Cat#: 100-0461) 包含

- 3 plates with anchor spot
- Capture Tag
- Biotin
- Antibody controls
- Reagents and Buffers for Conjugations and Assays

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
- **Human Angiogenesis Panel 1**
(ANG-2, FGFb, HB-EGF, HGF, PLGF, VEGF, VEGF-C)

The SP-X and Planar Array Assays

Ultra-sensitive cytokine profiling for broad monitoring of CAR-T cell activation and risk of cytokine release syndrome

- **Human Chemokine Panel 1** (IP-10, ITAC, MCP-1, MIP-3B)
- **Homebrew** for in-house custom assay development

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- 01.Hillier RJ, et al.. Aqueous Humor Cytokine Levels and Anatomic Response to Intravitreal Ranibizumab in Diabetic Macular Edema. *JAMA ophthalmology*. 2018;136:382-388.
- 02.Innocenti F, et al.Genetic variation determines VEGF-A plasma levels in cancer patients. *Sci Rep*. 2018;8:16332.
- 03.ayson GC, et al. Plasma Tie2 is a tumor vascular response biomarker for VEGF inhibitors in metastatic colorectal cancer. *Nature communications*. 2018;9:4672.
- 04.Kudo M, et al. Regional Differences in Efficacy, Safety, and Biomarkers for Second-Line Axitinib in Patients with Advanced Hepatocellular Carcinoma: From a Randomized Phase II Study. *Liver cancer*. 2018;7:148-164.
- 05.Salem A, et al. Cell Death, Inflammation, Tumor Burden, and Proliferation Blood Biomarkers Predict Lung Cancer Radiotherapy Response and Correlate With Tumor Volume and Proliferation Imaging. *Clinical lung cancer*. 2018;19:239-248.e7.
- 06.Young PA, et al. Activity of Anti-CD19 Chimeric Antigen Receptor T Cells Against B Cell Lymphoma Is Enhanced by Antibody-Targeted Interferon-Alpha. *Journal of interferon & cytokine research : the official journal of the International Society for Interferon and Cytokine Research*. 2018;38:239-254.
- 07.Clarke JM, et al. A phase Ib study of the combination regorafenib with PF-03446962 in patients with refractory metastatic colorectal cancer (REGAL-1 trial). *Cancer chemotherapy and pharmacology*. 2019.
- 08.Felfeli T, et al. Aqueous Humor Cytokines and Long-Term Response to Anti-Vascular Endothelial Growth Factor Therapy in Diabetic Macular Edema. *American journal of ophthalmology*. 2019.
- 09.George DJ, et al. Phase 1b trial of docetaxel, prednisone, and pazopanib in men with metastatic castration-resistant prostate cancer. *The Prostate*. 2019;79:1752-1761.
- 10.Gochuico BR, et al. A Comprehensive, Multidisciplinary, Precision Medicine Approach to Discover Effective Therapy for an Undiagnosed, Progressive, Fibro-inflammatory Disease. *Translational Research*. 2019.
- 11.Marti FEM, et al. Novel phase I trial design to evaluate the addition of cediranib or selumetinib to preoperative chemoradiotherapy for locally advanced rectal cancer: the DREAMtherapy trial. *European journal of cancer (Oxford, England : 1990)*. 2019;117:48-59.
- 12.Scheiman J et al. Meta-omics analysis of elite athletes identifies a performance-enhancing microbe that functions via lactate metabolism. *Nature medicine*. 2019.**
- 13.Tobos C, et al. Sensitivity and binding kinetics of an ultra-sensitive chemiluminescent enzyme-linked immunosorbent assay at arrays of antibodies. *Journal of Immunological Methods*. 2019:112643-undefined.**
- 14.Edwards MR, et al. Low-dose 17 α -ethinyl estradiol (EE) exposure exacerbates lupus renal disease and modulates immune responses to TLR7/9 agonists in genetically autoimmune-prone mice. *Scientific Reports*. 2020;10:5210.
- 15.Fu W, et al. An open-label, randomized trial of the combination of IFN- γ plus TFF2 with standard care in the treatment of patients with moderate COVID-19. *EClinicalMedicine*. 2020.**
- 16.Fu W, et al.. A clinical pilot study on the safety and efficacy of aerosol inhalation treatment of IFN- γ plus TFF2 in patients with moderate COVID-19. *EClinicalMedicine*. 2020.**
- 17.Tobos CI, Sheehan AJ, Duffy DC and Rissin DM. A Customizable Multiplex Antibody Array Immunoassay with Attomolar Sensitivities. *Anal Chem*. 2020.**
- 18.Vessillier S, et al. Development of the first reference antibody panel for qualification and validation of cytokine release assay platforms- Report of an international collaborative study. *Cytokine: X*. 2020:100042.
- 19.Vignoli A, et al. Effect of Estrogen Receptor Status on Circulatory Immune and Metabolomics Profiles of HER2-Positive Breast Cancer Patients Enrolled for Neoadjuvant Targeted Chemotherapy. *Cancers*. 2020;12:314.
- 20.Dorgham K, et al. Distinct cytokine profiles associated with COVID-19 severity and mortality. *The Journal of allergy and clinical immunology*. 2021.
- 21.Quintanilha JCF, et al. Plasma levels of angiopoietin-2, VEGF-A, and VCAM-1 as markers of bevacizumab-induced hypertension: CALGB 80303 and 90401 (Alliance). *Angiogenesis*. 2021.
- 22.Tan L, et al. Dietary saturated fatty acid palmitate promotes cartilage lesions and activates the unfolded protein response pathway in mouse knee joints. *PLoS One*. 2021;16:e0247237.

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