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Sterility Test System

工作原理 Work Principle

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

SM86 Sterility test system is based on the principle of Membrane Filtration by the action of directional peristaltic pressure, filling testing sample into the closed canister and filtration by the positive pressure of the device. Rinsing membrane to remove the bacteriostatic ingredient if samples included. Appropriate medium is pumped into each canister separately to promote the growth of the captured organisms. The canisters are incubated and examined for contamination in accordance with the relevant pharmacopoeia. Many types of sterile products can be tested, includes antibiotics and drugs containing antibacterial ingredients, infusion, injection, sterilization medical equipment, sterile water for injection, etc. First method of Sterility test in Chp2015, Chinese veterinary pharmacopoeia.

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1. Whole mirror polished stainless steel process, easy to clean and disinfection.
2. 7inch capacitive touch screen, supply good user experiences for information display and input.
3. The unique label print function, connect with the printer, can print the relevant information of the experiment to the sticker and stick it on the canister body directly, which avoids the traditional method of hand-writing on the canister body, improves the standardization of the labeling of experimental data.
4. Two operation mode: Program operation and Direct operation.
5. Equipped with medical waterproof footswitch, intermittent operation and continuous operation of 2 modes.
6. Run is quickly and speed is stable.
7. Drain tray can be hung and lay on the table. When it is hung, more than 180° angle for choice.
8. Side installed with pushbutton power switch with light.
9. Bottle holder can meet the demand on hanging for the pouch.
10. English and Chinese operation language for choice

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- ◆ 小型外接标签打印机, 单线与主设备连接, 轻松自动精准的标签打印。
 - ◆ 脚踏开关与操作系统界面泵启停的配合使用, 操作起来得心应手。脚踏开关的间断运行模式和连续运行模式, 大容量和小容量样品都能容易控制泵的运行。
 - ◆ 防水航空插头, 接触性好, 接口更安全。
- ◆ External min. label printer, single line connect with device, print easily and accurately.
◆ Footswitch works with operation system interface to control the operation mode: Intermittent operation and continuous operation. Large or small capacity samples can be controlled pump running well.
◆ Waterproof aerial plug, good contact, more secure interface.

独创培养器标签打印功能 Unique label print function

- ◆ 原先无菌检查结束后, 培养器上手写标注相关信息, 标注不规范, 而且费时。
 - ◆ 独创培养器标签打印功能, 一键完成, 打印所需标注信息分别贴到相应培养器杯壁上。
 - ◆ 培养时相应培养器读取信息方便。
- ◆ As traditional method after finish sterility test, canister body will be written with information by hand.
◆ GEEVO unique label print function, will print the relevant information of the experiment to the sticker and stick it on the canister body directly,
◆ Be convenient for reading the information of canister during incubator.

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
- ◆ Extra-large (7") industrial-grade IPS anti-glare and anti-uv capacitive touch screen, more convenient for data entry and more sensitive for operation.
◆ Bottle holder can 360° degree swing and up & down move to choose suitable position. also meet the demand of hanging pouch.
◆ Drain tray can be disassemble to be drain tray container and drain tray support, is convenient for clean and disinfection, avoiding pollution. It can be hung or lay on the table. When it is hung, more than 180 degree angle for choice to meet requirement for different operation habit of operator.
◆ Side installed with waterproof pushbutton switch with light.
◆ Rising design on front of the device to avoid the flow of liquid and protect the screen.

技术参数 Technology Parameter

| | |
|--------------------------------------|--|
| 电源 Voltage | 220V±110V±10%, 50/60Hz |
| 总功率 Total power | 120W |
| 电源开关 Pushbutton power switch | 防水按钮式开关, 带亮灯 / Waterproof button switch with light |
| 打印机插口 Printer inlet | 5 芯防水航空插口 / 5 core waterproof aerial plug |
| 脚踏开关插口 Footswitch inlet | 3 芯防水航空插口 / 3 core waterproof aerial plug |
| 转速范围 Rotation speed | 最高转速 240 转/分钟 / Up to 240rpm |
| 预存样品 SOP 信息 Prestore SOP information | 50 个 pcs |
| 显示屏 Touch screen | 7 英寸高清触摸式电容屏 / 7inch high definition touch capacitive screen |
| 操作语言 Operating language | 简体中文、英文 / Chinese \ English |
| 机箱 Pump housing | 304 不锈钢镜面抛光 / 304 stainless steel mirror polished |
| 培养器托板 Drain tray support | 316L 不锈钢镜面抛光 / 316L stainless steel mirror polished |
| 排液槽底槽 Drain tray container | 聚苯砜 (PPSU) |
| 主机总量 Weight of pump device | 14.8Kg |
| 主机尺寸 (mm) Pump device size | 610×380×420 (W×D×H) |
| 打印机尺寸 (mm) Printer device size | 133×186×128 (W×D×H) |
| 标签尺寸 (mm) Print paper size | 50×25 (W×L) |

集菌培养器 Sterility Test Closed Canister

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It works together with SM86 sterility test device.

1. Type choice depends on the testing sample package(LVP glass bottle, plastic bottle, soft pouch, vial and ampoule) and feature of testing sample(Water agent, oil or powder, no bacteriostasis or bacteriostasis, soluble or insoluble, etc)
2. Double and triple for choice
3. Blister box package, avoid closed canister being damaged during transportation effectively.
4. Sterilization mode: EO and Gamma.

培养器专用托盘 Pallet for Closed Canister

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PALLET-10 (10 孔) PALLET-20 (20 孔)

During handle the closed canister after the sterility test, place it on the worktable and incubator for culturing always have the trouble of overturn. And this one can solve it and be easy for operator to check.

1. 304 stainless steel
2. Rounded rectangular openings for portable, easy to take it from incubator.
3. Two types: Pallet-10 and pallet-20
Pallet-10 (10 holes) Pallet-20 (20holes)

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ADD.: No.509 Hongfeng Road, Yuhang Economic & Technical Development Area, Hangzhou, Zhejiang, China

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Hangzhou Geevo Technology Co., Ltd.