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

135V; 10ms; 2mm

参考文献 (REFERENCE)

1. Kubinieć RT, Liang H, Hui SW: Effect of pulse length and pulse strength on transfection by electroporation. *BioTechniques* 1990, 8:16-20.
2. Watanabe SY, Albsoul-Younes AM, Kawano T, Itoh H, Kaziro Y, Nakajima S, Nakajima Y: Calcium phosphate-mediated transfection of primary cultured brain neurons using GFP expression as a marker: application for single neuron electrophysiology. *Neurosci Res* 1999, 33:71-78.
3. Teissie J, Golzio M, Rols MP: Mechanisms of cell membrane electroporation: a minireview of our present knowledge. *Biochim Biophys Acta* 2005, 1724:270-280.
4. Mir LM, Bureau MF, Gehl J, Rangara R, Rouy D, Caillaud JM, Delaere P, Branellec D, Schwartz B, Scherman D: High-efficiency gene transfer into skeletal muscle mediated by electric pulses. *Proc Natl Acad Sci USA* 1999, 96:4262-4267.
5. Colleoni S, Donofrio G, Lagutina I, Duchi R, Galli C, Lazzari G: Establishment, differentiation, electroporation, viral transduction, and nuclear transfer of bovine and porcine mesenchymal stem cells. *Cloning Stem Cells* 2005, 7:154-165.
6. Peister A, Mellad JA, Wang M, Tucker HA, Prockop DJ: Stable transfection of MSCs by electroporation. *Gene Therapy* 2004, 11:224-228.
7. Ahmed S, Reynolds BA, Weiss S: BDNF enhances the differentiation but not the survival of CNS stem cell-derived neuronal precursors. *J Neurosci* 1995, 15:5765-5778.
8. Lu P, Jones LL, Tuszynski MH: BDNF-expressing marrow stromal cells support extensive axonal growth at sites of spinal cord injury. *Exp Neurol* 2005, 191:344-360.

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V-LEADER BIOTECHNOLOGY(BEIJING)CO.,LTD