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WAVE BIOREACTOR

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Web: www.jichen-biotech.com.cn/ www.jichen-biotech.com

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
The wave bioreactor is designed for such as cell therapy, antibody drug seed amplification, and vaccine production. It is suitable for various types of cell cultures, including CHO, NSO, hybridoma, HEK293, insect/baculovirus, adenovirus, T cells, plant cells, and primary human cell lines. By controlling parameters such as swing speed, swing angle, gas flow rate, pH, DO, temperature, etc., using a disposable reaction bag for cell culture results in stable performance, safety, and simplicity.

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The single-use cell culture bag is mounted on a rocking platform of the Wave Bioreactor and inflated for use. The culture medium and cells are delivered to the single-use cell culture bag through tins tubing. The motor movement in the rocking platform causes the cell culture fluid to produce a wavy rocking motion, thus achieving efficient gas transport and mixing. The motion speed and angle of the rocking platform are adjusted according to different cell types and culture conditions.

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1. suitable for most cell culture applications
2. Applicable to batch culture, fed-batch culture, and perfusion culture
3. suitable for suspension culture, microcarrier culture and Fibr disk culture
4. Single-use technology ensures sterility and greatly reduces the risk of contamination
5. 3 tray models to meet the different culture volume
6. swing mode enhances gas exchange and reduces shear force
7. single-use cell culture bags allow flexible switching between products, avoiding cross contamination
8. customizable sensors, software, and single-use cell culture bags to meet specific needs
9. powerful and simple to learn control system provide accurate control of process parameters, process parameter record, comparison and other functions, helping users to optimize the process

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Adapted single-use cell culture bags

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Composed of bag body, pipe, connecting parts and filter, which can be used by customers after irradiation sterilization.

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The type and size of the culture bag specified by the customer, including the length and diameter of the pipeline, the composition of the parts, the number of interfaces, etc.

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Composed of bag body, pipeline, connecting parts, built-in fibra disk and filter, provided for customers after irradiation sterilization;

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/Integrated perfusion filter membrane version

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It is composed of bag body, pipe, connecting piece, retaining film and filter. It can be used by customers after irradiation sterilization.

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
/Equipped with pH\DO electrode (high standard version)

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Composed of bag body, pipe, connector, ReadyMate sterile quick connection, filter and single-use electrode. It is provided to customers after irradiation sterilization

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
Use 50L WAVE to continuously carry out 3 batches of seed cell culture, the 1st batch culture volume is 10L, the 2nd-3rd batch culture is 25L, set up and process control according to the process parameters.

The cell parameters are stable (viable cell density, viability, specific growth rate) and metabolic parameters are stable (YLac/Gluc) after 72 hours of seed cell culture, which can meet the needs of seed cell culture.

细胞参数 Cell parameter	接种密度 / Inoculation density	10 ⁶ cells/mL	0.3
	最大工作体积 / Maximum working volume	L	25
	工作体积 / Working volume	L	10/25
反应器参数 Reactor parameter	温度 / Temperature	设定值 / Set value	37 °C
	转速 / 角度 / Speed / Angle	设定值 / Set value	15rpm/6° / 20 rpm/8°
	pH	设定值 / Set value	7.05±0.2
	DO	设定值 / Set value	40%
	Air	设定值 / Set value	20-50mL/min / 40-70mL/min
	O ₂	设定值 / Set value	0-200mL/min / 0-400mL/min
CO ₂	设定值 / Set value	0-100mL/min	