guava[®] easyCyte

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Guava Caspase Analysis

Unit Control
Sample Information
Analysis Sample List

0001:A01
0002:A02
0003:A03
0004:A04
0005:A05
0006:A06
0007:A07
0008:A08
0009:A09
0010:A10
0011:A11
0012:A12

SR with 7-AAD
FAM with 7-AAD
SR and FAM with 7-AAD

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Y-Axis
SR-Peptide Fluor (YLW-HLog)
7-AAD Fluor (RED-HLog)

Enable Gating

Gating Dot Plot

SR-Peptide Fluor (YLW-HLog) vs Forward Scatter (FSC-HLog)

Marker Position
X: 22
Y: 27
Angle X: 0
Angle Y: 0
Set

Apoptosis Analysis Dot Plot

7-AAD Fluor (RED-HLog) vs SR-Peptide Fluor (YLW-HLog)

Analysis Results

	Count	Cells/mL	% of Total	YLW MFI	RED MFI
Live (LL)	1032	6.72e04	34.40%	6.97	9.66
Mid Apoptotic (LR)	878	5.72e04	29.27%	49.97	12.91
Necrotic (UL)	24	1.56e03	0.80%	16.10	68.21
Late Apoptotic/Dead (UR)	1066	6.94e04	35.53%	123.64	221.78
Caspase +	1944	1.27e05	64.80%	90.37	127.44
7-AAD +	1090	7.10e04	36.33%	121.27	218.40

Analysis Operations

Open Data Set
Log Comment
Export to FCS 2.0
View Event Log
Export to Spreadsheet
Print Preview
Print

Apply Current Settings to Selected Samples

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Title	Source	Year
Tunable Signal Processing in Synthetic MAP Kinase Cascades	Cell	2011
A Rapid and Scalable System for Studying Gene Function in Mice Using Conditional RNA Interference	Cell	2011
Zscan4 regulates telomere elongation and genomic stability in ES cells	Nature	2010
HIF-1 antagonizes p53-mediated apoptosis through a secreted neuronal tyrosinase	Nature	2010
Toolkit for evaluating genes required for proliferation and survival using tetracycline-regulated RNAi	Nature Biotechnology	2010
PI3 kinase signals BCR-dependent mature B cell survival	Cell	2009
TRF2 functions as a protein hub and regulates telomere maintenance by recognizing specific peptide motifs	NAT STRUCT MOL BIOL	2009
Nanog and Oct4 associate with unique transcriptional repression complexes in embryonic stem cells	Nature Cell Biology	2008
SAR11 marine bacteria require exogenous reduced sulphur for growth	Nature	2008
An improved zinc-finger nuclease architecture for highly specific genome editing	Nature Biotechnology	2007
MiR-150 controls B cell differentiation by targeting the transcription factor c-Myb	Cell	2007
Heat shock factor 1 is a powerful multifaceted modifier of carcinogenesis	Cell	2007
Inhibition of CDK1 as a potential therapy for tumors over-expressing MYC	Nature Medicine	2007
A single nucleotide polymorphism in the MDM2 promoter attenuates the p53 tumor suppressor pathway and accelerates tumor formation in humans	Cell	2004
The codon 72 polymorphic variants of p53 have markedly different apoptotic potential	Nature Genetics	2003
Direct recognition of cytomegalovirus by activating and inhibitory NK cell receptors	Science	2002
NF- κ B p100 is a pro-apoptotic protein with anti-oncogenic function	Nature Cell Biology	2002

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guava® HT 系列流式细胞仪	目录编号
guava easyCyte 8HT, Base System	0500-4008B
guava easyCyte 6HT-2L, Base System	0500-4007B
guava easyCyte 6HT, Base System	0500-4005/0500-4006
guava easyCyte 5HT, Base System	0500-4005B
Flow Cell for HT series	0500-2260
Power Conditioner, International	1000-3340
guavaSoft Software Package for guava easyCyte Systems includes InCyte, Express Pro and guavaSuite Modules	0500-4115
Guava® Instrument Cleaning Fluid (ICF) (100 mL)	4200-0140
guava® Easy Check Kit	4500-0025

guava® easyCyte 系列流式细胞仪	目录编号
guava easyCyte 8, Base System	0500-5008B
guava easyCyte 6-2L, Base System	0500-5007B
guava easyCyte 6, Base System	0500-5005/0500-5006
guava easyCyte 5, Base System	0500-5005B
Flow Cell for easyCyte series	0500-2270
Power Conditioner, International	1000-3340
guavaSoft Software Package for guava easyCyte Systems includes InCyte, Express Pro and guavaSuite Modules	0500-5115
Guava® Instrument Cleaning Fluid (ICF) (100 mL)	4200-0140
guava® Easy Check Kit	4500-0025

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