BUEHLER® OmniMet® MHT MOTORIZED HARDNESS TESTING SYSTEM

- Fully or Semi-Automated Indentation Hardness Testing Systems
- Micro or Semi-Macro/ Macro Indentation Hardness Testing
- Supplied Complete with a Hardness Tester or as a Retrofit Kit to an Existing Supported Tester
- Provides Faster
 Answers to Routine
 Tasks Such As
 Effective Case Depth
 and Weld Profile
 Assessment





PRELIMINARY DRAFT



Automated Hardness Testing System Solutions



OmniMet® MHT shown here with twin indenter MicroMet® 5104 and 200mm x 80mm stage

Fully Automated OmniMet® MHT OmniMetMHT provides full automation of the hardness testing process on supported testers with minimal operator involvement for improved productivity.

Using Buehler microindentation hardness tester models MicroMet[®] 2104 or MicroMet[®] 5104, which offer a high degree of motorization, the operater has only to load the specimen, select the required profile of indentation pattern and initiate the test. The automated operation leaves the user free to undertake other work. Upon starting the test, all operations are automated, including rotation of the motorized turret from the objective to the selected indenter, moving the motorized stage to the indent positions, initiating the indent cycle for each indent, rotating the turret back to the measuring objective, focusing and measuring all of the indents, and finally, providing a report of the hardness profile with a determination of the case hardness depth if required. Following this automated operation, the test data may be saved to the MHT database for future review or shared with

co-workers using the supplied sample data viewer software.

Semi-Automated OmniMet[®] MHT

When used on testers with less motorization, such as the MicroMet[®] 2103, MicroMet[®] 5103 or the Buehler semi macro and macro Vickers tester models, the OmniMet MHT may be used in a semi-automated mode with operator involvement to turn a non-motorized turret or to focus the indentations during the measurement process.

Retrofitting to an Existing Tester

The OmniMet MHT system may be supplied complete with a preconfigured tester or as a retrofit kit to an existing tester that is supported.*

* Please contact your Buehler Sales Engineer for the current description, system specification and complete listing of supported accessory options, and hardness testers.

For Micro and Macro Indentation Needs



MicroMet 5104 tester with added motorized stage, automated focus and joystick control for use with the OmniMet MHT Fully Automated Microindentation Hardness Testing System.

Automated Microindentation Hardness Testing

The OmniMet MHT can be supplied for fully automated operation with a single indenter and two objectives or with three objectives and dual indenters for the added advantage of being able to switch quickly between Vickers and Knoop as well as having an intermediate magnification available for larger indents in softer materials.

With the motorized turret and focusing available on these testers, the whole indentation and hardness reading operation is automated. Microindentation hardness tester model MicroMet 5104 is available with the OmniMet MHT as a complete system for fully automated operation.

For semi-automated operation, the MHT system can be supplied with microindentation hardness tester model MicroMet 5103. In the semi automated mode, the user is prompted to turn the non-motorized turret and manually set the focus, all other functions remain automated. A number of additional microindentation



MacroVickers model 1900-5114 with motorized stage and joystick control for semi-automated operation.

testers^{*} are supported by OmniMet MHT for retrofitting for semi-automated operation.

Semi-Automated Macroindentation Hardness Testing

For applications demanding higher test loads, the OmniMet MHT can be supplied complete with a 5112 model semi-macro or 5114 model macro Vickers hardness tester. While most of the hardness testing operation can be automated with these testers, the user is prompted to set the focus manually for semi- automated operation.

A customer's existing macro Vickers or semi-macro Vickers tester that is supported^{*} may be retrofitted with the OmniMet MHT automation package, with the degree of automation being dependent on the model of the tester.

* Please contact your Buehler Sales Engineer for the current description, system specification and complete listing of supported accessory options, and hardness testers.

Generating Accurate Results, Rapidly & Consistently



The OmniMet MHT Graphic User Interface (GUI) is designed for ease of use - with good work flow using a few simply understood icons to run the hardness test. Here the GUI is shown with the alignment tool for setting the direction of the profile perpendicular to the specimen surface.





Once all the indentations have been made and measured, it is then only necessary to select the report style to be used. At left, a case hardness report is shown with a single hardness profile, table of results and two effective case depth limits. Above, a multiple profile report is shown with effective case hardness measurements for both profiles. In all cases, the user may select the unit of measure for depth and for hardness scale as well as selecting the hardness limit for the effective case depth determination.

Secure Data, Easily Retrieved and Exported



🖾 Mie	crosoft Excel - B	ook1								JÞ
E E	jle <u>E</u> dit ⊻jew	Insert Format	Iools Data Window Help					Type a question for	help - I	8 3
		₿ 🖪 💞 🐰 🛛	🖹 🛍 • 🚿 🗠 • 🖓 Σ	- 2↓ 150%	• 🕐 🐥 A	vrial	• 10 • B		1 :00 :00	
	A/ 🗸	∱× 82	C	<u> </u>	1	1	V		ħЛ	-
1	A Sample ID	D Drofile no	Liced profile	Depth	D1	D2	Tectload	Hardnoss	IVI	
2	B 2	1	vickers 10 indent	0.1	25.3	34.5	500	761.5		
2	D 2 D 2	1	vickers 10 indent	0.1	35.9	35.2	500	701.5		
3	D 2 D 2	1	vickers 10 indent	0.3	27.6	26.2	500	679.9		
4	D 2 D 2	1	vickers 10 indent	0.4	37.0 20 E	30.3	500	070.0 629.5		
5	D 2 D 2	1	vickers 10 indent	0.5	30.5	37.7	500	636.5		
0	B2	1	vickers 10 indent	0.6	39.9	40.5	500	5/4.5		
1	B2		vickers 10 indent	0.8	42.4	40.5	500	540.0		
8	B Z	1	vickers 10 Indent	0.9	42.2	42.8	500	514.0		
9	B 2	1	vickers 10 indent	1.0	42.6	43.7	500	497.8		
10	B 2	1	vickers 10 indent	1.1	44.5	44.6	500	467.5		
11	B 2	1	vickers 10 indent	1.3	45.4	45.3	500	451.0		
12	B 2	2	vickers 10 indent	0.1	34.6	37.2	500	718.2		
13	B 2	2	vickers 10 indent	0.3	36.7	37.2	500	678.7		
14	B 2	2	vickers 10 indent	0.4	38.5	36.6	500	658.2		
15	B 2	2	vickers 10 indent	0.5	38.5	37.0	500	650.2		
16	B 2	2	vickers 10 indent	0.6	36.9	39.5	500	634.6		
17	B 2	2	vickers 10 indent	0.8	39.0	39.5	500	601.7		
18	B 2	2	vickers 10 indent	0.9	39.9	39.8	500	584.5		
19	B 2	2	vickers 10 indent	1.0	41.5	42.1	500	531.1		
20	B 2	2	vickers 10 indent	1.1	42.2	43.9	500	500.4		
21	B 2	2	vickers 10 indent	1.3	44.9	43.0	500	479.8		-
22			ļ							

Each OmniMet MHT, fully or semi-automated, comes complete with an OmniMet MHT sample data viewer application. This allows viewing and report generation for the hardness data of any specimen previously stored in the database. The viewer application may be used on the local PC of the OmniMet MHT or on any LAN PC, when the MHT is connected to a network .

Shown at left is a typical view of a populated OmniMet MHT database using the viewer application. The displayed hardness profile graph is from the highlighted data entry for the sample.

Records may be displayed for the whole database, the last months data, or for the last 7 days. Data can also be displayed according to the sample name using the included search tool.

A full report for the highlighted entry may be generated according to various international standards from a choice of supplied templates in Microsoft[®] Access.

Alternatively, the data from the highlighted sample may be exported to Microsoft Excel for the users own data analysis (shown at left).

Automated Hardness System Models and Options Models:

88-7000 OmniMet MHT Fully Automated Microindentation Hardness Testing System. Consisting of an OmniMet PC automation package with choice of: motorized stage and focus, PC accessories, microindentation hardness tester, automated focusing accessory, motorized stage and sample vise.

88-7500 OmniMet MHT Semi-Automated Indentation Hardness Testing System. Consisting of an OmniMet PC semi-automation package with choice of: motorized stage, PC accessories, microindentation hardness tester or macro/semi-macro Vickers hardness tester, motorized stage and sample vise.

PC Accessories:	indenters				
86-1100 Memory Upgrade	Tector Chaice for 88 7500 (MacrolComi Macro);				
86-0160 CD R/W Drive	1900-5112A* Semi-Macro Vickers				
86-0180 DVD R/W Drive	1900-5114A * Macro Vickers				
85-0140 22" Monitor					
85-0141 19" Monitor	Motorized Stage Choice for 88-7000 and 88-7500 Systems:				
85-0142 20" Flat Screen LCD Monitor	88-7050 Motorized stage with 50mm x 50mm travel				
85-6010-160 Laser Printer (Monochrome)	88-7100 Motorized stage with 100mm x 80mm travel				
85-6020-160 Inkjet Printer	88-7200 Motorized stage with 200mm x 80mm travel				
85-6030-160 Laser Printer (Color) Tester Choice for 88-7000:	Add-On to OmniMet Imaging System Options: 88-7001 OmniMet [®] Enterprise				
1600-5122VAMHT * MicroMet 5104 with Vickers indenter	88-7002 OmniMet [®] Image Capture and Report System				
1600-5122KAMHT* MicroMet 5104 with Knoop	88-7003 OmniMet® Express				
indenter	88-7004 OmniMet [®] Archive				
1600-5122VKAMHT * MicroMet 5104 with Vickers and Knoop indenters	*"A" in the catalog number denotes 110V. For 220V, substitute A with D. For 240V, substitute A with E.				

Tester Choice for 88-7500 (Microindentation):

1600-5112VA* MicroMet 5103 with Vickers indenter

1600-5112KA* MicroMet 5103 with Knoop indenter

1600-5112VKA* MicroMet 5103 with Vickers and Knoop

System Specifications

88-7000 OmniMet MHT Fully Automated Microindentation Hardness Testing System and 88-7500 OmniMet Semi-Automated Indentation Hardness Testing System - Please contact your Buehler Sales Engineer for the current description, system specification and complete listing of supported accessory options and hardness testers.

Buehler continuously makes product improvements; therefore, technical specifications are subject to change without notice. For a complete listing of Buehler equipment and consumables please refer to our Buehler Consumables and Equipment Buyers Guides.

© 2004 BUEHLER LTD. Printed in U.S.A. 25M0404 FN0XXXX	*Fo
	BUEHLE In der St Postfach Telefon: Web Site Email: in
BUEHLER BUEHLER BUEHLS - Worldwide Headquarters 41 Waukegan Road • P.O. Box 1	BUEHLE Minipara 3, Chem Télépho Télécopi Web Sita Email: in
Lake Bluff, Illinois 60044-1699 USA Tel: 847/295-6500 • Fax: 847/295-7979 Sales: 1/800/BUEHLER • 1/800/283-4537 Web Site: http://www.buehler.com Email: info@buehler.com	BUEHLE 101, Loc Telephor Web Site Email: s

R GMBH eele 2 • 40599 Düsseldorf 16 03 55 • 40566 Düsseldorf (+49) 0211/974100 • Telefax: (+49) 0211/9741079 e: http://www.buehler-met.de nfo@buehler-met.de R SARL rc de Dardilly, Bát. 0 nin du Jubin • 69570 Dardilly-France ne: (+33) (04) 37 59 81 20 ie: (+33) (04) 37 59 81 29 e: http://www.buehler.fr nfo@buehler.fr ER UK -khurst Lane • Coventry • CV6 5SE ne: (+44) (0) 24 7658 2158 • Fax: (+44) (0) 24 7658 2159 : http://www.buehler.co.uk ales@buebler.co.uk

metallurgical equipment produced by BUEHLER LTD. USA

BUEHLER ANALYST[®] SECTION



BUEHLER CANADA 9999 Highway 48 Markham, Ontario L3P 3|3

Telephone: (905) 201-4686 • Fax: (905) 201-4683 Sales Telephone: 1-800-268-3593 Fmail: info@buehler.ca

BUEHLER ASIA 5/F Vogue Centre 696 Castle Peak Road Lai Chi Kok, Kowloon Hong Kong, SAR, China Telephone: (852) 2307 0909 Fax: (852) 2307 0233