



www.stereozoom.com



# Leica StereoZoom®

Technical information

*Leica*  
MICROSYSTEMS

## Accessories

### Stands, illuminators

Incident-light stands	Antistatic and Terminator
Focusing drives	Coarse, fine, tiltable for OEM applications (bonders), adjustable ease of movement
Focusing range	135mm
Microscope carrier	Optics carrier 360° rotatable
Swingarm stand	Antistatic
FlexArm	For fitting to bench or wall
Transmitted-light stands	 Sub base for incident-light stand, with movable mirror for inclined illumination; transmitted-light stands for bright field and for bright/dark field; high-performance base
Stages	Gliding stage, cup stage, polarizer on glass stage plate
Illuminations	 Leica L2 antistatic cold-light source for linking to stand; various fibre-optic light guides; accessories for coaxial, vertical and transmitted light, Nicholas, Leica LED1000 (Laser-Emitting-Diode illumination), Leica L5 FL fluorescent system



### Accessories

Measuring graticules	For measuring lengths and crosshair, inch and mm
Image rotator for S4 E and S6 models	Vertical/oblique observation all-round 45° bird's-eye view from above
Protective glass for objective	Protects against dirt and damage
Digital imaging Systems (S6 D/S8 APO)	Divers Leica Digital cameras
Image archiving and processing software	Leica Image Manager, Leica QWin, FW4000, Materials Workstation
Video (S6 D/S8 APO)	Divers video objectives and adapters for commercial CCD cameras
Microphoto systems (S6 D/S8 APO)	Leica MPS30 and MPS60, fully automatic, with Databack
SLR camera system (S6 D/S8 APO)	Ricoh XR-X3000D reflex camera with Databack

**Link into the future: [www.stereozoom.com](http://www.stereozoom.com)**

## Technical data, performance features

 Only the Leica StereoZoom® line offers

StereoZoom®	Leica S4 E	Leica S6
Optical system	12°-Greenough, using best-corrected (central) part of objectives	12°-Greenough, using best-corrected (central) part of objectives
Zoom	4.8:1	6.3:1
Viewing angle	38°	60°
ESD protection	antistatic	antistatic
Specific surface resistivity	2×10 <sup>11</sup> ohm/square, discharge time from 1000V to 100V: <2 sec	2×10 <sup>11</sup> ohm/square, discharge time from 1000V to 100V: <2 sec
Magnification range (basic outfit)	6.3×–30×	6.3×–40×
Maximum resolution	372lp/mm	432lp/mm
Maximum numerical aperture	0.124	0.144
Working distance (basic outfit)	110mm	110mm
Field diameter (basic outfit)	 36.5mm	 36.5mm
Adjustable zoom stops	2	2
Video-/photo outlet, switchable		
Photography with coaxial lighting		
Objectives, lead-free	Achromats 0.32×, 0.5×, 0.63×, 0.75×, 1.6×, 2×	Achromats 0.32×, 0.5×, 0.63×, 0.75×, 1.6×, 2×
ErgoObjectives™	0.6×–0.75×/77–137mm 0.7×–1×/48–98mm	0.6×–0.75×/77–137mm 0.7×–1×/48–98mm
Adjustable objective	0.3×–0.4×/200–350mm	0.3×–0.4×/200–350mm
Ergonomic eyepieces, fixed and adjustable with cups	10×/23, 16×/16, 20×/12	10×/23, 16×/16, 20×/12
Ergonomic eyepieces fixed and adjustable for eyeglass wearers, with cups	10×/23, 16×/15, 25×/9.5, 40×/6	10×/23, 16×/15, 25×/9.5, 40×/6
Eyebase	55–75mm	55–75mm

StereoZoom® is a trademark registered in the Principal Register of the US Patent and Trademark Office.

Leica S6 E	Leica S6 T	Leica S6 D	Leica S8 APO
12°-Greenough, using best-corrected (central) part of objectives	12°-Greenough, using best-corrected (central) part of objectives	12°-Greenough, using best-corrected (central) part of objectives	<i>L</i> apochromatic, 12°-Greenough, using best-corrected (central) part of objectives
6.3:1	6.3:1	6.3:1	<i>L</i> 8:1, apochromatic
38°	38°	38°	38°
antistatic	<i>L</i> Terminator, conductive	antistatic	antistatic
2×10 <sup>11</sup> ohm/square, discharge time from 1000V to 100V: <2 sec	10 <sup>2</sup> -10 <sup>6</sup> ohm/square, discharge time from 1000V to zero: <0.1 sec	2×10 <sup>11</sup> ohm/square, discharge time from 1000V to 100V: <2 sec	2×10 <sup>11</sup> ohm/square, discharge time from 1000V to 100V: <2 sec
6.3×–40×	6.3×–40×	6.3×–40×	10×–80×
432lp/mm	432lp/mm	432lp/mm	600lp/mm
0.144	0.144	0.144	0.2
110mm	110mm	110mm	75mm
<i>L</i> 36.5mm	<i>L</i> 36.5mm	<i>L</i> 36.5mm	<i>L</i> 23mm
2	2	2	
		100% visual or 100% video/photo and 100% visual in left eyepiece	100% visual or 100% video/photo and 100% visual in left eyepiece
		<i>L</i> yes	<i>L</i> yes
Achromats 0.32×, 0.5×, 0.63×, 0.75×, 1.6×, 2×	Achromats 0.32×, 0.5×, 0.63×, 0.75×, 1.6×, 2×	Achromats 0.32×, 0.5×, 0.63×, 0.75×, 1.6×, 2×	<i>L</i> Apochromats 0.63×, 1.6×, 2× Achromat 0.32×
0.6×–0.75×/77–137mm 0.7×–1×/48–98mm	0.6×–0.75×/77–137mm 0.7×–1×/48–98mm	0.6×–0.75×/77–137mm 0.7×–1×/48–98mm	
0.3×–0.4×/200–350mm	0.3×–0.4×/200–350mm	0.3×–0.4×/200–350mm	
10×/23, 16×/16, 20×/12	10×/23, 16×/16, 20×/12	10×/23, 16×/16, 20×/12	10×/23, 16×/16, 20×/12
10×/23, 16×/15, 25×/9.5, 40×/6	10×/23, 16×/15, 25×/9.5, 40×/6	10×/23, 16×/15, 25×/9.5, 40×/6	10×/23, 16×/15, 25×/9.5, 40×/6
55–75mm	55–75mm	55–75mm	55–75mm

# Your Line of Sight



Leica S8 APO



Leica S6 D



Leica S6 T



Leica S6




Leica S6 E



Leica S4 E

# System diagram Imaging



## Catalogue references

---

### Stereomicroscopes

---

10 446 293	Leica S4 E StereoZoom® 0.63×–3×, zoom range 4.8:1, 38° viewing angle
10 446 294	Leica S6 E StereoZoom® 0.63×–4×, zoom range 6.3:1, 38° viewing angle
10 446 295	Leica S6 StereoZoom® 0.63×–4×, zoom range 6.3:1, 60° viewing angle
10 446 296	Leica S6 T StereoZoom® 0.63×–4×, conductive, zoom range 6.3:1, 38° viewing angle
10 446 297	Leica S6 D Stereozoom® 0.63×–4×, zoom range 6.3:1, viewing angle 38°, video/photo tube
10 446 298	Leica S8 APO Stereozoom® 1×–8× with apochromatically corrected zoom optics, zoom range 8:1, viewing angle 38°, video/photo tube

---

### Eyepieces, standard

---

10 446 332	Eyepiece 10×/23
10 446 333	Eyepiece 10×/23, adjustable*
10 446 354	Eyepiece 16×/16
10 446 355	Eyepiece 16×/16, adjustable*
10 446 356	Eyepiece 20×/12
10 446 357	Eyepiece 20×/12, adjustable*
10 446 326	Eyepiece 10×/23B f/eyeglasses
10 446 329	Eyepiece 10×/23B adjust. f/eyeglasses*
10 447 059	Eyecups, one pair

---

### Eyepieces, ergonomic

---

	Eyecups included
10 447 130	Eyepiece 10×/23, fixed
10 447 132	Eyepiece 16×/16, fixed
10 447 134	Eyepiece 20×/12, fixed
10 447 131	Eyepiece 10×/23, adjustable
10 447 133	Eyepiece 16×/16, adjustable
10 447 135	Eyepiece 20×/12, adjustable
10 447 149	Eyecups, one pair for eyepieces
10 447 136	Eyepiece for eyeglasses 10×/23B, fixed
10 447 138	Eyepiece for eyeglasses 16×/15B, fixed
10 447 137	Eyepiece for eyeglasses 10×/23B, adjustable
10 447 139	Eyepiece for eyeglasses 16×/15B, adjustable
10 447 150	Eyecups, one pair for eyepiece for eyeglasses
10 445 302	Eyepiece 25×/9.5B adjust. f/eyeglasses*
10 445 303	Eyepiece 40×/6B adjust. f/eyeglasses*
10 447 036	Spacing ring 4mm, pair, required w/ eyep. 10 445 302/303 * adjustable eyepieces accept reticles

---

### S4 E, S6 E, S6, S6 T, S6 D Objectives

---

10 446 316	0.32× Objective, WD 300mm
10 446 318	0.5× Objective, WD 200mm
10 446 319	0.63× Objective, WD 155mm
10 446 320	0.75× Objective, WD 130mm
10 446 321	1.6× Objective, WD 55mm
10 446 322	2× Objective, WD 35mm
10 446 325	Adjustable objective 0.3×–0.4×, WD 200–350mm
10 446 323	ErgoLens 0.6×–0.75×, WD 77–137mm
10 446 317	ErgoLens 0.7×–1×, WD 48–98mm
10 446 324	Lens shield for S4/S6/S8APO

## Low Eyepoint Eyepieces

### Eyepieces, fixed


	Standard	ergonomic
 10x/23	10 446 332	10 447 130
 16x/16	10 446 354	10 447 132
 20x/12	10 446 356	10 447 134

### Eyepieces, adjustable\*

	Standard	ergonomic
 10x/23	10 446 333	10 447 131
 16x/16	10 446 355	10 447 133
 20x/12	10 446 357	10 447 135

\* Adjustable eyepieces accept reticles

## Leica StereoZoom® System diagram





# High Eyepiece Eyepieces

## Eyepieces, fixed

	Standard	ergonomic
10×/23	10 446 326	10 447 136
16×/15		10 447 138

## Eyepieces, adjustable\*

	Standard	ergonomic
10×/23	10 446 329	10 447 137
16×/15		10 447 139
25×/9.5B		10 445 302
40×/6B		10 445 303


\* Adjustable eyepieces accept reticles

10 447 036  
Spacing Ring required with eyepieces  
10 445 302, 10 445 303


Photo / Video documentation  
Leica DC cameras



**Leica S6 T**  
10 446 296



**Leica S6 D**  
10 446 297



**Leica S8 APO**  
10 446 298

## Objectives for S8 APO

0.32×	10 446 334
APO 0.63×	10 446 335
APO 1.6×	10 446 336
APO 2×	10 446 337

## Objectives

### for S4 E, S6 E, S6 T, S6 D

0.32×	10 446 316
0.5×	10 446 318
0.63×	10 446 319
0.75×	10 446 320
1.6×	10 446 321
2×	10 446 322


Lens shield	10 446 324
-------------	------------

## Adjustable Lens


0.3×–0.4×	10 446 325
-----------	------------

## Ergolenses


0.6×–0.75×	10 446 323
0.7×–1×	10 446 317




10 446 345




10 446 344




10 446 343




10 446 342




10 376 070




13 312 611




13 312 613




13 312 614



13 312 610



13 312 714



10 446 299

## Catalogue references

---

### S8 APO Objectives

10 446 335	APO 0.63× objective, WD 100mm, for S8 APO
10 446 336	APO 1.6× objective, WD 37mm, for S8 APO
10 446 337	APO 2× objective, WD 25mm, for S8 APO
10 446 334	0.32× objective, WD 200mm, for S8 APO

---

### Focus columns/arms/carriers

10 446 339	Focus column
10 446 347	Focus column T
10 447 106	Coarse/fine focus drive
10 446 343	Focus arm
10 446 344	Focus arm for 25mm post
10 446 345	Mountable focus arm
10 446 342	STZ adapter for S4/S6/S8 APO
10 446 395	Extension for 0.5× objective for S4/S6
10 446 394	Carrier for c/f focusing drive, S-series
10 376 070	Carrier rod 25mm, inclinable

---

### Bases

10 446 340	Base for S4/S6/S8 APO
10 446 348	Base T
10 446 341	Sub base with reflector
10 445 367	Transmitted-light base HL*
10 446 359	Transmitted-light base HL RC™*
10 445 387	Brightfield base, 20W
10 445 363	Brightfield/darkfield base*

\* The transmitted-light base HL, HL RC™ and BF/DF require an external fiber optic light source (e.g. CLS150) and a fiber optic light guide (e.g. 30150703)

---

### Stands

13 312 714	S-stand
10 446 299	Swing arm stand ESD
13 312 610	Flex Arm Stand
13 312 611	Flex Arm horiz. mount
13 312 613	Flex Arm vert. mount
13 312 614	Flex Arm clamp mount

---

### Stages/Stage plates

10 446 301	Gliding stage
10 446 303	Cup stage
10 447 058	Clear glass stage plate
10 447 153	Frosted glass stage plate
10 447 057	Opaque B&W stage plate
10 447 065	Stage plate S6 T, dissipating

---

### Polarizing accessories

10 446 228	Polarizer on glass stage plate
10 315 306	Analyzer in rotatable mount

---

### Reticles/Stage micrometers

10 446 447	Reticle 10mm/0.1mm
10 446 448	Reticle 5mm/0.1mm
10 446 449	Reticle 5mm/0.05mm
10 447 000	Reticle 100 Div./0.002"
10 447 001	Reticle 100 Div./0.001"
10 447 002	Reticle 150 Div./0.0005"
10 446 450	Crosshair
10 310 345	Stage micrometer 50mm, 0.1/0.01mm div.
10 447 037	Stage micrometer 1", 0.001" divisions

---

### Accessories

10 447 040	Image rotator S4 E/S6
10 447 039	Dust cover, antistatic

## Catalogue references

---

### **Video, photomicrography**

---

For detailed information on Leica camera systems, accessories and image management software, go to our Internet site at [www.leica-microsystems.com](http://www.leica-microsystems.com)

10 445 928	Video objective 0.32× with C-Mount for video-/phototubes
10 445 929	Video objective 0.5× with C-Mount for video-/phototubes
10 446 261	Video objective 0.63× with C-Mount for video-/phototubes
10 446 307	Video objective 0.8× with C-Mount for video-/phototubes
10 445 930	Video-/photoobjective 1× for video-/phototubes
10 445 932	Eyepiece tube 37mm
10 446 120	Photo eyepiece 8×
10 445 304	Photo eyepiece 10×
10 445 305	Photo eyepiece 16×
10 447 003	Framing reticle 10× S6 D/S8 APO

---

### **Illuminations**

---

13 313 305	Nicholas illuminator lamp house
13 410 310	Extendable linkage
13 313 534	Adapter
13 410 311	Nicholas transformer 120V
13 410 312	Nicholas transformer 230V
13 313 175	Bulb for Nicholas lamp

---

### **Leica L5 FL fluorescent system**


---

10 446 431	Leica L5 FL fluorescent system blue (Ex. 450/64 nm / Em. 550/54 nm), 230 V, S Series
10 446 433	Leica L5 FL fluorescent system blue (Ex. 450/64 nm / Em. 550/54 nm), 120 V, S Series
10 446 434	Leica L5 FL fluorescent system green (Ex. 534/40 nm / Em. 622/64 nm), 230 V, S Series
10 446 435	Leica L5 FL fluorescent system green (Ex. 534/40 nm / Em. 622/64 nm), 120 V, S Series

For eye protection an UV protection is absolutely recommended:

10 446 154	UV protection
10 399 211	Arm, clampable
10 446 421	Adapter

# Leica L2 Assembly diagram



## Catalogue references

### Leica L2

	For further details please refer to brochure M1-288-0.
10 446 385	Leica L2 fiber optic light source
10 447 015	Power supply Leica L2
10 446 386	Single flexible light guide, 550mm
10 446 387	Double flexible light guide, 750mm
10 446 388	Single gooseneck, 500mm
10 446 389	Double gooseneck, 500mm
10 446 390	6-point ringlight, 58mm I.D., 750mm
10 447 038	Transmitted light stage
10 446 391	Focusing lens
10 447 055	Day light conversion filter for focusing lens
10 446 392	Universal light guide
10 447 152	Double universal light guide
10 446 374	L2 adapter for focus column
10 446 375	L2 adapter for S-stand
10 446 376	L2 adapter for focus drive 300mm
10 446 377	L2 base
10 446 378	Flexible light guide mount
10 447 009	Universal light guide mount
10 447 152	Double universal light guide
10 446 372	Near vertical illuminator
10 446 373	Coaxial illuminator for S4/S6/S8 APO
10 280 636	Power cable 2m, CH
10 445 661	Power cable 2m, US
10 445 662	Power cable 2m, EURO
10 445 663	Power cable 2m, BS
10 447 056	Bulb 8V/20W for L2

## Information material


Leica StereoZoom® Greenough stereomicroscopes	M1-188-0
Leica L2 cold-light source	M1-288-0
Leica L5 FL cold-light fluorescence system	M1-205-1
Leica Z6 APO and Z16 APO zoom systems	M1-416-0
Leica Z6 APO A and Z16 APO A zoom systems	M1-417-0
Leica MZ16 and MZ16 A stereomicroscope	M1-116-0
Leica MZ125 stereomicroscope	M1-125-0
Leica MZ95 stereomicroscope	M1-195-0
Leica MZ75 stereomicroscope	M1-175-0
Leica MS5, MZ6 stereomicroscope	M1-141-0
Leica MZ FLIII fluorescence stereomicroscope	M1-160-0
Leica Fluo Combi™	M1-166-1
Leica MZ16 FA motorized fluorescence stereomicroscope	M1-116-5
Leica colposcope	M1-280-0
Leica IC A video module for Leica M stereomicroscopes	M1-393-1
Leica ICC A video module for Leica microscopes	M1-393-2
Leica MPS30 photoautomat	M1-330-0
Leica MPS60 photoautomat	M1-360-0
Leica ErgoModules®	M1-215-2
Leica heating stage	M1-227-0

Please also visit our homepage:

**[www.stereomicroscopy.com](http://www.stereomicroscopy.com)**

It provides the latest information and updates and numerous examples for the practical application of our stereomicroscopes in industry and science. You may also view, print and save any desired brochure and the latest manuals in 12 languages.

# Leica LED1000 Assembly diagram



- 30 210 001 Operating unit with power pack
- 30 220 001 LED Spot (19 LEDs)
- 30 221 005 Gooseneck for Spot 85mm
- 30 221 006 Gooseneck for Spot 200mm
- 30 221 007 Gooseneck for Spot 300mm
- 30 221 001 Blue color filter for Spot
- 30 221 002 Green color filter for Spot
- 30 221 003 Yellow color filter for Spot
- 30 221 004 Red color filter for Spot
- 30 221 005 Gooseneck for Spot 85mm
- 30 221 006 Gooseneck for Spot 200mm
- 30 221 007 Gooseneck for Spot 300mm

- 30 221 008 Base plate with M5 adapter for gooseneck
- 30 210 002 LED ring illuminator (40 LEDs) with objective adapter 58mm/66mm
- 30 123 107 Objective adapter 58mm/66mm
- 30 231 001 Polarization filter set for LED ring illuminator
- 30 117 102 Power cable, 3-pole IEC320 for J/USA/CAN
- 30 117 103 Power cable, 3-pole IEC320 for IT
- 30 117 104 Power cable, 3-pole IEC320 for EU
- 30 117 105 Power cable, 3-pole IEC320 for CH
- 30 117 106 Power cable, 3-pole IEC320 for UK
- 30 117 107 Power cable, 3-pole IEC320 for AUS

# Optical Data Leica StereoZoom® S8 APO

				Apochromats						Achromat		
				0.63×		1.6×		2×		0.32×		
Working distance		75mm		101mm		37mm		25mm		200mm		
Eyepieces	Zoom position	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	
<b>10 447 130</b>	1.0	10.0	23.0	6.3	36.5	16.0	14.4	20.0	11.5	3.2	71.9	
<b>10 447 131</b>	1.25	12.5	18.4	7.9	29.1	20.0	11.5	25.0	9.2	4.0	57.5	
<b>10 446 332</b>	1.6	16.0	14.4	10.1	22.8	25.6	9.0	32.0	7.2	5.1	45.1	
<b>10 446 333</b>	2.0	20.0	11.5	12.6	18.3	32.0	7.2	40.0	5.8	6.4	35.9	
<b>10×/23</b>	2.5	25.0	9.2	15.8	14.6	40.0	5.8	50.0	4.6	8.0	28.8	
	3.2	32.0	7.2	20.2	11.4	51.2	4.5	64.0	3.6	10.2	22.5	
<b>10 447 136</b>	4.0	40.0	5.8	25.2	9.1	64.0	3.6	80.0	2.9	12.8	18.0	
<b>10 447 137</b>	5.0	50.0	4.6	31.5	7.3	80.0	2.9	100.0	2.3	16.0	14.4	
<b>10 446 326</b>	6.3	63.0	3.7	39.7	5.8	100.8	2.3	126.0	1.8	20.2	11.4	
<b>10 446 329</b>	8.0	80.0	2.9	50.4	4.6	128.0	1.8	160.0	1.4	25.6	9.0	
<b>10 447 132</b>	1.0	16.0	16.0	10.1	25.3	25.6	10.0	32.0	8.0	5.1	43.9	
<b>10 447 133</b>	1.25	20.0	12.8	12.6	20.3	32.0	8.0	40.0	6.4	6.4	35.0	
<b>16×/16</b>	1.6	25.6	10.0	16.1	15.9	41.0	6.2	51.2	5.0	8.2	27.3	
	2.0	32.0	8.0	20.2	12.7	51.2	5.0	64.0	4.0	10.2	22.0	
	2.5	40.0	6.4	25.2	10.2	64.0	4.0	80.0	3.2	12.8	17.5	
	3.2	51.2	5.0	32.3	7.9	81.9	3.1	102.4	2.5	16.4	13.7	
	4.0	64.0	4.0	40.3	6.4	102.4	2.5	128.0	2.0	20.5	10.9	
	5.0	80.0	3.2	50.4	5.1	128.0	2.0	160.0	1.6	25.6	8.8	
<b>10 446 354</b>	6.3	100.8	2.5	63.5	4.0	161.3	1.6	201.6	1.3	32.3	6.9	
<b>10 446 355</b>	8.0	128.0	2.0	80.6	3.2	204.8	1.3	256.0	1.0	41.0	5.5	
<b>16×/15</b>	1.0	16.0	15.0	10.1	23.7	25.6	9.4	32.0	7.5	5.1	41.2	
	1.25	20.0	12.0	12.6	19.0	32.0	7.5	40.0	6.0	6.4	32.8	
	1.6	25.6	9.4	16.1	14.9	41.0	5.8	51.2	4.7	8.2	25.6	
	2.0	32.0	7.5	20.2	11.9	51.2	4.7	64.0	3.8	10.2	20.6	
	2.5	40.0	6.0	25.2	9.6	64.0	3.8	80.0	3.0	12.8	16.4	
	3.2	51.2	4.7	32.3	7.4	81.9	2.9	102.4	2.3	16.4	12.8	
	4.0	64.0	3.8	40.3	6.0	102.4	2.3	128.0	1.9	20.5	10.2	
	5.0	80.0	3.0	50.4	4.8	128.0	1.9	160.0	1.5	25.6	8.3	
	<b>10 447 138</b>	6.3	100.8	2.3	63.5	3.8	161.3	1.5	201.6	1.2	32.3	6.5
	<b>10 447 139</b>	8.0	128.0	1.9	80.6	3.0	204.8	1.2	256.0	0.9	41.0	5.2
<b>10 447 134</b>	1.0	20.0	12.0	12.6	19.0	32.0	7.5	40.0	6.0	6.4	37.5	
<b>10 447 135</b>	1.25	25.0	9.6	15.8	15.2	40.0	6.0	50.0	4.8	8.0	30.0	
<b>20×/12</b>	1.6	32.0	7.5	20.2	11.9	51.2	4.7	64.0	3.8	10.2	23.5	
	2.0	40.0	6.0	25.2	9.5	64.0	3.8	80.0	3.0	12.8	18.8	
	2.5	50.0	4.8	31.5	7.6	80.0	3.0	100.0	2.4	16.0	15.0	
	3.2	64.0	3.8	40.3	6.0	102.4	2.3	128.0	1.9	20.5	11.7	
	4.0	80.0	3.0	50.4	4.8	128.0	1.9	160.0	1.5	25.6	9.4	
	5.0	100.0	2.4	63.0	3.8	160.0	1.5	200.0	1.2	32.0	7.5	
	<b>10 446 356</b>	6.3	126.0	1.9	79.4	3.0	201.6	1.2	252.0	1.0	40.3	6.0
<b>10 446 357</b>	8.0	160.0	1.5	100.8	2.4	256.0	0.9	320.0	0.8	51.2	4.7	
<b>25×/9.5B</b>	1.0	25.0	9.5	15.8	15.0	40.0	5.9	50.0	4.8	8.0	29.7	
	1.25	31.3	7.6	19.7	12.1	50.0	4.8	62.5	3.8	10.0	23.8	
	1.6	40.0	5.9	25.2	9.4	64.0	3.7	80.0	3.0	12.8	18.6	
	2.0	50.0	4.8	31.5	7.5	80.0	3.0	100.0	2.4	16.0	14.8	
	2.5	62.5	3.8	39.4	6.0	100.0	2.4	125.0	1.9	20.0	11.9	
	3.2	80.0	3.0	50.4	4.7	128.0	1.9	160.0	1.5	25.6	9.3	
	4.0	100.0	2.4	63.0	3.8	160.0	1.5	200.0	1.2	32.0	7.4	
	5.0	125.0	1.9	78.8	3.0	200.0	1.2	250.0	1.0	40.0	5.9	
	6.3	157.5	1.5	99.2	2.4	252.0	0.9	315.0	0.8	50.4	4.7	
	<b>10 445 302</b>	8.0	200.0	1.2	126.0	1.9	320.0	0.7	400.0	0.6	64.0	3.7
<b>40×/6B</b>	1.0	40.0	6.0	25.2	9.5	64.0	3.8	80.0	3.0	12.8	18.8	
	1.25	50.0	4.8	31.5	7.6	80.0	3.0	100.0	2.4	16.0	15.0	
	1.6	64.0	3.8	40.3	6.0	102.4	2.3	128.0	1.9	20.5	11.7	
	2.0	80.0	3.0	50.4	4.8	128.0	1.9	160.0	1.5	25.6	9.4	
	2.5	100.0	2.4	63.0	3.8	160.0	1.5	200.0	1.2	32.0	7.5	
	3.2	128.0	1.9	80.6	3.0	204.8	1.2	256.0	0.9	41.0	5.9	
	4.0	160.0	1.5	100.8	2.4	256.0	0.9	320.0	0.8	51.2	4.7	
	5.0	200.0	1.2	126.0	1.9	320.0	0.8	400.0	0.6	64.0	3.8	
	6.3	252.0	1.0	158.8	1.5	403.2	0.6	504.0	0.5	80.6	3.0	
<b>10 445 303</b>	8.0	320.0	0.8	201.6	1.2	512.0	0.5	640.0	0.4	102.4	2.3	

# Optical Data StereoZoom® Leica S4 E, S6 E, S6, S6 T, S6 D

				with additional objectives										
				0.32×		0.5×		0.63×		0.75×		1.6×		
Working distance		110mm		300mm		200mm		155mm		130mm		55mm		
Eyepieces	Zoom position *4.0: S6	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	
<b>10 447 130</b>	0.63	6.3	36.5	2.0	115.0	3.2	71.9	4.0	57.5	4.7	48.9	10.1	22.8	
<b>10 447 131</b>	0.8	8.0	28.8	2.6	88.5	4.0	57.5	5.0	46.0	6.0	38.3	12.8	18.0	
<b>10 446 332</b>	1.0	10.0	23.0	3.2	71.9	5.0	46.0	6.3	36.5	7.5	30.7	16.0	14.4	
<b>10 446 333</b>	1.25	12.5	18.4	4.0	57.5	6.3	36.5	7.9	29.1	9.4	24.5	20.0	11.5	
<b>10×/23</b>	1.6	16.0	14.4	5.1	45.1	8.0	28.8	10.1	22.8	12.0	19.2	25.6	9.0	
<b>10 447 136</b>	2.0	20.0	11.5	6.4	35.9	10.0	23.0	12.6	18.3	15.0	15.3	32.0	7.2	
<b>10 447 137</b>	2.5	25.0	9.2	8.0	28.8	12.5	18.4	15.8	14.6	18.8	12.2	40.0	5.8	
<b>10 446 326</b>	3.2	32.0	7.2	10.2	22.5	16.0	14.4	20.2	11.4	24.0	9.6	51.2	4.5	
<b>10 446 329</b>	4.0*	40.0	5.8	12.8	18.0	20.0	11.5	25.2	9.1	30.0	7.7	64.0	3.6	
<b>10 447 132</b>	0.63	10.1	25.3	3.2	80.0	5.0	51.2	6.4	40.0	7.6	33.7	16.1	15.9	
<b>10 447 133</b>	0.8	12.8	20.0	4.1	62.4	6.4	40.0	8.1	31.6	9.6	26.7	20.5	12.5	
<b>16×/16</b>	1.0	16.0	16.0	5.1	50.2	8.0	32.0	10.1	25.3	12.0	21.3	25.6	10.0	
	1.25	20.0	12.8	6.4	40.0	10.0	25.6	12.6	20.3	15.0	17.1	32.0	8.0	
	1.6	25.6	10.0	8.2	31.2	12.8	20.0	16.1	15.9	19.2	13.3	41.0	6.2	
	2.0	32.0	8.0	10.2	25.1	16.0	16.0	20.2	12.7	24.0	10.7	51.2	5.0	
	2.5	40.0	6.4	12.8	20.0	20.0	12.8	25.2	10.2	30.0	8.5	64.0	4.0	
<b>10 446 354</b>	3.2	51.2	5.0	16.4	15.6	25.6	10.0	32.3	7.9	38.4	6.7	81.9	3.1	
<b>10 446 355</b>	4.0*	64.0	4.0	20.5	12.5	32.0	8.0	40.3	6.4	48.0	5.3	102.4	2.5	
<b>16×/15</b>	0.63	10.1	23.8	3.2	74.0	5.0	47.6	6.4	37.8	7.6	31.7	16.1	14.9	
	0.8	12.8	18.8	4.1	58.6	6.4	37.5	8.1	29.8	9.6	25.0	20.5	11.7	
	1.0	16.0	15.0	5.1	46.9	8.0	30.0	10.1	23.8	12.0	20.0	25.6	9.4	
	1.25	20.0	12.0	6.4	37.5	10.0	24.0	12.6	19.0	15.0	16.0	32.0	7.5	
	1.6	25.6	9.4	8.2	29.3	12.8	18.8	16.1	14.9	19.2	12.5	41.0	5.8	
	2.0	32.0	7.5	10.2	23.4	16.0	15.0	20.2	11.9	24.0	10.0	51.2	4.7	
	2.5	40.0	6.0	12.8	18.8	20.0	12.0	25.2	9.5	30.0	8.0	64.0	3.8	
	<b>10 447 138</b>	3.2	51.2	4.7	16.4	14.6	25.6	9.4	32.3	7.4	38.4	6.2	81.9	2.9
	<b>10 447 139</b>	4.0*	64.0	3.8	20.5	11.7	32.0	7.5	40.3	5.9	48.0	5.0	102.4	2.3
<b>10 447 134</b>	0.63	12.6	19.0	4.0	60.0	6.3	38.1	7.9	30.4	9.5	25.3	20.2	11.9	
<b>10 447 135</b>	0.8	16.0	15.0	5.1	47.1	8.0	30.0	10.1	23.8	12.0	20.0	25.6	9.4	
<b>20×/12</b>	1.0	20.0	12.0	6.4	37.5	10.0	24.0	12.6	19.0	15.0	16.0	32.0	7.5	
	1.25	25.0	9.6	8.0	30.0	12.5	19.2	15.8	15.2	18.8	12.8	40.0	6.0	
	1.6	32.0	7.5	10.2	23.5	16.0	15.0	20.2	11.9	24.0	10.0	51.2	4.7	
	2.0	40.0	6.0	12.8	18.8	20.0	12.0	25.2	9.5	30.0	8.0	64.0	3.8	
	2.5	50.0	4.8	16.0	15.0	25.0	9.6	31.5	7.6	37.5	6.4	80.0	3.0	
	3.2	64.0	3.8	20.5	11.7	32.0	7.5	40.3	6.0	48.0	5.0	102.4	2.3	
	<b>10 446 356</b>	4.0*	80.0	3.0	25.6	9.4	40.0	6.0	50.4	4.8	60.0	4.0	128.0	1.9
<b>25×/9.5B</b>	0.63	15.8	15.0	5.0	47.5	7.9	30.1	9.9	24.0	11.8	20.1	25.2	9.4	
	0.8	20.0	11.9	6.4	37.1	10.0	23.8	12.6	18.8	15.0	15.8	32.0	7.4	
	1.0	25.0	9.5	8.0	29.7	12.5	19.0	15.8	15.0	18.8	12.6	40.0	5.9	
	1.25	31.3	7.6	10.0	23.8	15.6	15.2	19.7	12.1	23.4	10.1	50.0	4.8	
	1.6	40.0	5.9	12.8	18.6	20.0	11.9	25.2	9.4	30.0	7.9	64.0	3.7	
	2.0	50.0	4.8	16.0	14.8	25.0	9.5	31.5	7.5	37.5	6.3	80.0	3.0	
	2.5	62.5	3.8	20.0	11.9	31.3	7.6	39.4	6.0	46.9	5.1	100.0	2.4	
	3.2	80.0	3.0	25.6	9.3	40.0	5.9	50.4	4.7	60.0	4.0	128.0	1.9	
	<b>10 445 302</b>	4.0*	100.0	2.4	32.0	7.4	50.0	4.8	63.0	3.8	75.0	3.2	160.0	1.5
<b>40×/6B</b>	0.63	25.2	9.5	8.1	29.6	12.6	19.0	15.9	15.1	18.9	12.7	40.3	6.0	
	0.8	32.0	7.5	10.2	23.5	16.0	15.0	20.2	11.9	24.0	10.0	51.2	4.7	
	1.0	40.0	6.0	12.8	18.8	20.0	12.0	25.2	9.5	30.0	8.0	64.0	3.8	
	1.25	50.0	4.8	16.0	15.0	25.0	9.6	31.5	7.6	37.5	6.4	80.0	3.0	
	1.6	64.0	3.8	20.5	11.7	32.0	7.5	40.3	6.0	48.0	5.0	102.4	2.3	
	2.0	80.0	3.0	25.6	9.4	40.0	6.0	50.4	4.8	60.0	4.0	128.0	1.9	
	2.5	100.0	2.4	32.0	7.5	50.0	4.8	63.0	3.8	75.0	3.2	160.0	1.5	
	3.2	128.0	1.9	41.0	5.9	64.0	3.8	80.6	3.0	96.0	2.5	204.8	1.2	
	<b>10 445 303</b>	4.0*	160.0	1.5	51.2	4.7	80.0	3.0	100.8	2.4	120.0	2.0	256.0	0.9









with additional objectives													
2x		0.3x-0.4x		0.3x-0.4x		0.6x-0.75x		0.6x-0.75x		0.7x-1x		0.7x-1x	
35mm		200mm		350mm		77mm		137mm		48mm		98mm	
Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm	Total magnification	Field diameter mm
12.6	18.3	2.5	92.0	1.8	127.8	4.7	48.9	3.5	65.7	6.2	37.1	4.5	51.1
16.0	14.4	3.1	74.2	2.2	104.5	6.0	38.3	4.5	51.1	7.9	29.1	5.7	40.4
20.0	11.5	3.9	59.0	2.8	82.1	7.5	30.7	5.6	41.1	9.9	23.2	7.1	32.4
25.0	9.2	4.9	46.9	3.5	65.7	9.4	24.5	7.0	32.9	12.4	18.5	8.9	25.8
32.0	7.2	6.2	37.1	4.5	51.1	12.0	19.2	9.0	25.6	15.8	14.6	11.4	20.2
40.0	5.8	7.8	29.5	5.6	41.1	15.0	15.3	11.2	20.5	19.8	11.6	14.2	16.2
50.0	4.6	9.8	23.5	7.0	32.9	18.8	12.2	14.0	16.4	24.8	9.3	17.8	12.9
64.0	3.6	12.5	18.4	9.0	25.6	24.0	9.6	17.9	12.8	31.7	7.3	22.7	10.1
80.0	2.9	15.6	14.7	11.2	20.5	30.0	7.7	22.4	10.3	39.6	5.8	28.4	8.1
20.2	12.7	3.9	65.6	2.8	91.4	7.6	33.7	5.6	45.7	10.0	25.6	7.2	35.6
25.6	10.0	5.0	51.2	3.6	71.1	9.6	26.7	7.2	35.6	12.7	20.2	9.1	28.1
32.0	8.0	6.2	41.3	4.5	56.9	12.0	21.3	9.0	28.4	15.8	16.2	11.4	22.5
40.0	6.4	7.8	32.8	5.6	45.7	15.0	17.1	11.2	22.9	19.8	12.9	14.2	18.0
51.2	5.0	10.0	25.6	7.2	35.6	19.2	13.3	14.3	17.9	25.3	10.1	18.2	14.1
64.0	4.0	12.5	20.5	9.0	28.4	24.0	10.7	17.9	14.3	31.7	8.1	22.7	11.3
80.0	3.2	15.6	16.4	11.2	22.9	30.0	8.5	22.4	11.4	39.6	6.5	28.4	9.0
102.4	2.5	20.0	12.8	14.3	17.9	38.4	6.7	28.7	8.9	50.7	5.0	36.4	7.0
128.0	2.0	25.0	10.2	17.9	14.3	48.0	5.3	35.8	7.2	63.4	4.0	45.4	5.6
20.2	11.9	3.9	61.5	2.8	85.7	7.6	31.6	5.6	42.9	10.0	24.0	7.2	33.3
25.6	9.4	5.0	48.0	3.6	66.6	9.6	25.0	7.2	33.3	12.7	18.9	9.1	26.4
32.0	7.5	6.2	38.7	4.5	53.4	12.0	20.0	9.0	26.7	15.8	15.2	11.4	21.0
40.0	6.0	7.8	30.8	5.6	42.9	15.0	16.0	11.2	21.4	19.8	12.1	14.2	16.9
51.2	4.7	10.0	24.0	7.2	33.3	19.2	12.5	14.3	16.8	25.3	9.5	18.2	13.2
64.0	3.8	12.5	19.2	9.0	26.7	24.0	10.0	17.9	13.4	31.7	7.6	22.7	10.6
80.0	3.0	15.6	15.4	11.2	21.4	30.0	8.0	22.4	10.7	39.6	6.1	28.4	8.5
102.4	2.3	20.0	12.0	14.3	16.8	38.4	6.2	28.7	8.4	50.7	4.7	36.4	6.6
128.0	1.9	25.0	9.6	17.9	13.4	48.0	5.0	35.8	6.8	63.4	3.8	45.4	5.3
25.2	9.5	4.9	49.0	3.5	68.6	9.5	25.3	7.1	33.8	12.5	19.2	8.9	27.0
32.0	7.5	6.2	38.7	4.5	53.3	12.0	20.0	9.0	26.7	15.8	15.2	11.4	21.1
40.0	6.0	7.8	30.8	5.6	42.9	15.0	16.0	11.2	21.4	19.8	12.1	14.2	16.9
50.0	4.8	9.8	24.5	7.0	34.3	18.8	12.8	14.0	17.1	24.8	9.7	17.8	13.5
64.0	3.8	12.5	19.2	9.0	26.7	24.0	10.0	17.9	13.4	31.7	7.6	22.7	10.6
80.0	3.0	15.6	15.4	11.2	21.4	30.0	8.0	22.4	10.7	39.6	6.1	28.4	8.5
100.0	2.4	19.5	12.3	14.0	17.1	37.5	6.4	28.0	8.6	49.5	4.8	35.5	6.8
128.0	1.9	25.0	9.6	17.9	13.4	48.0	5.0	35.8	6.7	63.4	3.8	45.4	5.3
160.0	1.5	31.2	7.7	22.4	10.7	60.0	4.0	44.8	5.4	79.2	3.0	56.8	4.2
31.5	7.5	6.1	38.9	4.4	54.0	11.8	20.1	8.8	27.0	15.6	15.2	11.2	21.2
40.0	5.9	7.8	30.4	5.6	42.4	15.0	15.8	11.2	21.2	19.8	12.0	14.2	16.7
50.0	4.8	9.8	24.2	7.0	33.9	18.8	12.6	14.0	17.0	24.8	9.6	17.8	13.3
62.5	3.8	12.2	19.5	8.8	27.0	23.4	10.1	17.5	13.6	30.9	7.7	22.2	10.7
80.0	3.0	15.6	15.2	11.2	21.2	30.0	7.9	22.4	10.6	39.6	6.0	28.4	8.4
100.0	2.4	19.5	12.2	14.0	17.0	37.5	6.3	28.0	8.5	49.5	4.8	35.5	6.7
125.0	1.9	24.4	9.7	17.5	13.6	46.9	5.1	35.0	6.8	61.9	3.8	44.4	5.3
160.0	1.5	31.2	7.6	22.4	10.6	60.0	4.0	44.8	5.3	79.2	3.0	56.8	4.2
200.0	1.2	39.0	6.1	28.0	8.5	75.0	3.2	56.0	4.2	99.0	2.4	71.0	3.3
50.4	4.8	9.8	24.5	7.1	33.8	18.9	12.7	14.1	17.0	24.9	9.6	17.9	13.4
64.0	3.8	12.5	19.2	9.0	26.7	24.0	10.0	17.9	13.4	31.7	7.6	22.7	10.6
80.0	3.0	15.6	15.4	11.2	21.4	30.0	8.0	22.4	10.7	39.6	6.1	28.4	8.5
100.0	2.4	19.5	12.3	14.0	17.1	37.5	6.4	28.0	8.6	49.5	4.8	35.5	6.8
128.0	1.9	25.0	9.6	17.9	13.4	48.0	5.0	35.8	6.7	63.4	3.8	45.4	5.3
160.0	1.5	31.2	7.7	22.4	10.7	60.0	4.0	44.8	5.4	79.2	3.0	56.8	4.2
200.0	1.2	39.0	6.2	28.0	8.6	75.0	3.2	56.0	4.3	99.0	2.4	71.0	3.4
256.0	0.9	49.9	4.8	35.8	6.7	96.0	2.5	71.7	3.3	126.7	1.9	90.9	2.6
320.0	0.8	62.4	3.8	44.8	5.4	120.0	2.0	89.6	2.7	158.4	1.5	113.6	2.1

# Dimensions Leica S4 E, S6 E, S6 T

with incident-light stand



with transmitted-light stand



Dimensions in mm

# Dimensions Leica S6




with mountable focus arm, tiltable






Dimensions in mm

# Dimensions Leica S6 D

## with incident-light stand






## with transmitted-light stand






Dimensions in mm

# Dimensions Leica S8 APO

## with incident-light stand



## with transmitted-light stand



Dimensions in mm

# Leica Microsystems – the brand for outstanding products

Leica Microsystems' Mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement, lithography and analysis of microstructures.

Leica, the leading brand for microscopes and scientific instruments, has developed from five brand names, all with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Leica symbolizes not only tradition, but also innovation.

## Leica Microsystems – an international company with a strong network of customer services

Australia:	Gladesville, NSW	Tel. +1 800 625 286	Fax +61 2 9817 8358
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 20 00	Fax +1 905 762 89 37
China:	Hong Kong	Tel. +8522 564 6699	Fax +8522 564 4163
Denmark:	Herlev	Tel. +45 44 5401 01	Fax +45 44 5401 11
France:	Rueil-Malmaison Cédex	Tel. +33 1 4732 8585	Fax +33 1 4732 8586
Germany:	Bensheim	Tel. +49 6251 1360	Fax +49 6251 136 155
Italy:	Milan	Tel. +39 02 57 486 1	Fax +39 02 5740 3273
Japan:	Tokyo	Tel. +81 3 543 596 09	Fax +81 3 543 596 15
Korea:	Seoul	Tel. +82 2 514 6543	Fax +82 2 514 6548
Netherlands:	Rijswijk	Tel. +31 70 41 32 130	Fax +31 70 41 32 109
Portugal:	Lisbon	Tel. +35 1 213 889 112	Fax +35 1 213 854 668
Singapore:		Tel. +65 6 77 97 823	Fax +65 6 77 30 628
Spain:	Barcelona	Tel. +34 93 494 9530	Fax +34 93 494 9532
Sweden:	Sollentuna	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Glattbrugg	Tel. +41 44 809 34 34	Fax +41 44 809 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 666 663	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 800 248 0123	Fax +1 847 405 0164

## and representatives of Leica Microsystems in more than 100 countries.

In accordance with the ISO 9001 certificate, Leica Microsystems (Switzerland) Ltd., Business Unit Stereo & Microscope Systems has at its disposal a management system that meets the requirements of the international standard for quality management. In addition, production meets the requirements of the international standard ISO 14001 for environmental management.

The companies of the Leica Microsystems Group operate internationally in five business segments, where we rank with the market leaders.

### Microscopy

Our expertise in microscopy is the basis for all our solutions for visualization, measurement and analysis of microstructures in life sciences and industry.

### Specimen Preparation

We specialize in supplying complete solutions for histology and cytopathology.

### Imaging Systems

With confocal laser technology and image analysis systems, we provide three-dimensional viewing facilities and offer new solutions for cytogenetics, pathology and material sciences.

### Medical Equipment

Innovative technologies in our surgical microscopes offer new therapeutic approaches in microsurgery. With automated instruments for ophthalmology, we enable new diagnostic methods to be applied.

### Semiconductor Equipment

Our automated, leading-edge measurement and inspection systems and our E-beam lithography systems make us the first choice supplier for semiconductor manufacturers all over the world.

Leica Microsystems (Switzerland) Ltd.  
Stereo & Microscope Systems  
CH-9435 Heerbrugg

Telephone +41 71 726 33 33  
Fax +41 71 726 33 99  
[www.leica-microsystems.com](http://www.leica-microsystems.com)  
[www.stereomicroscopy.com](http://www.stereomicroscopy.com)

**Leica**  
MICROSYSTEMS