

In-Vivo Electroporator CUY21VIVO-SQ

Revolutionary and innovative technology for in vivo electroporation



Applications

- **In ovo electroporation**
Neural tube / Lim bud / Mesencephalon / Optic vesicle etc.
- **In utero electroporation**
Spinal cord / Hippocampus / Cortex / Telencephalon etc.
- **Ex vivo electroporation**
Mouse's brain slice / Mouse embryo / Mouse embryo's gonad / Mouse's small intestine etc.
- **In vivo electroporation**
Mouse : Muscle / Brain / Skin
Rat : Knee joint / Retina / Kidney
Zebrafish : Fin / Retina
Honey bee : Brain
Planarian : Eye
Xenopus : Embryo
- **Electrochemotherapy**
Chondrosarcoma

Features

Impedance measurement

Impedance measurement function is a unique feature and enables a researcher to measure the impedance prior to electroporation. Under in vivo experiment, many factors change every experiment.

The fluctuation in physical factors such as distance between electrodes and sample volume etc. affects resistance and current. The current is highly related to the reproducibility of in vivo electroporation. Since the resistance and current are interrelated, the impedance value (resistance) is an important indication in in vivo electroporation. With this function, one is allowed to check the impedance (resistance) before running an electroporation program and adjust it by moving the position of electrodes slightly if it is far away from previous value.

Low current measurement

CUY21VIVO-SQ has the current measurement function that measures and displays the actual current immediately after electroporation. The current value is crucial data and helps a researcher to verify the process of electroporation in an electrical aspect. The current measurement function is improved and the minimum current to be measured is extended to 1mA.

The capability enables the in vivo electroporation with very sensitive targets such as adult mouse's brain in which the current must be below 10mA. More applications can be done without any extra devices.

Low voltage control

In in-vivo electroporation, the required voltage is below 100V. Since CUY21VIVO-SQ is designed for in-vivo electroporation, the accuracy of the voltage below 100V is very high due to the narrow voltage range. The voltage of CUY21Vivo-SQ can be set from 0.1V to 99.9V in 0.1V increments. The low voltage control is required for the electroporation with a single cell.

Safety function

As the conductive part of an electrode for in vivo electroporation is exposed and a researcher often grabs an electrode by hand in in vivo electroporation, there is a possibility that a researcher touches it during electroporation accidentally and gets injured seriously. In order to avoid such a tragedy, the safety function is built in CUY21Vivo-SQ and turns off the electroporation program automatically when the current gets beyond the limit.

Easy operation

CUY21VIVO-SQ carries seven independent LED displays and numerical keys. Each electroporation parameter is shown at its own LED display. The program settings and edited measured values can be checked at a glance. The numerical keys also make it easy for a researcher to make and edit a program.

Electrodes

In ovo electroporation

Chopstick electrode

For ectoderm (neural tube)

How to read a catalog number



LF611P3-1
 Length of conductive tip
 Length of bent part
 Platinum coating

Electrode list

LF611P3-1 LF611P7-5
 LF611P7-2 LF611P8-2
 LF611P7-4 (picture)

Z-shape chopstick electrode

For mesoderm and endoderm

How to read a catalog number



LF613P3
 Length of conductive tip
 Platinum coating

Electrode list

LF613P1 LF613P3 (picture)
 LF613P2 LF613P5

Tungsten needle electrode & holder

How to read a catalog number



LF614T
 T : tungsten
 S : Stainless steel

Electrode list

LF614T (picture) LF615C Electrode holder
 LF614S

Electrode holder



LF580 (for LF611 and LF613 series)

In vivo electroporation

Tweezers with disk electrodes

For mouse embryo, muscle, testis, retina, skin and zebra fish's fin etc.

How to read a catalog number



LF650P5
 Diameter of a disk electrode
 P : Platinum coating
 S : Stainless steel

Electrode list

LF650P1 LF650P7 LF650S5
 LF650P3 LF650P10 LF650S7
 LF650P5 (picture) LF650S10

Parallel fixed needle electrodes

For muscle

How to read a catalog number



LF560S5
 Distance between electrodes
 Stainless steel

Electrode list

LF560S5 (picture)
 LF560S10
 LF560S15

Ex vivo electroporation

Petri dish type electrode & Rod type electrode

For sliced tissue

How to read a catalog number



LF701P10E
 E : Petri dish type electrode
 L : Rod type electrode
 Diameter or square
 Platinum
 1 : Square electrode
 0 : Disk electrode

Electrode list

LF700P2E LF700P1L LF701P2E LF701P2L
 LF700P3E LF700P2L LF701P5E LF701P3L
 LF700P4E LF700P3L LF701P7E LF701P5L
 LF700P5E LF700P5L LF701P10E LF701P7L
 LF700P7E LF700P7L LF701P20E LF701P10L
 LF700P10E LF70010L LF701P20L
 LF700P20E LF700P20L

Bath type electrode

For mouse embryo and extracted organ

How to read a catalog number



LF520P5
 Distance between plate electrodes
 Platinum

Electrode list

LF520P5 (picture, L8 x H3mm plate)
 LF520P15 (L10 x H5mm plate)
 LF520P20 (L20 x H8mm plate)
 LF520P25 (L20 x H8mm plate)

There are more kinds of electrodes available and an electrode can be customized. Please kindly contact us.

Specification

DC wave form	Square	Current limit *	1.60A
Impedance measurement range	30ohms - 35Kohms	Voltage measurement range	0.1 - 99.9V
Voltage range	0.1 - 99.9V (0.1V increment)	Current measurement range	0.001 - 1.60A
Pulse length	0.05 - 99.9msec (0.01msec increment)		When measured current < 1.00A, 1mA resolution When measured current ≥ 1.00A, 10mA resolution
Pulse interval	0.1 - 999msec (0.1msec increment)	Power	115V or 220V 250VA 50/60Hz
No. of pulse	1 - 99 pulses	Dimensions	W360cm x L380cm x H180cm
Storage	99 programs	Weight	12.3Kg

* Safety function. When current exceeds the limit, the electroporation program is turned off automatically.

BEX

Manufactured by
BEX Co., Ltd.
 2-61-14 Itabashi Itabashi-ku
 Tokyo 173-0004, Japan

Distributed by
Fukuoka Bio-instrument Trading Co., Ltd.
 2-7-20 Kawarada Onojo
 Fukuoka 816-0932 Japan
 TEL: +81-92-404-0545 FAX: +81-92-404-0546
 Email: info@fbt-corp.jp