"Crushing" Lysis Efficiency - Nothing Resists It!

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电邮：tracy.meng@mpbio.com

FastPrep®-24 Instrument Ultra-Rapid, Thorough Homogenization

Lyse Any Tough or Frozen Sample In 40 Seconds

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
Genomic DNA from human ovarian tissue lysed with the FastPrep®-24 for 20 sec. Courtesy of Dr. David Smith, Oncotech Inc.

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
Overview of the FastPrep®-24 Procedure

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
NEW EASY-TO-USE PROGRAMMABLE INSTRUMENT!

Save up to 5 pre-sets that you use most, with specific combinations of speed and time.

Cat. #	Description	Price
6004500	FastPrep®-24 Instrument	

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

MP Biomedicals guarantees the BEST performance from your FastPrep®-24 Instrument when used in combination with FastPrep Lysing Matrix Tubes

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New Easy Dispenser Box for all Matrix Kits of 100 and 500 Tubes

Cat. #	Description	Pack size
6910-050	Lysing Matrix A	50 x 2 ml Tubes
6910-100	Lysing Matrix A	100 x 2 ml Tubes
6910-500	Lysing Matrix A	500 x 2 ml Tubes
6911-050	Lysing Matrix B	50 x 2 ml Tubes
6911-100	Lysing Matrix B	100 x 2 ml Tubes
6911-500	Lysing Matrix B	500 x 2 ml Tubes
6912-050	Lysing Matrix C	50 x 2 ml Tubes
6912-100	Lysing Matrix C	100 x 2 ml Tubes
6912-500	Lysing Matrix C	500 x 2 ml Tubes
6913-050	Lysing Matrix D	50 x 2 ml Tubes
6913-100	Lysing Matrix D	100 x 2 ml Tubes
6913-500	Lysing Matrix D	500 x 2 ml Tubes
6914-050	Lysing Matrix E	50 x 2 ml Tubes
6914-100	Lysing Matrix E	100 x 2 ml Tubes
6914-500	Lysing Matrix E	500 x 2 ml Tubes
6750-200	BioPulverizer™ System I (10 tubes each of 6910-6914)	50 x 2 ml Tubes

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

50ml Lysing Matrix Tubes

15ml Lysing Matrix Tubes


Cat. #	Description	Pack size
6950-010	BigA - Lysing Matrix Tubes	10 x 50ml Tubes
6950-050	BigA - Lysing Matrix Tubes	50 x 50ml Tubes
6951-010	BigB - Lysing Matrix Tubes	10 x 50ml Tubes
6951-050	BigB - Lysing Matrix Tubes	50 x 50ml Tubes
6953-010	BigD - Lysing Matrix Tubes	10 x 50ml Tubes
6953-050	BigD - Lysing Matrix Tubes	50 x 50ml Tubes
6960-010	BigClean - Lysing Matrix Tubes with stainless steel beads	10 x 50ml Tubes
6960-050	BigClean - Lysing Matrix Tubes with stainless steel beads	50 x 50ml Tubes
6954-010	Big E-Lysing Matrix Tubes	10 x 50ml Tubes
6954-050	Big E-Lysing Matrix Tubes	50 x 50ml Tubes

Cat. #	Description	Pack size
6930-050	TeenA - Lysing Matrix Tubes	50 x 15ml Tubes
6931-050	TeenB - Lysing Matrix Tubes	50 x 15ml Tubes
6932-050	TeenC - Lysing Matrix Tubes	50 x 15ml Tubes
6933-050	TeenD - Lysing Matrix Tubes	50 x 15ml Tubes
6934-050	TeenE - Lysing Matrix Tubes	50 x 15ml Tubes

RNA Extraction From Basil Leaves

Line 1. Marker
 Line 2. 1x40sec. with FastPrep® FP120 Instrument
 Line 3. 2x40sec. with FastPrep® FP120 Instrument
 Line 4. 3x40sec. with FastPrep® FP120 Instrument
 Line 5. 1x40sec. with FastPrep® 24 Instrument
 Line 6. 2x40sec. with FastPrep® 24 Instrument
 Line 7. 3x40sec. with FastPrep® 24 Instrument
 Line 8. 1x40sec. with FastPrep® FP120 Instrument
 Line 9. 1x40sec. with FastPrep® 24 Instrument

50mg of basil leaves were homogenized with the FastPrep®-24 or FastPrep® FP120 Instrument. One, two or three runs of 40 seconds were performed at a speed setting of 6.0 m/s and RNA was extracted with the FastRNA® Pro Green Kit.

Ready-to-use Protocols for DNA, RNA and Protein Isolation from Any Sample!

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FastRNA® Pro Soil-Direct Kit FastRNA® Pro Soil-Indirect Kit

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
DNA from various soil samples extracted with the FastDNA® SPIN Kit for Soil. 20% of the DNA isolated from 500mg soil was loaded on a 1.2% agarose gel (0.5X TAE). Soil was taken from: Lane 1: tomato pot; Lane 2 : sludge; Lane 3 : sandy soil; Lane 4 : under pine tree; Lane 5 : under palm tree; Lane 6 : green garden; Lane 7 : Nile Lilly pot; Lane 8 : lawn grass; Lane 9 : citrus tree; Lane 10 : avocado tree. DNA ranges from 4-20 kb.

References

1. Selesi D. et al (2005). Appl. Envir. Microbiol. Vol 71 : 175-184
2. Alexandrino M. et al (2004). Water Research. Vol 38 : 1340 - 1346
3. Mumy K.L. et al (2004). J. of Microbiological Methods. Vol 57 : 259 – 268


Cat #	Designation	Pack Size
6560-200	FastDNA® SPIN Kit for Soil	50 Preps

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Rat total RNA extracted with the FastRNA® Pro Green Kit. Approximately 2% of the total RNA isolated from 100 mg frozen tissue was loaded on to a 1.2% denaturing agarose gel (1X MOPS).

Cat #	Designation	Pack Size
6025-050	FastRNA® Pro Blue Kit (bacteria)	50 Preps
6035-050	FastRNA® Pro Red Kit (yeast and fungi)	50 Preps
6045-050	FastRNA® Pro Green Kit (plants and animals)	50 Preps

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
Figure 1
12% SDS PAGE of lysate of BL21 cells expressing the GST protein resulting from homogenization with the FastProtein™ Blue Matrix


Figure 2
12% SDS PAGE of lysate of yeast cells resulting from homogenization with FastProtein™ Red.

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Cat. #	Description	Size
6550-400	FastProtein™ Blue Matrix	50 x 2ml
6550-500	FastProtein™ Blue Matrix	100 x 2ml
6550-600	FastProtein™Red Matrix	50 x 2ml
6550-700	FastProtein™Red Matrix	100 x 2ml

Table of Typical FastPrep® Systems Settings for Optimal Sample Lysis

Sample Name	Sample Type	Quantity	Lysing Matrix	FastPrep® speed	FastPrep® time
HUMAN AND ANIMAL					
Human	Lung	50 mg	Lysing Matrix D	6.0	4x 30 sec.
Human	Breast	80 mg	Lysing Matrix D	6.0	2x 30 sec.
Human	Kidney	50 mg	Lysing Matrix D	6.0	40 sec.
Human	Thyroid Tumors	100 mg	Lysing Matrix A	6.0	3x 30 sec.
Mouse	Eye	10 mg	Lysing Matrix D	6.0	4x 30 sec.
Mouse	Heart	70 mg	Lysing Matrix D	6.0	4x 30 sec.
Mouse	Kidney	50 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Femur	40 mg	Lysing Matrix A	6.0	4x 30 sec.
Mouse	Leg Muscle	50 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Intestine	50 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Ear	45 mg	Lysing Matrix D	6.0	4x 30 sec.
Mouse	Tail	100 mg	Lysing Matrix A	6.0	4x 30 sec.
Mouse	Spleen	70 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Lung	50 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Liver	50 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Brain	50 mg	Lysing Matrix D	6.0	40 sec.
Mouse	Pancreatic cells (bHC9)	10 ⁷ cells	Lysing Matrix D	6.0	40 sec.

PLANT					
Alpowa Wheat	Leaf Tissue	75 mg	Lysing Matrix D	6.0	40 sec.
Alpowa Wheat	Seed	100 mg	Lysing Matrix A	6.0	40 sec.
Arabidopsis thaliana	Fresh Leaves	50 mg	Lysing Matrix D	6.0	40 sec.
Arabidopsis thaliana	Fresh Leaves	200 mg	Lysing Matrix D	6.0	2x 40 sec.
Bartlett Pear	Leaf Tissue	50 mg	Lysing Matrix D	6.0	40 sec.
Classic Oat	Leaf Tissue	75 mg	Lysing Matrix D	6.0	40 sec.
Classic Oat	Seed	100 mg	Lysing Matrix A	6.0	40 sec.
Corn	Leaf Tissue	100 mg	Lysing Matrix D	6.0	40 sec.
Crest Barley	Leaf Tissue	100 mg	Lysing Matrix D	6.0	40 sec.
Crest Barley	Root	300 mg	Lysing Matrix A	6.0	40 sec.
Kaybonnet Rice	Leaf Tissue	100 mg	Lysing Matrix D	6.0	40 sec.
Kaybonnet Rice	Seed	100 mg	Lysing Matrix A	6.0	40 sec.
Klages Barley	Root	300 mg	Lysing Matrix A	6.0	40 sec.
Klages Barley 70 mg Leaf Tissue 6.0 40 seconds	Leaf Tissue	70 mg	Lysing Matrix D	6.0	40 sec.
Tobacco	Leaf Tissue	75 mg	Lysing Matrix D	6.0	40 sec.

Table of Typical FastPrep® Systems Settings (continued)

Sample Name	Sample Type	Quantity	Lysing Matrix	FastPrep® speed	FastPrep® time
PLANT					
Lafitte Rice	Leaf Tissue	75 mg	Lysing Matrix D	6.0	40 sec.
Lafitte Rice	Sprout Leaf	100 mg	Lysing Matrix D	6.0	2x 30 sec.
Soybean	Seed	100 mg	Lysing Matrix A	6.0	40 sec.
Corn	Seed	100 mg	Lysing Matrix A	6.0	40 sec.
Oat FL 502	Leaf Tissue	75 mg	Lysing Matrix D	6.0	40 sec.
Oat FL 502	Seed	100 mg	Lysing Matrix A	6.0	40 sec.
Riser Oat	Leaf Tissue	70 mg	Lysing Matrix D	6.0	40 sec.
Richland Soybean	Leaf Tissue	100 mg	Lysing Matrix D	6.0	40 sec.
Tam Wheat	Leaf Tissue	75 mg	Lysing Matrix D	6.0	40 sec.
Tam Wheat	Root	80 mg	Lysing Matrix A	6.0	40 sec.
Tomato, Early Girl	Leaf Tissue	75 mg	Lysing Matrix D	6.0	4 x 30 sec.
Williams 82 Soybean	Leaf Tissue	70 mg	Lysing Matrix D	6.0	40 sec.
Wrens Rye	Leaf Tissue	100 mg	Lysing Matrix D	6.0	40 sec.
Pine	Needle	100 mg	Lysing Matrix A	6.0	30 sec.

BACTERIA					
Listeria monocytogenes	Cells	10 ⁹ cells	Lysing Matrix B	6.0	3x 30 sec.
Streptococcus pyogenes	Cells	10 ⁹ cells	Lysing Matrix B	6.0	20 sec.
Streptococcus mutans	Cells	10 ⁹ cells	Lysing Matrix B	6.0	30 sec.
Staphylococcus aureus	Cells	10 ⁸ cells	Lysing Matrix B	6.0	2x 40 sec.
Photobacterium luminescens	Cells	10 ⁹ cells	Lysing Matrix B	6.0	2x 30 sec.
Escherichia coli	Cells	10 ⁸ cells	Lysing Matrix B	6.0	30 sec.
Mycobacterium tuberculosis	Cells	10 ⁸ cells	Lysing Matrix B	6.0	2x 45 sec.
Lactococcus lactis	Cells	10 ⁸ cells	Lysing Matrix B	6.0	3x 30 sec.

YEAST AND FUNGI					
Saccharomyces cerevisiae	Cells	2x 10 ⁸ cells	Lysing Matrix C	6.0	40 sec.
Schizosaccharomyces pombe	Cells	10 ⁸ cells	Lysing Matrix C	5.0	4x 15 sec.
Candida albicans	Cells	10 ⁸ cells	Lysing Matrix C	6.0	2x 30 sec.
Cryptococcus neoformans	Cells	10 ⁸ cells	Lysing Matrix C	6.0	4x 30 sec.
Aspergillus fumigatus	Cells	10 ⁸ cells	Lysing Matrix C	6.0	2x 30 sec.
Fusarium solani	Cells	10 ⁸ cells	Lysing Matrix C	6.0	2x 30 sec.

ENVIRONMENTAL SAMPLES					
Sediments	Soil/rocks	50 mg	Lysing Matrix E	5.5	2x 30 sec.
Soil	Sandy sample	50 mg	Lysing Matrix E	4.0	4x 30 sec.
Soil	Litter	50 mg	Lysing Matrix E	5.5	30 sec.
Feces	Turd	300 mg	Lysing Matrix E	6.0	40 sec.

FastPrep®-24 Applications Questions and Answers

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Coloured scanning electron micrograph (SEM) of *Bacillus anthracis* spores

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
Quartz filter folded into the 15 ml tube in a TeenPrep adapter

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Quartz filter lysate after FastPrep lysing with no additional lysing matrix, for 30 seconds at 6.5 m/s speed settings

Yes!

- Please send me more information on FastPrep!
 I am ready for a demonstration, call me (_____) _____

Simply fill out and fax back, or submit online, and RECEIVE A FREE GIFT

North America: 1.800.334.6999
Europe: 00800.666.888
Japan: 03.3808.2401
Singapore: 65.6775.4536
Australia: 02.9838.7390

Fax:

First Name: _____

Last Name: _____

Phone Number: _____

Fax Number: _____

E-Mail Address: _____

Your email address is never shared and always maintained in confidence

I wish to receive email updates: **YES!** no

Company/Institution Name: _____

Internal Mail Code (if applicable): _____

Street Address: _____

Division/Department: _____

Building/Room: _____

City: _____ State/Province: _____

Zip / Postal Code: _____ Country: _____

I work for:

- Private Industry University/College
 Government Hospital/Medical Center
 Independent Analytical Lab Law Enforcement Agency
 Non for Profit Research Lab/Institute/Foundation
 Other (explain) _____

My primary field of work:

- Molecular Biology Biotechnology
 Pharmaceuticals Pathology
 Clinical Chemistry Environmental
 Medical/Biological/Clinical Forensics/Narcotics
 Plastics/Polymers/Rubber Agricultural/Foods/Beverage
 Energy/Petroleum Food and Beverages
 Inorganic Chemicals Organics Chemicals
 Other (explain) _____

Organization Type:

- University Scientific Research Institute
 Research Hospital Hospital
 Biotechnology Pharmaceutical
 Drugs Discovery Biopharmaceutical
 Chemical Manufacturer In Vitro Diagnostics Manufacturer
 Contract Research Organization Government Research
 Government Enforcement Agency Armed Services
 Group Purchasing Integrated Health Network
 Physician Office Lab Reference Lab
 Other (explain) _____

In addition to a FastPrep, I would like to receive information about the following products (please check all that apply):

- 2008/2009 Master Catalog Growth Media, Sera, Diets
 Molecular Biology Products Biochemicals
 Adsorbents RIA/EIA/ELISA assays
 In Vitro Diagnostic Products Radiochemicals
 Rare & Fine Chemicals Solvents
 Rapid Tests Alzheimer Research Kit
 Antibodies & other Immunobiologicals PCR Enzymes and Reagents

FAX or Submit online at www.mpbio.com/fastprep

Select One Free Gift:

Hurry! (Quantities are limited,
1 Gift per lab, Gifts may be substituted)

Clock

Travel mug

Laser Pointer

Mouse Pad

www.mpbio.com

Cat #02RI00254_FP24_0408_20M

MP Biomedicals North America, Tel: 1.800.854.0530 • MP Biomedicals Europe, Tel: 00800 7777 9999

MP Biomedicals Japan, Tel: 03.3808.2102 • MP Biomedicals Singapore, Tel: 65.6775.0008 • MP Biomedicals Australia, Tel: 02.9838.7422