

迷你数显折射计

ATAGO爱宕新的迷你PAL系列数显折射计 主要性能特点:

- 超凡的防尘防水性能 (IP65), 具备可完全流水冲洗的性能
- 革命性的E.L.I (外部光线阻止) 功能, 保证仪器在室外测量时的精度
- 清水归零, 操作方便
- 小巧轻便, 仅重100g, 尺寸5.5×3.1×10.9cm
- 具自动温度补偿功能 (ATC)
- 设计简洁, 可单手操作
- 配便携式挂钩, 便于携带



LAA699



LAA700



西域订货号	制造商型号	测量参数	测量范围	分辨率	精度	销售价
LAA699	PAL-1	糖度	0.0 ~ 53.0%	0.1%	± 0.2%	
LAA700	PAL-03S	融雪剂 (Nacl)	0.0 ~ 28.0%	0.1%	± 0.2%	
LAA701	PAL-10S	尿比重	1.000 ~ 1.060	0.001	± 0.001	
LAA702	PAL-22S	蜂蜜水分	12.0 ~ 30.0%	0.1%	± 0.2%	
LAA703	PAL-27S	豆浆浓度	0.0 ~ 20.0%	0.1%	± 0.2%	
LAA704	PAL-34S	乙醇浓度	0.0 ~ 45.0%	0.5%	± 1.0%	
LAA705	PAL-36S	甲醇浓度	0.0 ~ 35.0%	1.00%	± 2%	
LAA706	PAL-37S	异丙醇浓度	0.0 ~ 60.0%	0.2%	± 0.6% (0 ~ 30%) ~ ± 1.6% (60%)	
LAA707	PAL-38S	二甲基酰胺浓度	0.0 ~ 45.0%	0.2%	± 0.4%	
LAA708	PAL-39S	双氧水浓度	0.0 ~ 50.0%	0.2%	± 0.6%	
LAA709	PAL-40S	NaOH浓度	0.0 ~ 38.0%	0.1%	± 0.2%	
LAA710	PAL-41S	融雪剂 (Nacl)	0.0 ~ 30.0%	0.1%	± 0.2%	
LAA711	PAL-43S	Nacl	0.0 ~ 30.0%	0.1%	± 0.2%	
LAA712	PAL-79S	葡萄酒 T.A. 1990	0.0 ~ 26.0%	0.1%	± 0.2%	
LAA713	PAL-80S	葡萄酒 T.A. 1971	0.0 ~ 25.0%	0.1%	± 0.2%	
LAA714	PAL-83S	葡萄酒 KMW (or babo)	0.0 ~ 40.0°	0.1°	± 0.2°	
LAA715	PAL-84S	葡萄酒 baume	0.0 ~ 21.0°	0.1°	± 0.2°	
LAA716	PAL-85S	聚乙烯醇	0.0 ~ 20.0%	0.1%	± 0.2%	
LAA717	PAL-86S	葡萄酒	Oe (Ger): 0 ~ 240° Brix: 0.0 ~ 53.0%	1° 0.1%	± 1° ± 0.2%	
LAA718	PAL-87S	葡萄酒	Oe (Oechsle): 0 ~ 240° Brix: 0.0 ~ 53.0%	1° 0.1%	± 1° ± 0.2%	
LAA719	PAL-88S	丙二醇	丙二醇浓度 (V/V): 0.0 ~ 90.0% 丙二醇冷冻温度: 0 ~ -50 °C	0.2% (V/V) 1 °C	± 0.4% ± 1 °C	
LAA720	PAL-91S	乙二醇	乙二醇浓度 (V/V) 0.0 ~ 90.0% 乙二醇冷冻温度 0 ~ -50 °C	0.2% (V/V) 1 °C	± 0.4% ± 1 °C	
LAA721	PAL-96S	面汤浓度	0.0 ~ 20.0%	0.1%	± 0.2%	
LAA722	PAL-102S	切削油浓度	0.0 ~ 70.0%	0.1%	± 0.2%	
LAA723	PAL-Pâtissier	糕点糖度	Brix: 0.0 ~ 85.0 % Baume: 0 ~ 45°	0.1 % 1°	± 0.2% ± 1°	
LAA724	PAL-Plato	发酵麦芽汁浓度	0.0 ~ 30.0°P	0.1°P	± 0.2°P	

实验室检测分析仪器

数字式折射计

- Palette系列结构紧凑, 经多年使用验证, 功能强大, 易于操作
- 体积小, 精度高
- 方便的温度显示
- 具自动温度补偿功能
- 具外部光线阻止功能



LAA736



LAA737



LAA743



西域订货号	制造商型号	品牌	测量参数	测量范围	分辨率	精度	销售价
LAA736	PR-32 α	ATAGO	糖度	0.0 ~ 32.0%	0.1%	± 0.1%	
LAA737	PR-60PA	ATAGO	异丙醇浓度	0.0 ~ 60.0% (W/W)	0.1%	± 0.5% (0 ~ 30%)	
LAA738	PR-50HO	ATAGO	双氧水浓度	0.0 ~ 50.0% (W/W)	0.1%	± 0.5%	
LAA739	PET-109	ATAGO	乙醇浓度	0.0 ~ 45.0% (W/W)	0.1%	± 0.5% (0 ~ 30%)	
LAA740	PR-100SA	ATAGO	海水盐度计	0.0 ~ 100.0‰	1‰	± 2‰	
LAA741	PR-40DMF	ATAGO	二甲基酰胺浓度	0.0 ~ 40.0% (W/W)	0.1%	± 0.3%	
LAA742	UG- α	ATAGO	尿比重	1.0000 ~ 1.0600	0.000	± 0.0010	
LAA743	WM-7	ATAGO	葡萄酒	Brix: 0.0 ~ 45.0% T.A. 1990: 0.0 ~ 26.0% T.A. 1971: 0.0 ~ 25%	Brix: 0.1% T.A. 1990: 0.1% T.A. 1971: 0.1%	Brix: ± 0.1% T.A. 1990: ± 0.2% T.A. 1971: ± 0.2%	
LAA749	DBR35	-	糖度	0 ~ 35%	0.2%	± 0.1%	
LAA750	DBR45	-	糖度	0 ~ 45%	0.2%	± 0.1%	
LAA751	DSA1	-	盐度	0 ~ 28%	0.2%	± 0.1%	