

Leica EM AMW

Automatic Microwave Tissue Processor for Electron Microscopy



Bringing Speed and Automation to

With more than 130 years experience in the manufacturing of scientific instruments, Leica is proud to introduce the worlds first Automatic Microwave Tissue Processor for electron microscopy specimen preparation. The unique Leica EM AMW has been developed in close cooperation with the scientific community in order to meet the requirement's of today's modern EM facility:

- Reducing sample turn around time
- High level of automation with minimal user interface
- Accurate reagent temperature measurement & control
- Ease of set up and ease of use
- Reproducible results
- Safety

With the Leica EM AMW, tissue can be rapidly processed, embedded and polymerised into resin for subsequent analysis with unmatched automation, in hours rather than days. Automated microwave tissue processing for research, routine & diagnostic EM applications are possible. At the nucleus of the new Leica EM AMW resides the technology which makes it so unique

The Mono-mode microwave chamber design directs the microwave energy into a defined area, resulting in a homogenous field pattern surrounding the sample. Guiding the microwave energy in this manner enables the system to provide a uniform effect over the samples & reagents. 100 % of the microwave energy is absorbed by the processing reagent and sample without hot and cold spots and negating the need for water loads.

The other major advantage of the Mono-mode design is the energy density. Even with a low power setting, the energy density is actually much greater compared to a standard laboratory microwave tissue processor.

In addition, the user has the freedom to choose between continuous or pulsed power settings for various applications, such as pulsing for fixation steps and continuous power for dehydration. The software and non-contact infra-red temperature sensor ensures accurate processing parameters are maintained with maximum safety for user and specimens.



EM Tissue Processing

Load Samples ...

20 polypropylene vials of 20 ml are available to the user for one processing run. Reagents can be pre-loaded in a fume cupboard into the vials already placed onto the carousel. Vials are automatically closed via the sealing lid to avoid escape of toxic vapour when transporting the carousel to the instrument and during processing. The carousel with pre filled vials is slid into place and the Leica EM AMW is ready to start.

Up to 40 specimens are supported in baskets and stacked onto a stem. The microwave lid holds the basket assembly in place while processing.

According to the programme in use, the first vial is moved automatically into the microwave chamber where the specimen stack is positioned. A fully automatic run can then be performed.



Specimens ready to be inserted into the microwave chamber



Detail view, carousel in working position, sealing lid open and vial on the way into the MW chamber



Mouse liver (immersion fixed). Hepatocytes with or- Mouse kidney, cortex (immersion fixed); Parts of a ganelle-rich cytoplasm and round cell nuclei contain- proximal and of a distal renal tubule with blood capiling prominent nucleoli; in between a blood sinusoid laries (filled with erythrocytes) in between. with preserved space of Disse and bile canaliculi with numerous microvilli.

Mouse kidney, detail of a cell of the distal renal tubule; the cytoplasm displays convolutions of the cell membrane and numerous mitochondria.

Courtesy of Dr. Josef A. Schroeder, Uniklinikum Regensburg

... and Press START to Process

The colour screen and intuitive software makes it the easiest-to-use microwave tissue processor on the market, even for those who have little or no experience with microwave technology in the laboratory. However, the software also offers unmatched flexibility, allowing microwave parameters to be adjusted according to your requirements. The design of the user interface provides all necessary programme information at a glance.

All necessary process information is available with the click of a mouse.

- Vial in process
- Reagent details
- Temperature & time scale graphic
- Actual temperature progression
- Time to end



Schematic of the instrument status Graphic overview of the selected programme with all vials and their set temperatures.

Leica MICROSYSTEMS		time to end 02:43:09 h
	processing	temperature 37°C
Program 1	°C	
vial 1		
Karnovsky	60	
37°C	40	
00:05:00 h		min
Tissue Processin	g 0 1 2 3	4 5
Pause St	op Menu Loc	k 11:02:47 26.04.2007

Schematic of the processing run Showing the temperature profile for each step in use. Data transfer via USB-Stick possible.

... Press START to Polymerize ...

Polymerization of up to 20 specimens per run can also be performed under microwave irradiation within a short time. Samples are placed in re-useable polymerisation forms and stacked. After placing the stack onto the carousel a polymerisation programme can be started. If more samples need to be polymerized, a second or third stack can be placed onto the carousel for automatic polymerization.



	processing 56°C
Polymerisation	°C
vial 1	
Epon	60
80°C	40
02:54:50 h	
	20 20 00 400 450

Select a polymerisation programme from the memory and press START...

...hardened blocks are ready for ultrathin sectioning.

EM Tissue Processing

Specimen top loading system for easy access

The 6.5" mouse controlled colour screen

in combination with the intuitive user interface makes programming a matter of a mouse click. Storing of up to 99 programmes and unlimited temperature slope settings makes the instrument highly flexible.



The robotic reagent system

handles the 20 reagent vials sitting on the carousel and performs automatic reagent changes according to the programme. Up to 40 specimens located in special re-usable, flow through baskets, can be processed at once without any interaction from the user.

Leica Microsystems – the brand for outstanding products

Leica Microsystems' mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement and analysis of micro-structures.

Leica, the leading brand for microscopes and scientific instruments, developed from five brand names, all with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Yet Leica symbolizes innovation as well as tradition.

Leica Microsystems – an international company with a strong network of customer services

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Tel. +33 1 47 32 85 85	Fax +33 1 47 32 85 86
Germany:	Bensheim	Tel. +49 6251 136 0	Fax +49 6251 136 155
Italy:	Milan	Tel. +39 0257 486.1	Fax +39 0257 40 3475
Japan:	Tokyo	Tel. + 81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives of Leica Microsystems in more than 100 countries. The companies of the Leica Microsystems Group operate internationally in three business segments, where we rank with the market leaders.

• Microscopy Systems

Our expertise in microscopy is the basis for all our solutions for visualization, measurement and analysis of microstructures in life sciences and industry. With confocal laser technology and image analysis systems, we provide three-dimensional viewing facilities and offer new solutions for cytogenetics, pathology and materials sciences.

• Specimen Preparation

We provide comprehensive systems and services for clinical histo- and cytopathology applications, biomedical research and industrial quality assurance. Our product range includes instruments, systems and consumables for tissue infiltration and embedding, microtomes and cryostats as well as automated stainers and coverslippers.

• Medical Equipment

Innovative technologies in our surgical microscopes offer new therapeutic approaches in microsurgery.

