



## 球磨机

球磨机是利用滚筒内的金属球在伴随滚筒的滚动的过程中对样品进行冲击破碎的设备, 适合涂料, 矿业及陶瓷, 药物原料, 食品等行业的破碎需求。



高强度耐冲击金属滚筒, 内腔为技术材料, 外部经特殊处理, 适合滚动驱动, 保证稳定有效的破碎效果, 有不锈钢和铝两种滚筒可以选择, 也可以根据客户要求定制  
数字显示与控制, 滚筒滚动 50 to 300 RPM 可调, 采用反馈控制技术, 设备将根据滚筒内样品的装载量自动调整输出功率, 保证滚筒在设定转速下连续稳定运行

数字定时功能: 1min-99hr, 或者连续运行

滚筒最大长度可以到 700mm, 有五种不同体积可以选择: 1x15L; 1x11L; 2x5L; 3x3L; 4x1.5L

安全设置及执行标准:

1. 球磨机工作腔带安全盖, 避免意外情况对操作人员的伤害
2. 如果在研磨工作过程中, 滚筒盖意外打开, 球磨机将自动停机
3. Community directives: 2004/108/EC 2006/42/EC, 2006/95/EC
4. Standards: EN 61010-1 EN 61326

设备外形尺寸: 125 × 49 × 34(W × D × H);

重量: 72Kg;

功率: 150W

## 订货信息

订货号	产品描述
ML 007	BALLS MILL 主机, 220 V / 50-60 Hz
<b>不锈钢滚筒及破碎滚珠</b>	
PI 062	Jar of 1 liter of stainless steel
PI 063	Jar of 3 liters of stainless steel
PI 064	Jar of 5 liters of stainless steel
PI 226	Jar of 15 liters of stainless steel
PI 061	Balls of Ø9 mm of stainless steel (1 Kg)
PI 058	Balls of Ø15 mm of stainless steel (1 Kg)
PI 059	Balls of Ø20 mm of stainless steel (1 Kg)
PI 060	Balls de Ø30 mm of stainless steel (1 Kg)
<b>铝制滚筒及破碎滚珠</b>	
PV 035	Jar of alumina of 1 liter
PV 036	Jar of alumina 3 liters
PV 037	Jar of alumina 5 liters
PV 040	Ball of alumina of 20 mm (1Kg)
PV 042	Ball of alumina of 30 mm (1Kg)