



PIPETTING JUST GOT BETTER

美杏高德**Personal Pipettor+**无人值守

集聚多种移液需求功能于一身的移液工作站。仅需单击运行即可完成需要多次格式变化操作、不同类型的实验。

Personal Pipettor+无人值守是具有直观操作界面、在线自动更换吸头功能的自动化工作站，能帮助您取得最佳的实验结果 - 更加精确、可靠和可重复。

Personal Pipettor+无人值守降低工作量（如：梯度稀释），消除手工操作误差，在有限的空间内得出满意的实验结果。

简单易用的高性能移液工作站

应用领域

板复制
板格式重组
试剂添加
化合物添加
逐列、逐排梯度稀释
纯化
PCR体系构建
Elisa
NGS文库构建
细胞实验
更多.....

功能

5个直线型板位，加1个清洗板位/废液槽
1, 8, 12, 16, 24, 96, 384 通道移液
多样化移液锁头
独立于移液系统的高精度移液锁头
快速在线更换的EZ-load吸头技术
直观的iPad控制界面，简单易用
同时支持电脑软件操作
可整合微孔板叠板机和机械手
真正实现自动化无人值守
可整合钢针清洗模块
体积小巧，最大化节省操作空间，
适合放置于多数标准通风橱

PERSONAL PIPETTOR 无人值守移液工作站

参数	
尺寸	宽711mm x 深483mm x 高607mm
重量	约45kg
电源	100-240VAC, 50/60Hz
工作温度	10-37°C

www.apricotdesigns.com

Apricot Designs, Inc. 677 Arrow Grand Circle, Covina, CA 91722 USA

Telephone: 626-966-3299 | Fax: 626-966-3200

Apricot Designs (Shanghai) Co., Ltd. 上海美杏高德生物科技有限公司
上海市浦东新区张江路665号德宏大厦1008室 邮编: 201203

电话: (86)21 50796338
传真: (86)21 50796228

© 2017 Apricot Designs, Inc. The Apricot Designs logo is a trademark of Apricot Designs, Inc.
Specifications subject to change. All information is deemed accurate at the time of publication.
Apricot Designs is not responsible for any errors that may occur.