



The Homogenizer Company™

# Omni Mixer Series Homogenizers

Homogenizing tough samples for over 50 years



# **Omni Mixer Series** Homogenizers

Omni Mixer and Omni Macro Homogenizers offer Available accessories make them flexible enough to process a wide variety and volume of samples.







- > Bulk vaccine products
- > Processing fruits & vegetables
- > Plant homogenization for environmental testing
- > Petroleum by-products & waste
- > Mixing anesthetics for research
- > Lotions & oil emulsions
- > Culinary use featured in Modernist Cuisine
- > Animal tissue
- > Paper pulp processing
- > Paint, ink and dye application

#### < Omni Mixer

600 watt homogenizer

**ORDER NO. 17105** 

Processing Range: 0.25 mL to 10 L

Generator Probes: 7, 10, 20, 32, 35, 45, 55 mm

Blade Assemblies: 1" and 2" in sealed chambers Power Rating: 600 watt

Speed: Variable from 500 to 20,000 rpm Dimensions (with stand): 24" H x 10" W

Weight: 9.5 kg (21 lbs.)

Warranty: 3 year motor warranty



### Omni Macro 1,800 watt homogenizer

**ORDER NO. 17505** 

Processing Range: 0.25 mL to 30 L

Generator Probes: 7, 10, 20, 32, 35, 45, 55 mm Blade Assemblies: 1" and 2" in sealed chambers

Power Rating: 1,800 watt

Speed: Variable from 1,000 to 20,000 rpm Dimensions (with stand): 24" H x 10" W

Weight: 9.5 kg (21 lbs.) Sound Level: <68 db

Warranty: 3 year motor warranty



#### **Omni Macro ES**

1,800 watt digital, programmable, homogenizer

**ORDER NO. MES-115** 

Processing Range: 0.03 mL to 50 L

Generator Probes: 7, 10, 20, 32, 35, 45, 55 mm Blade Assemblies: 1" and 2" in sealed chambers

Power Rating: 1,800 watt

Speed: Variable from 1,000 to 20,000 rpm, soft

start, 0.5% accuracy

Programs: 12 programs with 12 time and speed

steps, including ramping function

Interface: RS-232 and windows software

Dimensions: 34.5" max Weight: 14.5 kg (32 lbs.) Sound Level: <68 db

Warranty: 3 year motor warranty

Compliance: CE approved



The Homogenizer Company™

The versatile **Omni Mixer Line** offers a wide variety of accessories for all of your processing needs.



#### Micro-Homogenizer Attachment (MHA)

- Converts your OmniMixer,
   OmniMacro or Digital MacroES
   into a Micro-Sample Homogenizer
- Sample volume ranges from 0.5 to 5 mL with your choice of either a Stainless Steel or Nylon Insert Chamber



## **Aerosol Sealed Stainless Steel Chamber Assemblies**

- > Fluted to enhance flow
- Polished interior to facilitate cleaning
- > Available in 316 L stainless steel
- Compatible with generator probes and cutting blade assemblies
- Available in 50 mL, 200 mL and 400 mL sizes



### **Aerosol Sealed Glass Chamber Assemblies**

- > Rectangular shape to enhance flow
- Not recommended for substances that may cause breakage
- Compatible with generator probes and cutting blade assemblies
- Available in 50 mL, 200 mL, 473 mL, and 946 mL sizes



## Aerosol Sealed Plastic Chamber Assemblies

- Process sample in a centrifuge tube
   without need to transfer a second container
- Compatible with generator probes and cutting blade assemblies
- Available in 50 mL, 200 mL, 473 mL and 946 mL sizes





### Stainless Steel Generator Probes

- > Wide variety of diameters and lengths
- > Flat bottom or saw tooth models for different sample types
- > Available in 7mm, 10mm, 20 mm, 35 mm, 45 mm and 55 mm diameter sizes