



P5

PLASMA PEN

等离子清洗笔

The Plasma Pen is a hand-held atmospheric cold plasma device that requires no external gas supply; just plug it in and start cleaning! Ideal for spot organic contaminant cleaning, hydrophilic treatment and surface modification of small samples such as SEM specimens in research labs.

P5等离子笔是一种便携手持大气冷等离子发生器，不需要外部气体供应，只需将其插电即可开始清洗！特别适用于实验室样品表面的有机污染物等离子清洗、亲水化处理和表面改性。比如扫描电镜SEM样品观测前处理等用途。



High Efficiency!

高效!

Soft Plasma!

温和等离子体!

Plug and Play!

即插即用!

SuPro Instruments LTD

速普仪器

Nanshan ,


Shenzhen, China

[Http://www.suproinst.com](http://www.suproinst.com)

Tel: 86-755-26642901

Fax: 86-755-26419205

技术规格	参数
工作模式	手持式大气等离子体清洗笔
工作气体	空气/氩气气氛, 无需外部气源气瓶
等离子温度	< 50 摄氏度
处理宽度	5-10 mm
工作距离	5-10 mm
清洗速度	5-20 mm/s
清洗功率	5 W
操作方式	手持按钮控制/ 即插即用
使用用途	材料表面活化、亲水处理、有机污染物清洁、等离子杀毒等用途
使用环境	严禁在易燃易爆气体环境中使用 严禁在易燃易爆粉尘环境中使用 环境温度~10-40 摄氏度, 相对湿度 < 80%
重量	~ 0.25 kg
尺寸	200 mm 长* 45 mm 直径
电源	外接电源 5 V DC/ 2 A
功耗	< 10 W
冷却方式	风冷
质量保证	一年免费/终身维护
备注	以上所列技术规格与参数更新恕不另行通知, 如有疑问请联系我们



深圳市速普仪器有限公司

地址: 深圳市南山区桂庙路22号向南瑞峰B2009

电话: 0755-26642901 传真: 0755-26419205

邮件: sales@suproinst.com

Http://www.suproinst.com

深圳市速普仪器有限公司 (北京办事处)

地址: 北京市海淀区太平路甲40号金玉元E座105室

电话: 010-63839091 传真: 010-63839091

邮件: zkchang@suproinst.com

淘宝店: 速普仪器