



www.labstac.com

AIR-JACKETED MULTI GAS INCUBATOR

Designed for a broad range of applications in the biomedical, pharmaceutical, medical and clinical research. Airjacketed Multi gas Incubator system help to facilitate a stable cell culture environment. Easy to clean and minimizes recovery time and the risk of contamination.

Used in Biomedical, Pharmaceutical, Medical, Clinical research, Laboratory.

Also known as Laboratory Air-jacketed Multi gas Incubator, Multi gas Incubator Air-jacketed.

IO31 AIR-JACKETED MULTI GAS INCUBATOR

151 liters internal volume
Zironuim oxide Sensor
IR CO2 sensor
Water Level Alarm
Integrated electric siphon pump
Auto-zero
Temperature deviation (at 37° C)
Moist-heat disinfection at 90°C
Relative humidity ≥95%



SPECIFICATIONS

Model	IO311	IO312	
Oxigen	High	Low	
Capacity	151 L		
Temperature Range	Amb. +2 to 55°C		
Temperature Accuracy	<±0.1°C		
Temperature recovery Time	≤7 mins(After 30 sec. Door opening)		
CO2 Range	0~20%		
CO2 Accuracy	±0.1		
CO2 Sensor	IR standard or TC optional		
O2 range(% CO2)	3%-20%, 22%-85%		
O2 accuracy(%CO2)	±0.2		
O2 sensor	zirconuim		
Humidity Range	≥95%±3%		
Interior	Stainless steel, type 304		
Temperature Sensors	PT1	PT1000	
Internal Dimention (WxDxH)	470×530	470×530×607 mm	
Overall Dimension (WxDxH)	637×768×869 mm		
Shelves/Trays	3 standard, 10 maximum		
Shelf Dimension(WxD)	423×445 mm		
Max. Load per Shelf	10 kg		
CO2 control system	Microprocessor PID		
Weight	80	80 kg	
Power	≤650 VA+10%		

www.labstac.com

2

Power Supply	220 V±10%/ 50 Hz
--------------	------------------

ACCESSORIES

Accessory Code	Name	Units
IX1007	Standard shelf	3 Sets
IX1026	CO2 steel cylinder auto-switch device	1 Set

OPTIONAL ACCESSORIES

Accessory Code	Name
IX1016	CO2 regulator
IX1018	Air connector
IX1027	Shelf
IX1028	Shelf support

Accessory Code	Name
IX1029	Small HEPA filter
IX1030	O2 regulator
IX1031	N2 regulator



LABSTAC LTD

Kemp House 152 – 160 City Road, London, EC1V 2NX, United Kingdom Email: contact@labstac.com | Website: www.labstac.com