

pH-meter OpH218


OrigaMeter PH 计



- ✓ **Sustainable and repairable**
- ✓ **Reliable and fast results**
- ✓ **Easy to use and simple connections**
- ✓ **Data transfers (Regressi, ExAo, Excel)**
- ✓ 可持续和可修复
- ✓ 可靠而快速的结果
- ✓ 易于使用且连接简单
- ✓ 数据传输 (Regressi, ExAo, Excel)

OrigaLys was founded by four R&D engineers from Radiometer company.


测量范围	pH: - 9 to 23 pH
	mV: \pm 2000 mV
	°C: - 10°C to 110°C
分辨率	pH: \pm 0,01 pH
	mV: \pm 0,1 mV
	°C: \pm 0,1°C
输入电极阻抗	$> 2 \times 10^{12}$ Ohms
电极极化	可以在BNC输入的连接电极上施加10A的电流
稳定性	3 mV/min (\approx 0,05 pH/min)
测量	2 种模式 : <ul style="list-style-type: none"> • 连续: 连续显示pH/mV和温度 • 自动: 当漂移低于稳定性标准时, 将记录并存储pH值结果.
校准	1 to 3 点
缓冲区选择	3种模式: <ul style="list-style-type: none"> • 自动识别 (Series: IUPAC or 4-7-10) • 参照手册选择范围 (IUPAC and 4-7-10) • 手动调节pH 值
校准范围	<ul style="list-style-type: none"> • 斜率: 95 to 102% • pH值: 5.80 to 7.50 pH
语言	法语, 英语, 西班牙语, 德语和意大利语
显示	128x64, OLED, 尺寸 60 x 30 mm
输入 / 输出	<ul style="list-style-type: none"> • 1 个玻璃电极或组合电极输入(BNC 接口) • 1 个参比电极输入(TAKE BANANE 4 mm) • 1 个温度传感器输入(take RCA / CINCH) • 1 个RS232 接口(SUB-D 9针口) • 1 个USB 2.0 接口(type B socket) • 1 模拟输出(take mini-DIN8)
外壳	坚固及防污(INOX - PC - PMMA)
尺寸 (H x W x D)	80 x 140 x 180 mm
重量	1 Kg
电源	2 种: <ul style="list-style-type: none"> • AC 交流适配器 12Vdc, 1A, 12W • USB 2.0
使用环境	<ul style="list-style-type: none"> • 温度范围: 5 to 40°C • 湿度范围: 20 to 80%



参比电极

其它

型号	① OGR005 Type REF321	② OGR004 Type REF421	③ OGR006 Type XR300	④ OGR003 Type XR110	⑤ D110GL008 Type AL120
尺寸	∅ 8 x 103mm	∅ 8 x 103mm	∅ 8 x 120mm	∅ 8 x 120mm	∅ 8 x 140mm
材质	Glass	Glass	Glass	Glass	Glass
参比系统	Ag/AgCl	Calomel	Ag/AgCl	Calomel	-
电解液	KCl 3M with saturated AgCl	Saturated KCl	KCl 3M with saturated AgCl	Saturated KCl	-



组合式 pH 电极

Non-combined pH electrodes

型号	① OGP201 Type pH2401-8	② OGP202 Type pH3001	③ OGP203 Type pH3005	④ OGP001 Type pHG301	⑤ OGP002 Type pHG311
pH 范围	0 - 12	0 - 12	0 - 12	0 - 12	0 - 14
温度范围	-5 to 80°C	-5 to 80°C	-5 to 80°C	-5 to 80°C	-5 to 80°C
尺寸	∅ 12 x 103mm	∅ 12 x 103mm	∅ 8 x 103mm	∅ 8 x 103mm	∅ 8 x 103mm
材质	Glass	Glass	Epoxy	Glass	Glass
参比系统	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl	Ag/AgCl
电解液	KCl 3M with saturated AgCl	KCl 3M with saturated AgCl	KCl 3M with saturated AgCl	Saturated KCl	Saturated KCl

See the full list of electrodes e-mail to
jeffzhou@chohao.com

pH-meter Packs

独立甘汞版

pH-meter OpH218



- ✓ 两根电缆
- ✓ 一个pH电极
- ✓ 一枚甘汞参比电极

独立 Ag/AgCl 版

pH-meter OpH218



- ✓ 两根电缆
- ✓ 一个pH电极
- ✓ 一个 Ag/AgCl 参比电极

组合式

pH-meter OpH218



- ✓ 一根电缆
- ✓ 一个组合式pH电极
- ✓ (玻璃或环氧树脂)

单机版

pH-meter OpH218



Susceptible
d'être modifié:
17-04-2020

中昊清远(北京)科技有限公司

北京市朝阳区延静里中街3号院 3号楼309室
☎ (010) 84786530

☎ (010) 84786530

E-mail: jeffzhou@chohao.com

www.chohao.com