

## Comparison vs. TM3000



TM3030  
External view




		TM3030	TM3000
Acc. Voltage	5kV	3 Stage (Low mag., Mid mag., High mag.)	1 Stage
	15kV	1 Stage	1 Stage
	15kV-HC	1 Stage	1 Stage
Magnification		15x to 30,000x	15x to 30,000x
Specimen Stage		Manual Motor Driven	Manual Motor Driven
Applicable OS		Windows7 (32/64bit)	Windows7 (32bit)

## @5kV surface detailed image quality enhanced

1. Thanks to improvement of an optical system, 5kV observation is further enhanced throughout high magnifications.
2. It is very effective function when low vacuum observation or for ion milled specimens.



TM3000



TM3030 (High mag. mode)

Specimen : Magnetic head  
Mag. : 10,000 x

Low mag., Mid. mag., & High mag. are selectable according to specimens or observation magnifications at 5kV mode.