

Beijing Byrd Hengyue Science & Trading Co., Ltd.

Room 13 west building No.1 Yard Saoziying
Haidian Dist. Beijing China(100091)
Tel:86-010-62883589 62869761
Fax:86-010-62871930
E-mail:great77dream@yahoo.com.cn

Densitometer KDM01

tap-density measurement according to DIN ISO 3953

Application / Features

tap-density measurement acc. to DIN ISO 3953 :

- in particular of ductile metal-flakes
- other flake-like powder and mixtures
- equiaxial powder and powder in general
- automatic operation with 3 tap-frequencies
- standard-laboratory test-cylinders (25-100ml)
- reproducibility, quality-control
- easy cleaning and maintenance

Technical Data, Dimensions



L x B x H (without testcyl.)	390 x 380 x 270
net weight	17 kg
test-cylinder	25/100 ml
tap-frequency	100 / 200 / 300 Hz
shift	3 mm
drive power	55 W
electrical control	24V-110/220VAC

Options

- alternative test-cylinders and adapters
- fine-scale 500 g



Densitometer KDM01 with test-cylinders 25 & 100ml



北京博德恒悦科贸有限公司

地 址: 北京市海淀区骚子营 1 号院西楼 13 室, 100091
电 话: 010-62883589 62869761 传 真: 010-62871930
联系人: 孟德强 邮 箱: great77dream@yahoo.com.cn



www.zoz.de - info@zoz.de