

Simethicone Shaker



Model 85-AA

The Simethicone Wrist Action® Shaker has been modified to meet the requirements for Simethicone testing as described in the current edition of the U.S. Pharmacopoeia. The test employs a radius of 13.3 (+/-) 0.4 cm, measured from the center of the shaft and shakes through an arc of 10 degrees (lever arm set at approximately 5) at a frequency of 300 (+/- 10%) strokes per minute. Burrell Scientific redesigned the shaking mechanism and supplies the Multi-SUREGRIP™ Clamps and Shafts which allow the Shaker to achieve the required performance levels for the Simethicone Foaming Solution Test. The Shaker meets these professional standards

Specifications:

- 220V, 50CYC, 300 Osc/min
- 2 clamps, 4 vessels held
- Max load tested - 10LBS
- Weight: 68.00lbs
- Dimensions (HxDxW): 20 x 19 x18)

ACCESSORIES

1. Multi SUREGRIP™ Clamp 8 Place (PAIR)



2. Simethicone Bottle 8 oz

