

# OMNI TH<sub>Q</sub> DIGITAL HANDHELD TISSUE HOMOGENIZER

The **OMNI TH<sub>Q</sub>** with Whisper Drive Technology™ is a powerful tissue homogenizer that can handle a wide range of samples from <200µL to 120mL.

It provides reliable, repeatable performance with an ultra-quiet brushless motor. An integrated digital speed control stores the last set speed for ultimate repeatability. Tissue samples are usually processed in less than 30 seconds.

The **Omni TH<sub>Q</sub>** operates with patented **OmniTip™ Plastic Generator Probes**, available in both soft tissue and hard tissue versions. Omni Tips are ideal for customers concerned about cross-contamination and carryover between samples. Stainless steel generator probes (5 mm to 10 mm diameter) are available for larger and tougher sample types. Specially designed hybrid steel/plastic probes are also available.

The **Omni TH<sub>Q</sub>** is an ideal fit for a wide range of applications from tissue homogenization and dissociation, PD/PK studies, tablet disruption, emulsions and suspensions, cell disruption, food samples and much more. Lightweight enough to be hand-held it may also be stand mounted for lengthier processing.

**Processing Range:** <200 µL to 120 mL

**Speed:** Variable from 8,000 to 32,000 rpm  
(increments of 500 rpm)

**Power Rating:** 40 watt max

**Dimensions:** 21.8cm (8.6 in) height,  
5.54 cm (2.18 in) max. diameter

**Weight:** 340 g (12.0 oz.)

**Sound Level:** <68 db

**Warranty:** 3 year motor warranty

**Standards Approval:** CE Approved

## Features

- Three digit speed readout accurate to +/- <5% at top speed
- Whisper Drive Technology™ provides ultra quiet operation
- On-board memory stores last set speed
- Bayonet probe mounting system for fast release
- Uses stainless steel probes or disposable Omni Tip™ plastic probes
- Motor stand mount included
- Light weight ergonomic design



# OMNI International Stainless Steel Generator Probes

OMNI's **Stainless Steel Generator Probes** come with PTFE upper and lower bearings. The entire probe can be autoclaved, or easily disassembled for cleaning. Optional Stainless steel upper bearings extend bearing life but are not autoclavable.



ORDER NO. B5-075

ORDER NO. B7-110ST

ORDER NO. B10-110ST



## Patented Omni Tip™ plastic disposable generator probes

Proprietary **Omni Tip™ Plastic Probes** are available in both soft tip (deaggregation) and hard tissue (frozen) versions. Engineered from durable plastics, Omni Tips™ can be cleaned and reused, yet are economical enough that they can be disposed of after each use to eliminate cross contamination.

*The Omni Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.*

### A selection of our best-selling plastic probes:



SOFT TISSUE PLASTIC GENERATOR PROBES		ORDER NO.
Omni Tip™ Plastic Generator Probes, Soft Tissue	25 pack	30750
Omni Tip™ Plastic Generator Probes, Soft Tissue	50 pack	32750
Omni Tip™ Plastic Generator Probes, Soft Tissue	100 pack	34750
HARD TISSUE PLASTIC GENERATOR PROBES		ORDER NO.
Omni Tip™ Plastic Generator Probes, Hard Tissue	25 pack	30750H
Omni Tip™ Plastic Generator Probes, Hard Tissue	50 pack	32750H
Omni Tip™ Plastic Generator Probes, Hard Tissue	100 pack	34750H

Omni TH<sub>Q</sub> package includes: TH<sub>Q</sub> motor, 3 year motor warranty, stand clamp, 6 Soft Tissue & 6 Hard Tissue Omni Tip™ Plastic Homogenizer Probes, tool kit, and instruction manual. **ORDER NO. 12-500**

**AORAN**  
**奥然科技**  
全国统一服务热线  
4009-670-570  
[www.aoran.cn](http://www.aoran.cn)

**奥然科技 上海**  
地址：上海市徐汇区漕溪路250号银海大厦A1012室  
电话：021-6451 6293 6451 6297 6451 6299

**奥然科技 北京**  
地址：北京市海淀区北三环西路甲30号双天大厦423室  
电话：010-6840 0176 6840 0177

**奥然科技 广州**  
地址：广州市天河区五山路141号尚德大厦A栋1707室  
电话：020-8551 6682 8551 6679

**奥然科技 南京**  
地址：南京市中央路417号先锋广场1036室  
电话：025-5850 1866 5850 1966