## JHD MICRO HOMOGENIZER

The **µHb Micro Homogenizer** is a hand-held battery powered rotor-stator homogenizer designed for use inside or outside of a laboratory setting. The µHb is powered by 12 volt lithium ion battery that provides operating speeds from 10,000 to 35,000 rpm for rapid homogenization of samples in volumes from <200 µL to 100 mL.

The  $\mu$ Hb is an ideal solution for field based homogenization procedures involving the processing of tissues, faeces or soft plant material. Compatible with stainless steel generator probes from 5 mm to 10 mm or disposable 7 mm Omni Tip<sup>TM</sup> plastic generator probes for customers concerned with cross contamination. Hybrid generator probes are also available for applications requiring the durability of stainless steel with the convenience and cleanliness of plastic probes.

## **Specifications**

Processing Range: <200 µL to 100 mL

Speed: Continuously variable from 10,000 to 35,000 rpm

Motor: High performance 12 Volt motor

Dimensions: 21.8cm (8.6 in) height, 5.54 cm (2.18 in) max. diameter

Warranty: 1 year motor warranty

ORDER NO. UHB-115 UHB-220

### **Features**

- Powerful cordless 12 Volt motor
- Light weight portable homogenizer
- 1/4 turn quick connect coupling
- Compatible with Omni Tip<sup>™</sup> plastic and stainless steel generator probes





## OMNI International Stainless Steel Generator Probes

#### OMNI's Stainless Steel Generator Probes

come with PTFE upper and lower bearings. The entire probe can be autoclaved, or easily disassembled for cleaning. Optional Stainless steel upper bearings extend bearing life but are not autoclavable.



The Omni Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.



<b>ORDER NO. B5-075</b>	
-------------------------	--

ORDER NO. B7-110ST

ORDER NO. B10-110ST

# Patented Omni Tip<sup>™</sup> plastic disposable generator probes

Proprietary **Omni Tip<sup>™</sup> Plastic Probes** are available in both soft tip (deaggregation) and hard tissue (frozen) versions. Engineered from durable plastics,

Omni Tips<sup>™</sup> can be cleaned and reused, yet are economical enough that they can be disposed of after each use to eliminate cross contamination.

#### A selection of our best-selling plastic probes:

SOFT TISSUE PLASTIC GENERATOR PROBES		ORDER NO.
Omni Tip™ Plastic Generator Probes, Soft Tissue	25 pack	30750
Omni Tip™ Plastic Generator Probes, Soft Tissue	50 pack	32750
Omni Tip™ Plastic Generator Probes, Soft Tissue	100 pack	34750
HARD TISSUE PLASTIC GENERATOR PROBES ORDER NO.		
Omni Tip™ Plastic Generator Probes, Hard Tissue	25 pack	30750H
Omni Tip™ Plastic Generator Probes, Hard Tissue	50 pack	32750H
Omni Tip™ Plastic Generator Probes, Hard Tissue		

#### A 家 和 生 国 統 一 服 务 热 线 4009-670-570

#### 奥然科技 上海

地址:上海市徐汇区漕溪路250号银海大厦A1012室 电话:021-6451 6293 6451 6297 6451 6299

#### 奥然科技 北京

地址:北京市海淀区北三环西路甲30号双天大厦423室 电话:010-6840 0176 6840 0177

#### 奥然科技 广州

地址:广州市天河区五山路141号尚德大厦A栋1707室 电话:020-8551 6682 8551 6679

#### 奥然科技 南京

地址:南京市中央路417号先锋广场1036室 电话:025-5850 1866 5850 1966